


Little Wave XP Aluminum Seat Pan - Growth



Part Number	Tape Size	Seat Pan Depth
004458	3" X 4"	10" - 11"
004458 & 004459	3" X 4" & 3" X 5"	12"
004459	3" X 5"	13"
004460	3" X 6"	14" - 19"

Screw (Item 2) Qty Per Chair	
Seat Depth	Screw Qty
10" - 12"	6
13" - 15"	8
16" - 18"	10
19"	12

Rigid Seat Pans may be quoted, but must be ordered through Customer Service. If your chair was ordered prior to February 2013 a custom seat pan may be required. The seat pan width measures approximately 0.75" and the seat pan depth measures approximately 0.50" less than the chair size.

Item Number	Part Number	Description	Quantity
1a, 3a	505952	Seat Pan Rigid with AB Tape 10 x 10	1
1b, 3a	505953	Seat Pan Rigid with AB Tape 10 x 11	1
1c, 3a, 3b	505954	Seat Pan Rigid with AB Tape 10 x 12	1

Item Number	Part Number	Description	Quantity
1d, 3b	505955	Seat Pan Rigid with AB Tape 10 x 13	1
1e, 3c	505956	Seat Pan Rigid with AB Tape 10 x 14	1
1f, 3c	505957	Seat Pan Rigid with AB Tape 10 x 15	1
1g, 3c	505958	Seat Pan Rigid with AB Tape 10 x 16	1
1h, 3c	505959	Seat Pan Rigid with AB Tape 10 x 17	1
1i, 3c	505960	Seat Pan Rigid with AB Tape 10 x 18	1
1j, 3c	505961	Seat Pan Rigid with AB Tape 10 x 19	1
1k, 3a	505962	Seat Pan Rigid with AB Tape 11 x 10	1
1l, 3a	505963	Seat Pan Rigid with AB Tape 11 x 11	1
1m, 3a, 3b	505964	Seat Pan Rigid with AB Tape 11 x 12	1
1n, 3b	505965	Seat Pan Rigid with AB Tape 11 x 13	1
1o, 3c	505966	Seat Pan Rigid with AB Tape 11 x 14	1
1p, 3c	505967	Seat Pan Rigid with AB Tape 11 x 15	1
1q, 3c	505968	Seat Pan Rigid with AB Tape 11 x 16	1
1r, 3c	505969	Seat Pan Rigid with AB Tape 11 x 17	1
1s, 3c	505970	Seat Pan Rigid with AB Tape 11 x 18	1
1t, 3c	505971	Seat Pan Rigid with AB Tape 11 x 19	1
1u, 3a	505972	Seat Pan Rigid with AB Tape 12 x 10	1
1v, 3a	505973	Seat Pan Rigid with AB Tape 12 x 11	1
1w, 3a, 3b	505974	Seat Pan Rigid with AB Tape 12 x 12	1
1x, 3b	505975	Seat Pan Rigid with AB Tape 12 x 13	1
1y, 3c	505976	Seat Pan Rigid with AB Tape 12 x 14	1
1z, 3c	505977	Seat Pan Rigid with AB Tape 12 x 15	1
1aa, 3c	505978	Seat Pan Rigid with AB Tape 12 x 16	1
1ab, 3c	505979	Seat Pan Rigid with AB Tape 12 x 17	1
1ac, 3c	505980	Seat Pan Rigid with AB Tape 12 x 18	1
1ad, 3c	505981	Seat Pan Rigid with AB Tape 12 x 19	1
1ae, 3a	505983	Seat Pan Rigid with AB Tape 13 x 10	1

Item Number	Part Number	Description	Quantity
1af, 3a	505984	Seat Pan Rigid with AB Tape 13 x 11	1
1ag, 3a, 3b	505985	Seat Pan Rigid with AB Tape 13 x 12	1
1ah, 3b	505986	Seat Pan Rigid with AB Tape 13 x 13	1
1ai, 3c	505987	Seat Pan Rigid with AB Tape 13 x 14	1
1aj, 3c	505988	Seat Pan Rigid with AB Tape 13 x 15	1
1ak, 3c	505989	Seat Pan Rigid with AB Tape 13 x 16	1
1al, 3c	505990	Seat Pan Rigid with AB Tape 13 x 17	1
1am, 3c	505991	Seat Pan Rigid with AB Tape 13 x 18	1
1an, 3c	505992	Seat Pan Rigid with AB Tape 13 x 19	1
1ao, 3a	505994	Seat Pan Rigid with AB Tape 14 x 10	1
1ap, 3a	505995	Seat Pan Rigid with AB Tape 14 x 11	1
1aq, 3a, 3b	505996	Seat Pan Rigid with AB Tape 14 x 12	1
1ar, 3b	505997	Seat Pan Rigid with AB Tape 14 x 13	1
1as, 3c	505859	Seat Pan Rigid with AB Tape 14 x 14	1
1at, 3c	505860	Seat Pan Rigid with AB Tape 14 x 15	1
1au, 3c	505861	Seat Pan Rigid with AB Tape 14 x 16	1
1av, 3c	505862	Seat Pan Rigid with AB Tape 14 x 17	1
1aw, 3c	505863	Seat Pan Rigid with AB Tape 14 x 18	1
1ax, 3c	505864	Seat Pan Rigid with AB Tape 14 x 19	1
1ay, 3a	505998	Seat Pan Rigid with AB Tape 15 x 10	1
1az, 3a	505999	Seat Pan Rigid with AB Tape 15 x 11	1
1ba, 3a, 3b	506000	Seat Pan Rigid with AB Tape 15 x 12	1
1bb, 3b	506001	Seat Pan Rigid with AB Tape 15 x 13	1
1bc, 3c	505868	Seat Pan Rigid with AB Tape 15 x 14	1
1bd, 3c	505869	Seat Pan Rigid with AB Tape 15 x 15	1
1be, 3c	505870	Seat Pan Rigid with AB Tape 15 x 16	1
1bf, 3c	505871	Seat Pan Rigid with AB Tape 15 x 17	1
1bg, 3c	505872	Seat Pan Rigid with AB Tape 15 x 18	1

Item Number	Part Number	Description	Quantity
1bh, 3c	505873	Seat Pan Rigid with AB Tape 15 x 19	1
1bi, 3a	506002	Seat Pan Rigid with AB Tape 16 x 10	1
1bj, 3a	506003	Seat Pan Rigid with AB Tape 16 x 11	1
1bk, 3a, 3b	506004	Seat Pan Rigid with AB Tape 16 x 12	1
1bl, 3b	506005	Seat Pan Rigid with AB Tape 16 x 13	1
1bm, 3c	505877	Seat Pan Rigid with AB Tape 16 x 14	1
1bn, 3c	505878	Seat Pan Rigid with AB Tape 16 x 15	1
1bo, 3c	505879	Seat Pan Rigid with AB Tape 16 x 16	1
1bp, 3c	505880	Seat Pan Rigid with AB Tape 16 x 17	1
1bq, 3c	505881	Seat Pan Rigid with AB Tape 16 x 18	1
1br, 3c	505882	Seat Pan Rigid with AB Tape 16 x 19	1
1bs, 3a, 3b	506006	Seat Pan Rigid with AB Tape 17 x 12	1
1bt, 3b	506007	Seat Pan Rigid with AB Tape 17 x 13	1
1bu, 3c	505886	Seat Pan Rigid with AB Tape 17 x 14	1
1bv, 3c	505887	Seat Pan Rigid with AB Tape 17 x 15	1
1bw, 3c	505888	Seat Pan Rigid with AB Tape 17 x 16	1
1bx, 3c	505889	Seat Pan Rigid with AB Tape 17 x 17	1
1by, 3c	505890	Seat Pan Rigid with AB Tape 17 x 18	1
1bz, 3c	505891	Seat Pan Rigid with AB Tape 17 x 19	1
1ca, 3b	506008	Seat Pan Rigid with AB Tape 18 x 13	1
1cb, 3c	505895	Seat Pan Rigid with AB Tape 18 x 14	1
1cc, 3c	505896	Seat Pan Rigid with AB Tape 18 x 15	1
1cd, 3c	505897	Seat Pan Rigid with AB Tape 18 x 16	1
1ce, 3c	505898	Seat Pan Rigid with AB Tape 18 x 17	1
1cf, 3c	505899	Seat Pan Rigid with AB Tape 18 x 18	1
1cg, 3c	505900	Seat Pan Rigid with AB Tape 18 x 19	1
1ch, 3c	505901	Seat Pan Rigid with AB Tape 18 x 20	1
1ci, 3c	505902	Seat Pan Rigid with AB Tape 18 x 21	1

Item Number	Part Number	Description	Quantity
1cj, 3c	505903	Seat Pan Rigid with AB Tape 18 x 22	1
1ck, 3c	505904	Seat Pan Rigid with AB Tape 19 x 14	1
1cl, 3c	505905	Seat Pan Rigid with AB Tape 19 x 15	1
1cm, 3c	505906	Seat Pan Rigid with AB Tape 19 x 16	1
1cn, 3c	505907	Seat Pan Rigid with AB Tape 19 x 17	1
1co, 3c	505908	Seat Pan Rigid with AB Tape 19 x 18	1
1cp, 3c	505909	Seat Pan Rigid with AB Tape 19 x 19	1
2	100734	M5 x 16 PHMS ZC	12
3a	004458	AB Tape, 3" x 4"	4
3b	004459	AB Tape, 3" x 5"	4
3c	004460	AB Tape, 3" x 6"	4