



ORDER FORM

To see QUICKIE Power Options
or Accessories, scan here or visit:
[https://www.sunrisemedical.com/
power-wheelchair-options-accessories](https://www.sunrisemedical.com/power-wheelchair-options-accessories)



Q700M


with SEDEO® PRO™, PRO™ ADVANCED
and JAY® Syncra™ Seating

EIPW25

Customer	Ship To
Account #	Name
Date	Address
PO#	City
Buyer	State Zip
Mark For	Phone

Client Measurements

a) Top of head to seat	in
b) Top of shoulder to seat	in
c) Underarm to seat	in
d) Elbow to seat	in
e) Top of leg to seat	in
f) Shoulder width	in
g) Chest width	in
i) Hip width	in
j) Actual sitting depth	in
m) Foot length	in
n) Knee to foot	in
o) Seat width	in



Notes

Chair Specifications

HCPC	K0848 / K0856 / K0861 / K0868 / K0877 / K0884
Typical Weight w/out Batteries	330 lbs.
Typical Weight with Batteries	410 lbs.
Power Seating	50° of CG Tilt and Power Recline with 4" mechanical sheer reduction / 12" seat elevator
Seat Type	SEDEO® PRO™ / JAY® Syncra™
Seat Size Range	Width 16" to 22" Depth 14" to 22"
Seat-to-Floor Height Range	16.5" to 19"

Mid-Wheel Drive Base Specifications

Drive Wheel Position	Mid-wheel drive
User Weight Capacity	300 lbs.
Base Size	Width 24.5", Length 36.5"
Minimum Turning Radius	19.8" (excludes footplates)
Drive Control	R-net PM 120A / R-net PM 90A
Driving Range	Up to 18 miles (w/ Group 24 Batteries)
Drive Motors	4-pole 6 mph, 8 mph
Drive Wheels	14" w/flat-free inserts (black)
Casters	6" front / 6" rear (black)
Charger	10A Dual Mode (Compact)
Battery Size	Group 24 / Group 34
Color	9 shroud colors and 7 color accents
Maximum Ground Clearance	3"

BASE			
OPTION	HCPC	DESCRIPTION	RETAIL
Group 3 Standard Duty Mid Wheel Drive Base End Item Number EIPW25 (PRO™ Seat)			\$
250PC197	K0848	Quickie Q700 M NPC No Power Option for SEDEO® PRO™ Seat Included in Base Pricing : Q700 M Base, R-net 90amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5".	\$8,979.00
250PC198	K0856	Quickie Q700 M SPC Single Power Option for SEDEO® PRO™ Seat Included in Base Pricing : Q700 M Base, R-net 90amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5".	\$8,979.00
250PC203	K0861	Quickie Q700 M MPC Tilt/Recline for SEDEO® PRO™ Seat Included in Base Pricing: Q700 M Base, R-net 90amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5".	\$8,979.00
Group 4 Standard Duty Mid Wheel Drive Base End Item Number EIPW25 (PRO™ Seat)			\$
250PC201	K0868	Quickie Q700 M NPD No Power Option for SEDEO® PRO™ Seat Included in Base Pricing : Q700 M Base, R-net 120amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 8mph, Weight Capacity = 300lbs., STFH = 16.5".	\$14,297.00
250PC202	K0877	Quickie Q700 M SPD Single Power Option for SEDEO® PRO™ Seat Included in Base Pricing : Q700 M Base, R-net 120amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 8mph, Weight Capacity = 300lbs., STFH = 16.5".	\$14,297.00
250PC205	K0884	Quickie Q700 M MPD Tilt/Recline for SEDEO® PRO™ Seat Included in Base Pricing: Q700 M Base, R-net 120amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 8mph, Weight Capacity = 300lbs., STFH = 16.5".	\$14,297.00

BASE			
OPTION	HCPC	DESCRIPTION	RETAIL
User Weight			\$
250WG23	—	Up to 300 lbs.	NC
Motor Package			\$
250PF99	—	4-pole 10 kph / 6 mph <i>(Standard for Group 3 K0848, K0856, K0861)</i>	NC
250PF100	—	4-pole 12.5 kph / 8 mph <i>(Standard for Group 4 K0868, K0877, K0884)</i>	NC
250PF101	—	4-pole 12.5 kph / 8 mph <i>(Upgrade for Group 3 K0848, K0856, K0861)</i>	\$1,316.00
Tracking			\$
250MI01	K0108	SureTrac® Module & Motor Encoders <i>(Improves drive tracking)</i>	\$1,596.00
Drive Wheels			\$
250DW58	—	Drive Wheel, 14", Dual Density, Knobby, Black	NC
Caster Wheels			\$
250CW55	—	Caster Wheels 6", Solid, Black	NC
Shroud Color			\$
250C119	—	Matte Black	NC
250C165	—	Candy Apple Red	NC
250C166	—	Candy Blue	NC
250C167	—	Candy Purple	NC
250C176	—	Green Apple	NC
250C178	—	Pearl Pink	NC
250C182	—	Silver	NC
250C125	—	Pearl White	NC
250C170	—	American Flag	\$590.00
Color Accents			\$
<i>Note: Color Accents and Shroud Colors may not match.</i>			
250EZ20	—	Green	NC
250EZ22	—	Black	NC
250EZ21	—	Blue	NC
250EZ23	—	Orange	NC
250EZ24	—	Red	NC
250EZ25	—	Silver	NC
250EZ19	—	White	NC

BASE							
OPTION	HCPC	DESCRIPTION		RETAIL			
Service Decals				\$			
250D10	—	Army		NC			
250D11	—	Air Force		NC			
250D13	—	Marines		NC			
250D14	—	Navy		NC			
250D15	—	Purple Heart		NC			
ISO Transit Tie Downs				\$			
250T050	—	ISO Transit Package		NC			
LED Lights				\$			
250LI04	—	LED Lights and Indicators		\$1,373.00			
250LI08	—	LED Fender Lights		\$555.00			
250LI09	—	LED Fender Lights (Available with Multi-Power only))		NC			
250LI02	—	Omit LED Lights		NC			
Seat Heights				\$			
<i>Note: Available seat-to-floor heights are dependent on 2 options: 1- Battery Size and 2- Powered Seating. Use the matrix below to determine available ranges for your desired options. Measured from front of seat pan to floor with 0° frame angle.</i>							
		Seat Height	Power Seating	No Power Seating	PRO ADVANCED Recline		
					w/o Elevate	w/ Elevate	
250HT103	—	16.5"	Group 34	Group 34	N/A	N/A	NC
250HT55	—	17"	Group 34	GP24 or GP34	GP34	GP34	NC
250HT104	—	17.5"	Group 34	GP24 or GP34	GP34	GP34	NC
250HT56	—	18"	GP24 or GP34	GP24 or GP34	GP24 / GP34	GP34	NC
250HT105	—	18.5"	GP24 or GP34	GP24 or GP34	GP24 / GP34	GP24 / GP34	NC
250HT57	—	19"	GP24 or GP34	GP24 or GP34	GP24 / GP34	GP24 / GP34	NC
Batteries and Chargers				\$			
Group 24 Battery Options <i>(Recommended for max daily range)</i>							
250AG25	E2363	GP24 Flat Top Pair w/Strap MK Gel				\$1,116.00	
250AG26	E2363	GP24 Flat Top Pair w/Strap Brightway AGM				\$840.00	
Group 34 Battery Options <i>(Required for lower seat-to-floor heights)</i>							
250AG27	E2359	GP34 Flat Top Pair w/Straps MK Gel				\$1,075.00	
250AG29	E2359	GP34 Flat Top Pair w/Straps Brightway AGM				\$812.00	
Omit Batteries							
250AG24	—	Battery Harnesses Only <i>(Omit Batteries)</i>				NC	
Chargers							
250BJ28	—	10A Dual Mode (Compact)				NC	

POWERED SEATING			
OPTION	HCPC	DESCRIPTION	RETAIL
Power Seating Options			\$
Seat Elevator			
250PI08	E2300	C-Me® Seat Elevator, Powered, 300 mm / 12" <i>(C-Me speed up to 3 mph)</i>	\$3,717.00
250PI07	E2300	Seat Elevator, Powered, 12" <i>(Required for Anterior Tilt)</i>	\$3,717.00
Seat Tilt			
250PA35	E1002	Posterior 50° Tilt	\$8,292.00
Seat Recline			
250PA37	E1004	150° Powered Recline w/ 4" Shear Reduction	\$8,292.00
Seat Tilt/Recline			
250PA41	E1007	Posterior 50° Tilt and 150° Recline Powered, w/ 4" Shear Reduction	\$13,440.00
250PA62	E1007	Posterior 50° Tilt and PRO ADVANCED 165° Recline Powered, w/ 7" Shear Reduction	\$13,745.00
Seat Anterior Tilt			
250TG02	K0108	Anterior Tilt 0° to 30°	\$2,963.00
Power Recline Backrest Preset Home Angle			
250QK1	—	90°	NC
250QK2	—	95° (Default Position)	NC
Pre Set Seat Frame Angle <i>(No Power Seating Options, Increases Seat to Floor Height)</i>			\$
250FA10	—	0° <i>(Adds 0" to STFH)</i>	NC
250FA3	—	3° <i>(Adds 1" to STFH)</i>	NC
Backrest Frame Type <i>(N/A with Powered Recline)</i>			\$
250BF10	—	Tool Angle Adjustable Back, 90°-120°	NC
250BF11	E1225	User Adjustable Manual Semi-Recline Back, 90°-130°, Right <i>(Adjustment lever mounted on opposite side of joystick)</i>	\$669.00
250BF12	E1225	User Adjustable Manual Semi-Recline Back, 90°-130°, Left <i>(Adjustment lever mounted on opposite side of joystick)</i>	\$669.00
Backrest Frame Option <i>(N/A with PRO ADVANCED)</i>			\$
250BF13	—	Flip Down, Quick Release	\$144.00
Back Cushion Type			\$
250BT134	—	SEDEO® Back Cushion <i>(Choose Seating in Section A) (N/A with PRO ADVANCED)</i>	NC
250BT138	—	JAY® Syncra™ <i>(Choose Seating in Section B)</i>	NC
250BT135	—	JAY® J3® <i>(Choose Seating in Section C)</i>	NC
250BT137	—	SEDEO® Aftermarket Back Adapter <i>(Height adjustable round tubes for traditional mounting)</i>	\$515.00
250BT133	—	OMIT <i>(For personalized mounting. No headrest mount)</i>	NC

POWERED SEATING

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

A - SEDEO® Backrest Dimensions (N/A with PRO ADVANCED)

\$

Back Width

250BN28	—	400mm / 16"	NC
250BN29	—	460mm / 18"	NC
250BN30	—	500mm / 20"	NC
250BN31	—	560mm / 22"	NC

Back Shell Height

			Tool Angle Adjust.	150° Recline	Installed Back Height	
250B255	—	460mm / 18"	22°-23°	23°-25°		NC
250B256	—	560mm / 22"	24°-27°	26°-29°		NC
250B257	—	660mm / 26"	28°-30°	30°		NC

Installed Back Height - (From Pan to Top of Back)

250JR32	—	22"	NC
250JR33	—	23"	NC
250JR34	—	24"	NC
250JR35	—	25"	NC
250JR36	—	26"	NC
250JR37	—	27"	NC
250JR38	—	28"	NC
250JR39	—	29"	NC
250JR40	—	30"	NC

POWERED SEATING			
OPTION	HCPC	DESCRIPTION	RETAIL
B - JAY® Syncra™ Seating			
B - JAY® Syncra™ Back  P. 22			\$
250JP14	E2620	JAY® Syncra™ Back 16-20" Wide with 3" Contour Depth	\$1,116.00
250JP15	E2621	JAY® Syncra™ Back 22" Wide with 3" Contour Depth	\$1,230.00
B - JAY® Syncra™ Back Width			\$
<i>Back can be ordered a different width than seat pan. Armrest will be adjusted to back width ordered.</i>			
250JQ07	—	16" Wide (Adjustment range 16 to 17" Wide)*	NC
250JQ08	—	18" Wide (Adjustment range 18 to 19" Wide)*	NC
250JQ09	—	20" Wide (Adjustment range 20 to 21" Wide)*	NC
250JQ10	—	22" Wide (Adjustment range 22 to 23" Wide)	NC
<i>*Additional +1" of shell width adjustment is achievable. Requires purchase of 1" larger foam and upholstery separately.</i>			
B - JAY® Syncra™ Back Shell Options and Height Ranges			\$
<i>JAY® Syncra™ Back shells are mounted with a minimim 5.5" gap between the bottom of the shell to the seat pan to accommodate the mechanical shear. Back height can be adjusted in the field within the ranges shown based on the back frame type.</i>			
JAY® Syncra™ Shell Height			
250JR42	—	16" Tall Shell (22" to 24" Back Height)	NC
250JR43	—	18" Tall Shell (24" to 26" Back Height)	NC
250JR44	—	20" Tall Shell (26" to 28" Back Height)	NC
250JR45	—	22" Tall Shell (28" to 30" Back Height)	NC
JAY® Syncra™ Back Installed Back Height (From Pan to Top of Back)			
250JR50	—	22" (16" Shell Height)	NC
250JR51	—	23" (16" Shell Height)	NC
250JR52	—	24" (16" Shell Height)	NC
250JR53	—	24" (18" Shell Height)	NC
250JR54	—	25" (18" Shell Height)	NC
250JR55	—	26" (18" Shell Height)	NC
250JR56	—	26" (20" Shell Height)	NC
250JR57	—	27" (20" Shell Height)	NC
250JR58	—	28" (20" Shell Height)	NC
250JR67	—	28" (22" Shell Height)	NC
250JR68	—	29" (22" Shell Height)	NC
250JR69	—	30" (22" Shell Height)	NC
B - JAY® Syncra™ Back Upholstery			\$
250BU96	—	Syncra™ Premium Leatherette	NC
250BU97	—	Syncra™ Stretch	NC

POWERED SEATING

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

C - JAY® J3® Back

C - JAY® J3® Back - Contour

\$

250JP01	E2613	J3® PA Back for 16"-20" Wide	\$613.00
250JP02	E2614	J3® PA Back for 22" Wide	\$910.00
250JP11	E2615	J3® PL Back for 16"-20" Wide	\$715.00
250JP13	E2620	J3® PD Back for 16"-20" Wide (18" and 21" shell height only)	\$858.00

C - JAY® J3® Back Width

\$

Back can be ordered a different width than seat pan. Armrest will be adjusted to back width ordered.

250JQ02	—	16"	NC
250JQ03	—	18"	NC
250JQ04	—	20"	NC
250JQ05	—	22"	NC

C - JAY® J3® Back Shell Options and Height Ranges

\$

JAY® J3® shells are mounted with a minimim 4" gap between the bottom of the shell to the seat pan to accommodate the mechanical shear. Back height can be adjusted in the field within the ranges shown, based on the back frame type and the amount of usable hardware adjustability.

JAY® J3® Shell Height

Tool Angle Adjust.

150° Recline

PRO ADVANCED

250JR41	—	15" tall shell (Upper Thoracic Med) for 22" Wide	22"-23"	22"-23"	21"-24"	Installed Back Height	NC
250JR14	—	16.5" tall shell (Upper Thoracic Short)	21"-23"	22"-24"	21"-24"		NC
250JR15	—	18" tall shell (Upper Thoracic Med)	22"-24"	23"-25"	22"-26"		NC
250JR16	—	19.5" tall shell (Upper Thoracic Tall)	24"-26"	24"-26"	23"-27"		NC
250JR18	—	18" tall shell (Shoulder Height) for 22" wide	24"	24"	21"-25"		NC
250JR17	—	21" tall shell (Shoulder Height)	24"-26"	25"-28"	25"-28"		NC
250JR27	—	22.5" tall shell (Shoulder Height)	26"-28"	27"-30"	26"-29"		NC
250JR28	—	24" tall shell (Shoulder Height)	27"-29"	28"-30"	28"-30"		NC

JAY® J3® Installed Back Height (From Pan to Top of Back)

250JR29	—	21"	NC
250JR19	—	22"	NC
250JR20	—	23"	NC
250JR21	—	24"	NC
250JR22	—	25"	NC
250JR23	—	26"	NC
250JR24	—	27"	NC
250JR25	—	28"	NC
250JR30	—	29"	NC
250JR31	—	30"	NC

POWERED SEATING			
OPTION	HCPC	DESCRIPTION	RETAIL
C - JAY® Anterior Thoracic Support			\$
<i>Classic-Structure Pull Down Harness (4 point strap, only available with JAY® backrest)</i>			
250K10	E0960	Medium (ATS8302)	\$174.00
250K11	E0960	Large (ATS8303)	\$174.00
JAY® Lateral Thoracic Swing-Away Supports (Requires Pad, Upholstery and Mount)			\$
JAY® Lateral Thoracic Support Pads			
250JW35	E0956	Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Right	\$144.00
250JW38	E0956	Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Left	\$144.00
250JW36	E0956	Lateral Thoracic (L x H) 100 mm x 152 mm / 4" x 6" Pad - Right	\$144.00
250JW39	E0956	Lateral Thoracic (L x H) 100 mm x 52 mm / 4" x 6" Pad - Left	\$144.00
250JW37	E0956	Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Right	\$144.00
250JW40	E0956	Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Left	\$144.00
JAY® Lateral Thoracic Support Pad Upholstery			
250NU01	—	Syncra™ Premium Leatherette	NC
250NU02	—	Syncra™ Stretch	NC
JAY® Lateral Thoracic Swing-Away Mount			
250JW41	E1028	Lateral Thoracic Swing Away Bracket - Right	\$286.00
250JW42	E1028	Lateral Thoracic Swing Away Bracket - Left	\$286.00
JAY® Syncra™ Fixed Lateral Support Accessory			\$
JAY® Fixed Lateral Support Accessory Hardware (JAY® Fixed Lateral Support is height adjustable)			
250JW33	E0956	Fixed Lateral Support Accessory (L x H) 90 mm x 140 mm / 3.5" x 8" - Right	\$144.00
250JW34	E0956	Fixed Lateral Support Accessory (L x H) 90 mm x 140 mm / 3.5" x 8" - Left	\$144.00
Headrests (Requires Pad, Upholstery, & Mount)			\$
Headrest Pad			
250HE36	E0955	SEDEO® Pad Small With Cover (LxH) 240 x 150 mm / 10" x 6" <i>Can only be used with Multi Position Headrest Mount with Quick Release and SEDEO® Upholstery</i>	\$184.00
250HE38	E0955	Whitmyer® Plush Headrest Pad 203 mm / 8"	\$234.00
250HE39	E0955	Whitmyer® Plush Headrest Pad 254 mm / 10"	\$234.00
250HE44	E0955	Whitmyer® Plush Headrest Pad 356 mm / 14"	\$245.00
250HE46	E0955	Whitmyer® Adjust-A-Plush Headrest Pad Medium (Adjustable from 7 to 16" Wide x 4" Tall)	\$360.00
250HE48	—	OMIT Headrest Pad	NC
Headrest Pad Upholstery			
250NU10	—	SEDEO® Upholstery	NC
250NU04	—	Syncra™ Premium Leatherette	NC
250NU05	—	Syncra™ Stretch	NC
250NU11	—	Whitmyer® Lycra	NC

POWERED SEATING			
OPTION	HCPC	DESCRIPTION	RETAIL
Headrest Mount			
250MH32	E1028	SEDEO® Multi-Position Headrest Mount with Quick Release <i>(Only available with SEDEO® Pad)</i>	\$229.00
250MH33	E1028	Whitmyer® Linx Headrest Mount with Quick Release	\$188.00
250MH34	E1028	Whitmyer® Linx-2 Headrest Mount with Quick Release	\$229.00
250MH35	E1028	Whitmyer® AXYS Headrest Mount with Quick Release	\$121.00
250MH17	—	OMIT Headrest Mount	NC
Push Handles			
250PH27	—	Height Adjustable Push Handles	\$343.00
Seat Frame Rails			
250SN14	—	Seat Rail Short <i>(For effective seat depth 14" to 18")</i>	NC
250SN15	—	Seat Rail Long <i>(For effective seat depth 17" to 22")</i>	NC
Seat Cushion Type			
250ST79	—	SEDEO® (See Section D) <i>(N/A with PRO ADVANCED)</i>	NC
250ST81	—	JAY® Syncra™ Cushion (See Section E)	NC
250ST80	—	JAY® (See Section F)	NC
250ST75	—	Omit	NC
C - Seat Size			
Seat Pan Width <i>(Odd sizes not available on SEDEO® cushions)</i>			
250R53	—	16"	NC
250R54	—	17"	NC
250R55	—	18"	NC
250R56	—	19"	NC
250R57	—	20"	NC
250R58	—	21" <i>(Not available with JAY® Syncra™ Cushions)</i>	NC
250R59	—	22"	NC
Seat Pan Depth <i>(Front of seat pan to front of back cushion – effective depth)</i>			
250RF52	—	14" <i>(N/A with PRO ADVANCED)</i>	NC
250RF53	—	15" <i>(N/A with PRO ADVANCED)</i>	NC
250RF51	—	16" <i>(N/A with PRO ADVANCED)</i>	NC
250RF45	—	17"	NC
250RF46	—	18"	NC
250RF47	—	19"	NC
250RF48	—	20"	NC
250RF49	—	21" <i>(Not available with JAY® Syncra™ Cushions)</i>	NC
250RF50	—	22"	NC

POWERED SEATING

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

D - SEDEO® Cushions (N/A with PRO ADVANCED) \$

250CU82	E2607	SEDEO® Cushion (Sizes 16", 18", & 20" Wide) (Not available with 14" Depth)	\$451.00
250CU83	E2608	SEDEO® Cushion (Size 22" Wide)	\$567.00

E - JAY® Syncra™ Cushion P. 22 (JAY Your Way modifications available through Customer Service (800) 333-4000) \$

250CU89	E2622	JAY® Syncra™ Cushion featuring JAY Flow® Fluid with Cryo™ Technology (16" to 20" wide)	\$618.00
250CU90	E2623	JAY® Syncra™ Cushion featuring JAY Flow® Fluid with Cryo™ Technology (22" wide)	\$715.00

JAY® Syncra™ Cushion Cover Upholstery

250YR73	—	Syncra™ Premium Leatherette	NC
250YR74	—	Syncra™ Stretch	NC

F - JAY® Seat Cushion \$

250CU59	E2607	JAY® Ion® (Sizes 16" - 21" Wide)	\$451.00
250CU60	E2608	JAY® Ion® (Size 22" Wide)	\$567.00
250CU61	E2607	JAY® Union® (Sizes 16" - 21" Wide)	\$510.00
250CU62	E2608	JAY® Union® (Size 22" Wide)	\$623.00
250CU64	E2622	JAY® Fusion® (Sizes 16" - 21" Wide)	\$504.00
250CU65	E2623	JAY® Fusion® (Size 22" Wide)	\$601.00

F - JAY® Solid Seat \$

250CU84	—	JAY® Solid Seat	\$98.00
---------	---	-----------------	---------

POSITIONING SUPPORTS & ACCESSORIES

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

Knee and Thigh Swing-Away Support Pads \$

Lateral Knee Support Pads and Hardware *(must pick pad and mount)*

250XY05	E0953	Lateral Knee Support Pad, Right (4" round)	\$201.00
250XY06	E0953	Lateral Knee Support Pad, Left (4" round)	\$201.00
250XX08	E1028	Lateral Swing-Away Bracket Knee - Right	\$263.00
250XX09	E1028	Lateral Swing-Away Bracket Knee - Left	\$263.00

Lateral Thigh Support Pads and Hardware *(must pick pad and mount)*

250KF22	E0953	Lateral Thigh Support Pad Short - Right (5"x5")	\$144.00
250KF23	E0953	Lateral Thigh Support Pad Short - Left (5"x5")	\$144.00
250KF24	E0953	Lateral Thigh Support Pad Long - Right (5"x7")	\$144.00
250KF25	E0953	Lateral Thigh Support Pad Long - Left (5"x7")	\$144.00
250KF26	E1028	Lateral Thigh Removable Mount - Right	\$286.00
250KF27	E1028	Lateral Thigh Removable Mount - Left	\$286.00

JAY® Lateral Thigh Supports (Must pick pad, upholstery, and mount) \$

JAY® Lateral Thigh Support Pads

250KF16	E0953	Lateral Thigh Support Pad Standard, Right (8.25" x 3.25")	\$144.00
250KF17	E0953	Lateral Thigh Support Pad Standard, Left (8.25" x 3.25")	\$144.00
250KF18	E0953	Lateral Thigh Support Pad Large, Right (12" x 3.25")	\$144.00
250KF19	E0953	Lateral Thigh Support Pad Large, Left (12" x 3.25")	\$144.00

JAY® Lateral Thigh Support Pad Upholstery

250NU07	—	Syncra™ Premium Leatherette	NC
250NU08	—	Syncra™ Stretch	NC

JAY® Lateral Thigh Removable Mount

250KF20	E1028	Lateral Thigh Removable Mount - Right	\$286.00
250KF21	E1028	Lateral Thigh Removable Mount - Left	\$286.00

Armrest Structure \$

250EG35	—	Reclining, Flip Back, Long (Height adjustable 9" to 15")	NC
250EG41	—	Reclining, Flip Back, Short (Best used with shorter seat depths) (Height adjustable 9" to 15")	NC
250EG36	—	Reclining, Flip Back with 2-point post locking arm (Height adjustable 11" to 15" x .625")	\$286.00
250EG51	—	Reclining, Flip Back, Dual Link Armrest (Required for PRO ADVANCED) (Height adjustable 8" to 15" including length and lateral adjustment)	NC
250EG32	—	Single Post, Flip Back (N/A with thigh support), (Height adjustable 10" to 15")	NC

Arm Pads \$

SEDEO® Armrest Pad

250AB43	—	300 mm / 12", Integral Foam, Right	NC
250AE43	—	300 mm / 12", Integral Foam, Left	NC
250AB44	—	400 mm / 16", Integral Foam, Right	NC
250AE44	—	400 mm / 16", Integral Foam, Left	NC

Comfort Arm Pad (Foam)

250AB46	—	2" Wide Right	\$86.00
250AE46	—	2" Wide Left	\$86.00

POSITIONING SUPPORTS & ACCESSORIES

OPTION	HCPC	DESCRIPTION	RETAIL
Comfort Arm Pad (Gel)			
250AB48	—	2" Wide Right	\$172.00
250AE48	—	2" Wide Left	\$172.00
Ottobock® Arm Troughs			
250AE06	—	Ottobock® Arm Trough, Med (12"), Left	\$133.00
250AB06	—	Ottobock® Arm Trough, Med (12"), Right	\$133.00
250AE07	—	Ottobock® Arm Trough, Large (14") Left,	\$133.00
250AB07	—	Ottobock® Arm Trough, Large (14"), Right	\$133.00
250OL01	—	Ottobock® Hand Pad, Med (7"), Left	\$58.00
250OR01	—	Ottobock® Hand Pad, Med (7"), Right	\$58.00
250OL02	—	Ottobock® Hand Pad, Large (9"), Left	\$58.00
250OR02	—	Ottobock® Hand Pad, Large (9"), Right	\$58.00
250AE78	—	Ottobock® One Piece Arm Trough, (18"), Left	\$186.00
250AB78	—	Ottobock® One Piece Arm Trough, (18"), Right	\$186.00
Side Guards (Pair)			\$
250SG44	E0953	Side Guards Integral Foam (11"L x 5.5"H) (N/A with Reclining Armrests)	\$234.00
Fixed Centermount Legrest			\$
250H183	—	SEDEO® Fixed Centermount Legrest (3" to 17.5" Leg Length)	\$286.00
250H184	—	Fixed Centermount Legrest, Short (Includes Calf Pads & Single Footplate) (11.5" to 14.5" Leg Length)	NC
250H185	—	Fixed Centermount Legrest, Long (Includes Calf Pads & Single Footplate) (13.5" to 18.5" Leg Length)	NC
Power Centermount Legrest			\$
250H175	E1012	Power Centermount Legrest, 90° - 0° Angle, 7.5" Extension (Mechanical Articulation) (5" to 17.5" Leg Length)	\$2,653.00
250H109	E1012	Power Centermount Legrest, 90° - 0° Angle, 7.5" Extension (Powered Articulation) (3" to 17.5" Leg Length) (Required for PRO ADVANCED)	\$3,717.00
Power Centermount Legrest Accessories			\$
Footplates			
250F127	K0040	Dual, Height, Width and Angle Adjustable, Standard (14"W x 10"D)	\$193.00
250F128	K0040	Dual, Height, Width and Angle Adjustable, Large (16"W x 12"D)	\$193.00
Calf Pad Options			
250DI34	—	Calf Pads, Pair, Standard	NC
250DI35	E0995	Extra Large Pair Calf Pads for Centermount Legrests. Fully adjustable (N/A w/ Std footplate)	\$400.00
Knee Block (Required for Anterior Tilt)			
250DI37	E0957	Knee Blocks: Dual, Height, Width and Angle Adjustable	\$406.00
250DI55	E1028	Knee Block Quick Release Removable Hardware	\$138.00
SEDEO® PRO™ Swing-Away Legrest (Pair)			\$
250H177	—	Legrests, Swing-Away, 70°, SEDEO® PRO™ (14" to 21" Leg Length)	NC
250H178	E0990	Legrests, Manual Elevating, Swing-Away, SEDEO® PRO™ (Includes Padded Calf Pads) (14".5 to 21.5" Leg Length)	\$400.00
250H179	E1010	Legrests, Power Elevating, Swing-Away, SEDEO® PRO™ (Includes Padded Calf Pads) (15.5" to 22.5" Leg Length)	\$2,459.00
250DI76	—	Omit Legrests	NC
Swing-Away Legrest Footplates			
250F369	K0040	Angle Adjustable Flip-Up Footplates	\$193.00

POSITIONING SUPPORTS & ACCESSORIES

OPTION	HCPC	DESCRIPTION	RETAIL
Swing-Away Legrest Calf Supports			
250DI42	—	Calf Support, Depth Adjustable, Padded <i>(6.5" to 8" Wide, Depending On Chair Width, x 9.5" Tall)</i>	\$451.00
Swing-Away Legrest Options			
250DI44	E0951	Heel Loops	\$57.00
250KG4	—	Knee Pads, Integral Foam <i>(2.5" Wide x 4.5" Tall)</i>	\$229.00
SEDEO® Lower Extremity Swing-Away Support			
250VD01	E1020	Lower Extremity Support, Padded, Swing-Away, SEDEO®, Right	\$451.00
250VD02	E1020	Lower Extremity Support, Padded, Swing-Away, SEDEO®, Left	\$451.00
Pelvic Positioning Belts			
250A837	—	Positioning Belt, Auto Style Buckle, 50 mm 1600 mm / 2"	NC
250A839	—	Positioning Belt, Aircraft Buckle, 50 mm x 1300 mm / 2"	NC
250A840	E0978	Positioning Belt, Aircraft Buckle, Padded 50 mm x 1300 mm / 2"	\$86.00
250A845	K0108	Bodypoint - Evoflex Pelvic Stabilizer Belt	\$199.00
250A846	K0108	Bodypoint - Evoflex Pelvic Stabilizer Belt with Extension	\$220.00
Chest Straps			
250ZR06	—	Chest Strap, 2-Point, Padded <i>(Required)</i>	NC
250ZR07	—	Additional Chest Strap, 2-Point, Padded	\$286.00
Chest Bars			
250ZR08	—	Chest Bar, Detachable, Round Pad	\$567.00
250ZR09	—	Chest Bar, Detachable, Flat Pad, 200 mm x 100 mm	\$603.00
Accessories			\$
Ventilator Option			
250VT12	E1029	Vent Tray Fixed <i>(Pulmonetic LTV Models)</i>	\$531.00
O2 Tank Option			
250A718	E2208	O2 Tank Holder Mounted to Backrest <i>(D and E Styles)</i>	\$372.00
Wheelchair Tray Table			
2507A09	E0950	Tray, Swing-Away	\$400.00
Backpack			
250A1	—	Black Utility	\$62.00
Medial Thigh Support (N/A with Centermount Legrests)			
2508A80	E0957	Medial Thigh Support	\$193.00
2508A81	E1028	Pommel Swing-Away Hardware	\$286.00
Misc. Accessories			
2508A59	—	Transfer Bar <i>(Pair)</i>	\$400.00
2508A52	—	Ram Mount Phone	\$98.00
2508A53	—	Ram Mount Tablet	\$400.00
2508A54	—	Ram Mount Cup Holder	\$98.00
2508A58	—	Rear View Mirror	\$98.00
250A611	—	Crutch Holder <i>(N/A w/Power Recline)</i>	\$144.00
2508A73	—	Bag for Armrest	\$98.00
2508A75	—	Reflectors <i>(Backrest shroud mounted)</i>	\$98.00
2509A10	—	Utility Hooks	\$59.00

JOYSTICK / ELECTRONICS			
OPTION	HCPC	DESCRIPTION	RETAIL
R-net Joystick <i>(Remote)</i>			\$
250J121	—	Standard LED Joystick <i>(w/o Screen Display)</i>	NC
250J116	—	Color Joystick <i>(2.3" Color LCD Screen Display, CJS1M1)</i>	NC
250J118	—	High Resolution Color Joystick,(CJS1M2+BT) <i>3.5" Color LCD Display, Paddle Switches, Programmable Shortcut Keys, ICONs, Stereo Jack, IR and dual (IOS/Android) Bluetooth control</i>	\$1,206.00
250J119	—	Omit Remote Control <i>(Required with arm mounted specialty controls)</i>	NC
Joystick Mounting			\$
250JM34	—	Fixed on Armrest	NC
250JM35	E1028	Parallel Swing-Away	\$338.00
250JM36	E1028	Link-It Mid Line Mount on Armrest <i>(Swing-Away)</i>	\$972.00
Joystick Mount (Location)			\$
250HM16	—	Right Hand Armrest	NC
250HM17	—	Left Hand Armrest	NC
Display Language Selection R-net			\$
250L1	—	English	NC
250L3	—	French	NC
250L7	—	Spanish	NC
Joystick Handle (Optional)			\$
250JH13	E2323	Ball Handle	\$144.00
250JH15	E2323	T-Handle	\$144.00
250JH24	E2323	Bodypoint – Rubber Dome Joystick Handle	\$187.00
250JH25	E2323	Bodypoint – U-Shaped Joystick Handle <i>(4 inch)</i>	\$187.00
250JH26	E2323	Bodypoint – Grooved Mushroom Joystick	\$187.00
Electronic Programmers			\$
250PR10	—	R-net PC Programmer w/ Dongle	\$601.00
R-net Power Modules			\$
250ER11	—	Power Module, R-net, 120A, Non-Expandable <i>(Only available with Group 4)</i>	NC
250ER10	—	Power Module, R-net, 120A, Expandable <i>(Only available with Group 4)</i> <i>Required for specialty controls or 3+ actuators</i>	\$1,425.00
250EX09	E2377	Expandable 120A Controller - R-net	\$749.00
250EX10	E2313	Harness for 120A Expandable Controller	\$676.00
250ER18	—	Power Module, R-net, 90A, Non-Expandable	NC
250ER17	—	Power Module, R-net, 90A, Expandable <i>Required for specialty controls or 3+ actuators</i>	\$1,425.00
250EX15	E2377	Expandable 90A Controller - R-net	\$749.00
250EX16	E2313	Harness for 90A Expandable Controller	\$676.00

JOYSTICK / ELECTRONICS

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

Seating Control \$

Powered Seating Control

250E040	—	Omit Thru-Drive <i>(Requires Switch Control)</i>	NC
250E037	E2310	Thru-Drive Control for 1 Actuator	\$2,196.00
250E038	E2311	Thru-Drive Control for 2+ Actuators	\$3,975.00
250E032	E2311	CxSM, Thru-Drive Control for 2+ Actuators <i>(Required for PRO ADVANCED)</i>	\$3,975.00

Switch Control *(Available with Omit Thru-Drive only) (1 actuator for Tilt or Elevate only)*

250E033	—	Dual Toggle for 1 Actuator	NC
250E034	—	2.5" Single Button Switch for 1 Actuator <i>(Button switches are not mounted)</i>	NC
250E041	—	Ctrl+5 Switch Box for 2 Actuators or More	NC

Additional Power Seating Control *(Available with Thru-Drive only)*

250QG11	—	Ctrl+5 Switch Box with Free Assignable Functions	\$1,430.00
250QG12	—	Ctrl+5 Switch Box with Free Assignable Functions <i>(Multi-power only)</i>	NC

Dual Toggle Mounting Location

250TN1	—	Toggle Mount, Right	NC
250TN2	—	Toggle Mount, Left	NC

Ctrl+5 Mounting Location

250EY19	—	Right Armrest	NC
250EY20	—	Left Armrest	NC
250EY21	—	Back Cane, Right	\$172.00
250EY24	—	Back Cane, Left	\$172.00

Ctrl+5 Switch Type

250EY22	—	Push Buttons	NC
250EY23	K0108	Toggle Switches	\$343.00

Attendant Controls \$

2500A05	E2331	Dual Control Module for Attendant	\$1,487.00
2500A06	E2331	R-net Color LCD (CJSM1), Attendant Control Only	\$1,487.00

Attendant Mounting Location \$

250DN5	—	Right Back Cane	NC
250DN6	—	Left Back Cane	NC

R-net Programming \$

Assignable Buttons

Assignable Button feature allows for personalization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.

250PM12	—	Include Assignable Buttons	NC
---------	---	----------------------------	----

Power Centermount Footplate Extension Option

250IE03	K0108	Include 7.5" Powered Footplate Extension for Transfer Assist <i>(Electronics programmed to allow footplates to be extended with or without legrest being raised.) Only available with Powered Articulation.</i>	\$1,430.00
---------	-------	---	------------

JOYSTICK / ELECTRONICS

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

R-net Omni2 Options *For Expandable Controls, choose all that apply.* **\$**

2500N14	—	R-net Omni2 Only Input Control Module and 3.5" LCD Display w/ Built-In Infrared & Bluetooth Choose Mount Below	\$2,031.00
2500N15	—	R-net Omni2 Bundled Input Control Module and 3.5" LCD Display w/ Built-In Infrared & Bluetooth Bundled Omni2 Available w/ Switch-It Input Devices	NC

Omni2 Mount *If used w/ Standard Joystick, it must be mounted opposite side.* **\$**

250DS14	E1028	Link-It with Parallel Swing-Away Omni2 Mount w/ Rear Mounted Interface, Right	\$567.00
250DS15	E1028	Link-It with Parallel Swing-Away Omni2 Mount w/ Rear Mounted Interface, Left	\$567.00

PRO ADVANCED Programming **\$**
Biometric Repositioning

250PM03	K0108	Include Biometric Repositioning <i>(Requires Memory Position #1, #2, #3 below)</i>	NC
---------	-------	--	----

Memory Seating

"Memory Seating" feature allows for up to six custom seating positions that are accessed either through the joystick, Ctrl+5, joystick shortcut keys or programmable jacks on Omni2 and joystick. Choose desired preset settings below. Default positions are dealer adjustable.

250PM06	—	Include #1 Memory Seating (Seated) <i>(Required)</i>	NC
250PM07	—	Include #2 Memory Seating (Supine)	NC
250PM08	—	Include #3 Memory Seating (Pressure Relief)	NC
250PM09	—	Include #4 Memory Seating (Feet Above Heart)	NC
250PM10	—	Include #5 Memory Seating (Transfer)	NC
250PM11	—	Include #6 Memory Seating (C-Me)	NC

PRAR Pressure Relief Alert and Reminders

PRAR Pressure Relief Alert and Record feature provides programmable alerts to be set into the R-net electronics that will prompt the user when it is time to perform a pressure relief. The R-net system will record how often the pressure relief is completed or ignored.

250PM14	K0108	Include PRAR <i>(Minimum requirement Memory Position #1 and #3)</i>	NC
---------	-------	---	----

SPECIALTY INPUT DEVICES

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

Hand Control Packages *Requires Standard Joystick Omit / Only one hand control device can be ordered.*

\$

Proportional Input Devices

250EU14	E2312	Switch-It MicroGuide	\$3,769.00
250EU15	E2312	Switch-It MicroPilot	\$3,769.00
250EU17	E2328	Switch-It VersaGuide	\$2,076.00
250EU18	E2328	Switch-It VersaGuide EZ	\$2,453.00

Mounting Position *(Required)*

250EM04	E1028	Link-It to the Backrest <i>(Not Available w/ Omit Back Cushion Type)</i>	\$1,087.00
250EM05	E1028	Link-It to the Fixed Joystick Mount	\$338.00
250EM06	E1028	Link-It to the Parallel Swing-Away Joystick Mount	\$567.00
250EM07	E1028	Link-It to the Midline Swing-Away	\$1,201.00
250EM08	E1028	Omit Mount	NC

Mounting Platforms *(Optional)* *Not available w/ Link-It to Backrest*

250MG16	E0950	SunDisk Hand Tray <i>(Only available w/ MicroSeries Joysticks)</i>	\$338.00
250MG15	E0950	Bullet Hand Tray	\$338.00

Mounting Location *(Required)*

250IM09	—	Right	NC
250IM10	—	Left	NC

Switch Options

—	—	See Switch Selections <i>(Maximum 2)</i>	—
---	---	--	---

Head Control Packages *(Select headrest mount)*

\$

Only one head control device can be ordered.

Proportional Input Devices

250EU25	E2328	Switch-It DualPro® Proportional Head Array - Facial Pads	\$6,537.00
250EU26	E2328	Switch-It DualPro® Proportional Head Array - Spot Pads	\$6,537.00

Switched Input Devices

250EU28	E2330	Switch-It 3-Switch Head Array - Short Curved Wings	\$5,026.00
250EU30	E2330	Switch-It 3-Switch Head Array - Long Curved Wings	\$5,026.00

Switch Options

—	—	See Switch Selections <i>(Maximum 2)</i>	—
---	---	--	---

Sip 'N Puff

\$

250ID73	E2330	Switch-It Sip 'N Puff Combo Head Array with Link-It Hardware <i>(Requires headrest mount)</i>	\$5,713.00
250ID74	K0108	Switch-It OptiStop with Sip 'N Puff	\$2,225.00

Headrest Mounts *(Applies to Head/Sip 'N Puff control packages)*

\$

250IJ13	E1028	Whitmyer® Linx Headrest Mount with Quick Release	\$188.00
250IJ14	E1028	Whitmyer® Linx-2 Headrest Mount with Quick Release	\$229.00
250IJ15	E1028	Whitmyer® AXYS Headrest Mount with Quick Release	\$121.00

SPECIALTY INPUT DEVICES

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

Chin Control Packages

\$

Proportional Input Devices

250EU34	E2312	Switch-It MicroGuide <i>(Switches mounted to the headrest only.)</i>	\$3,769.00
250EU35	E2312	Switch-It MicroPilot <i>(Switches mounted to the headrest only.)</i>	\$3,769.00
250EU36	E2328	Switch-It VersaGuide	\$2,076.00
250EU37	E2328	Switch-It VersaGuide EZ	\$2,453.00

Mounting Position

250EM09	E1028	Link-It To Chin Harness, Small (150mm / 6")	\$1,430.00
250EM10	E1028	Link-It To Chin Harness, Large (220mm / 8")	\$1,430.00
250EM11	E1028	Link-It To Headrest (Swing-Away)	\$1,430.00
250EM53	E1028	Link-It To Manual Chin Boom - Right Side	\$4,341.00
250EM54	E1028	Link-It To Manual Chin Boom - Left Side	\$4,341.00
250EM16	E1028	Omit Mount	NC

Input Device Mounting Location

250IM13	—	Right	NC
250IM14	—	Left	NC

Switch Options

—	—	See Switch Selections	—
---	---	-----------------------	---

Safety Devices

\$

250SA06	—	Switch-It Remote Stop Switch - XLR Connector	\$1,070.00
---------	---	--	------------

Switched Driving Packages

\$

250SW44	E2328	Switch-It™ Cool Cube with DB9 connector for proximity/proportional switch driving. <i>(Pro Spot proportional switches are required if proportional driving is desired.)</i> <i>(For non-proportional driving, single switches can be used.)</i>	\$3,969.00
---------	-------	---	------------

SPECIALTY INPUT DEVICES

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

Switches

Single switches can be programmed to activate any function of the wheelchair including: On/Off, Mode, Profile, Speeds, Seating, Memory Positions, Horn, Lights & Indicators, IR functions, Bluetooth mouse clicks, Environmental Controls when plugged into assignable port.

Switches can connect via: Joysticks (Standard, Micro series, Versa series) Omni2, SJX6 Hub, head controls

Switch ordering instructions:

1. Select switch type and quantity in Section "A" using columns below. (Max 6 switches. See limitations below.)
2. Select switch mounting position from Section "B" using corresponding switch column.
3. Select switch mounting location from Section "C" using corresponding switch column.

Section A – Switch Selections

Switch / Quantity

Single Function Switches - Mechanical Switches	Switch 1	Switch 2	Switch 3	Switch 4	Switch 5	Switch 6	\$
2.5" Single Switch - Black							\$223.00
2.5" Single Switch - Blue							\$223.00
2.5" Single Switch - Green							\$223.00
2.5" Single Switch - Red							\$223.00
2.5" Single Switch - Yellow							\$223.00
Egg Switch - Black							\$126.00
Egg Switch - Blue							\$126.00
Egg Switch - Green							\$126.00
Egg Switch - Red							\$126.00
Egg Switch - Yellow							\$126.00
MicroLite <i>(Includes color labels)</i>							\$188.00
Mini Button 1 - Blue <i>(Velcro)</i>							\$121.00
Mini Button 1 - Green <i>(Velcro)</i>							\$121.00
Mini Button 1 - Red <i>(Velcro)</i>							\$121.00
1" Single Switch - Black							\$315.00
1" Single Switch - Blue							\$315.00
1" Single Switch - Green							\$315.00
1" Single Switch - Orange							\$315.00
1" Single Switch - Red							\$315.00
1" Single Switch - Yellow							\$315.00

Electronic Switches

Flex Switch (Headrest mount only)							\$245.00
-----------------------------------	--	--	--	--	--	--	----------

ProSpot - Proportional Switches

ProSpot Proportional Switch - Red							\$543.00
ProSpot Proportional Switch - Green							\$543.00
ProSpot Proportional Switch - Blue							\$543.00
ProSpot Proportional Switch - Yellow							\$543.00
ProSpot Proportional Switch - Orange							\$543.00

Section B – Switch Mounting Position

Link-It to the Backrest <i>(Up to 1 per side)</i>							\$744.00
Link-It to the Headrest - Gooseneck & Rigid <i>(Up to 2 per side)</i>							\$744.00
Link-It to the Fixed Joystick Mount <i>(Up to 2 per side)</i>							\$567.00
Link-It to the Parallel Swing-Away Joystick Mount <i>(Up to 2 per side)</i>							\$567.00
Link-It to the Seat Rail <i>(Up to 2 per side)</i> <i>(Not avail w/knee & thigh supports)</i>							\$451.00
Omit Mount <i>(includes Velcro)</i>							NC

Section C – Switch Mounting Location

Right side of chair							NC
Left side of chair							NC

Joystick Integrated Switches (Pair)

250J003	–	Proximity Switches on Flex Mount, Pair <i>(Used with VersaGuide/VersaGuide EZ)</i>	\$1138.00
---------	---	--	-----------

JAY® SYNCRA™ SEATING MEASURING GUIDE

JAY® SYNCRA™ BACK



Backrest Width	16"	18"	20"	22"
Max User Trunk Width	13.5"	15.5"	17.5"	19.5"


JAY® Swing-Away Lateral Thoracic Adjustment Range

Syncra Back Width	Width Range	Height range from bottom of shell to seat pan
16"	10" to 24"	12" to 19.5"
18"	12" to 26"	12.5" to 21.5"
20"	14" to 28"	12.5" to 22"
22"	16" to 30"	13.5" to 23"

Syncra™ Fixed Lateral Support Accessory Range

Syncra Back Width	Width Range	Height range from bottom of shell to seat pan	Contour Depth Range
16"	16.4" to 16.8"	6" to 8"	3" to 5"
18"	18.4" to 18.8"	6" to 8"	3" to 5"
20"	20.4" to 20.8"	6" to 8"	3" to 5"
22"	22.4" to 22.8"	6" to 8"	3" to 5"

Syncra™ Back Shell Height Diagram



JAY® SYNCRA™ CUSHION



JAY® Lateral Thigh Support Adjustment Range

Syncra Cushion Width	Width Range
16"	8.5" to 16.0"
17"	9.5" to 17.0"
18"	10.5" to 18.0"
19"	11.5" to 19.0"
20"	12.5" to 20.0"
22"	14.5" to 22.0"

built4me

Q700 M with SEDEO® PRO™

\$

Base Price

Total of additional equipment

Packaging and transport

Total Price

Date

Signature

Month

Day

Year

