


# Flip for Leckey 1.25" Soft Roll (Poly) Casters



*These casters effective February 8th, 2021 If replacing discontinued casters, a quantity of 2 will be needed. New casters do not match the sizes of the discontinued casters. If changing size of casters, new forks may be needed.*

Item Number	Part Number	Description	Quantity
1a	100786	Caster 6 x 1.25" Soft Roll <i>Replaced part number 100521</i>	1
1b	111863	Caster 8 x 1.25" Soft Roll <i>Replaced part number 100525</i>	1