


## Arc Dual Hand Tilt Mechanism - Growth



Cable (Item 1)		
Backrest Style	Extra-Short (12-15D)	Short, Medium, and Long
Ht. Adj. 18-23" Back Ht.	112940	112942
Ht. Adj. 24-27" Back Ht.	112941	112943
22" Fxd. Ht. Ang. Adj.	112940	112942
26" Fxd. Ht. Ang. Adj.	112941	112943
20" Stroller Handle	112938	112940
24" Stroller Handle	112939	112941
22" Recline	112938	112940
26" Recline	112939	112941

*Is your depth changing? Will your existing tilt cables achieve your new depth? If no, please select the appropriate tilt cable based off your new configuration on the table above. Each cable kit contains two cables and two e-rings.*

Item Number	Part Number	Description	Quantity
1a, 2	112938	Cable, Tilt, Dual, Short, Arc	1
1b, 2	112939	Cable, Tilt, Dual, Medium, Arc	1
1c, 2	112940	Cable, Tilt, Dual, Long, Arc	1
1d, 2	112941	Cable, Tilt, Dual, XL, Arc	1
1e, 2	112942	Cable, Tilt, Dual, XXL, Arc	1
1f, 2	112943	Cable, Tilt, Dual, XXXL, Arc	1