


# Little Wave XP Seat Upholstery - Growth



Seat upholstery width measures approximately 1" less than the chair size Seat upholstery depth measures true to the chair size If long sling present, depth of upholstery is chair depth + adjustable depth

Item Number	Part Number	Description	Quantity
1a	101500	Seat Upholstery Rigid 1010	1
1b	101501	Seat Upholstery Rigid 1011	1
1c	101502	Seat Upholstery Rigid 1012	1
1d	101503	Seat Upholstery Rigid 1013	1
1e	101504	Seat Upholstery Rigid 1014	1
1f	101505	Seat Upholstery Rigid 1015	1
1g	101506	Seat Upholstery Rigid 1016	1
1h	101565	Seat Upholstery Rigid 1017	1
1i	101566	Seat Upholstery Rigid 1018	1
1j	101507	Seat Upholstery Rigid 1110	1

Item Number	Part Number	Description	Quantity
1k	101508	Seat Upholstery Rigid 1111	1
1l	101509	Seat Upholstery Rigid 1112	1
1m	101510	Seat Upholstery Rigid 1113	1
1n	101511	Seat Upholstery Rigid 1114	1
1o	101512	Seat Upholstery Rigid 1115	1
1p	101513	Seat Upholstery Rigid 1116	1
1q	101567	Seat Upholstery Rigid 1117	1
1r	101568	Seat Upholstery Rigid 1118	1
1s	101170	Seat Upholstery Rigid 1210	1
1t	101228	Seat Upholstery Rigid 1211	1
1u	101514	Seat Upholstery Rigid 1212	1
1v	101515	Seat Upholstery Rigid 1213	1
1w	101516	Seat Upholstery Rigid 1214	1
1x	101517	Seat Upholstery Rigid 1215	1
1y	101518	Seat Upholstery Rigid 1216	1
1z	101519	Seat Upholstery Rigid 1217	1
1aa	101520	Seat Upholstery Rigid 1218	1
1ab	101569	Seat Upholstery Rigid 1219	1
1ac	101552	Seat Upholstery Rigid 1310	1
1ad	101553	Seat Upholstery Rigid 1311	1
1ae	101521	Seat Upholstery Rigid 1312	1
1af	101522	Seat Upholstery Rigid 1313	1
1ag	101523	Seat Upholstery Rigid 1314	1
1ah	101524	Seat Upholstery Rigid 1315	1
1ai	101525	Seat Upholstery Rigid 1316	1
1aj	101526	Seat Upholstery Rigid 1317	1
1ak	101527	Seat Upholstery Rigid 1318	1
1al	101571	Seat Upholstery Rigid 1319	1

Item Number	Part Number	Description	Quantity
2a	000548	Seat Rail 10" Deep RIGID	2
2b	000549	Seat Rail 11" Deep RIGID	2
2c	000550	Seat Rail 12" Deep RIGID	2
2d	000551	Seat Rail 13" Deep RIGID	2
2e	000241	Seat Rail 14" Deep Rigid/Liberty	2
2f	000242	Seat Rail 15" Deep RIGID	2
2g	000243	Seat Rail 16" Deep Rigid/Liberty	2
2h	000244	Seat Rail 17" Deep RIGID	2
2i	000245	Seat Rail 18" Deep Rigid/Liberty	2
2j	000246	Seat Rail 19" Deep RIGID	2
3	100734	M5 x 16 PHMS ZC	12
4a	101570	Seat Upholstery Rigid 1410	1
4b	101077	Seat Upholstery Rigid 1411	1
4c	101078	Seat Upholstery Rigid 1412	1
4d	101079	Seat Upholstery Rigid 1413	1
4e	101000	Seat Upholstery Rigid/Liberty 1414	1
4f	101005	Seat Upholstery Rigid 1415	1
4g	101012	Seat Upholstery Rigid/Liberty 1416	1
4h	101019	Seat Upholstery Rigid 1417	1
4i	101026	Seat Upholstery Rigid/Liberty 1418	1
4j	101033	Seat Upholstery Rigid 1419	1
4k	101573	Seat Upholstery Rigid 1510	1
4l	101574	Seat Upholstery Rigid 1511	1
4m	101575	Seat Upholstery Rigid 1512	1
4n	101080	Seat Upholstery Rigid 1513	1
4o	101001	Seat Upholstery Rigid/Liberty 1514	1
4p	101006	Seat Upholstery Rigid 1515	1
4q	101013	Seat Upholstery Rigid/Liberty 1516	1

Item Number	Part Number	Description	Quantity
4r	101020	Seat Upholstery Rigid 1517	1
4s	101027	Seat Upholstery Rigid/Liberty 1518	1
4t	101034	Seat Upholstery Rigid 1519	1
4u	101576	Seat Upholstery Rigid 1610	1
4v	101577	Seat Upholstery Rigid 1611	1
4w	101578	Seat Upholstery Rigid 1612	1
4x	101579	Seat Upholstery Rigid 1613	1
4y	101002	Seat Upholstery Rigid/Liberty 1614	1
4z	101007	Seat Upholstery Rigid 1615	1
4aa	101014	Seat Upholstery Rigid/Liberty 1616	1
4ab	101021	Seat Upholstery Rigid 1617	1
4ac	101028	Seat Upholstery Rigid/Liberty 1618	1
4ad	101035	Seat Upholstery Rigid 1619	1
4ae	102454	Seat Upholstery Rigid 1712	1
4af	101580	Seat Upholstery Rigid 1713	1
4ag	101003	Seat Upholstery Rigid/Liberty 1714	1
4ah	101008	Seat Upholstery Rigid 1715	1
4ai	101015	Seat Upholstery Rigid/Liberty 1716	1
4aj	101022	Seat Upholstery Rigid 1717	1
4ak	101029	Seat Upholstery Rigid/Liberty 1718	1
4al	101036	Seat Upholstery Rigid 1719	1
4am	101004	Seat Upholstery Rigid/Liberty 1814	1
4an	101009	Seat Upholstery Rigid 1815	1
4ao	101016	Seat Upholstery Rigid/Liberty 1816	1
4ap	101023	Seat Upholstery Rigid 1817	1
4aq	101030	Seat Upholstery Rigid/Liberty 1818	1
4ar	101037	Seat Upholstery Rigid 1819	1
4as	101010	Seat Upholstery Rigid 1915	1

<b>Item Number</b>	<b>Part Number</b>	<b>Description</b>	<b>Quantity</b>
4at	101017	Seat Upholstery Rigid/Liberty 1916	1
4au	101024	Seat Upholstery Rigid 1917	1
4av	101031	Seat Upholstery Rigid/Liberty 1918	1
4aw	101038	Seat Upholstery Rigid 1919	1