


Rogue XP Seat Upholstery



SCREW (ITEM 3) QTY PER CHAIR	
SEAT DEPTH	SCREW QTY
10" - 12"	6
13" - 15"	8
16" - 18"	10
19" - 20"	12

Seat upholstery width measures approximately 1" less than the chair size. Seat upholstery depth measures true to the chair size.

Item Number	Part Number	Description	Quantity	Retail
1a	101500	Seat Upholstery Rigid 1010	1	\$94.55
1b	101501	Seat Upholstery Rigid 1011	1	\$94.55

Item Number	Part Number	Description	Quantity	Retail
1c	101502	Seat Upholstery Rigid 1012	1	\$94.55
1d	101503	Seat Upholstery Rigid 1013	1	\$94.55
1e	101504	Seat Upholstery Rigid 1014	1	\$94.55
1f	101505	Seat Upholstery Rigid 1015	1	\$94.55
1g	101506	Seat Upholstery Rigid 1016	1	\$94.55
1h	101565	Seat Upholstery Rigid 1017	1	\$94.55
1i	101566	Seat Upholstery Rigid 1018	1	\$94.55
1j	102528	Seat Upholstery Rigid 1019	1	\$94.55
1k	101507	Seat Upholstery Rigid 1110	1	\$94.55
1l	101508	Seat Upholstery Rigid 1111	1	\$94.55
1m	101509	Seat Upholstery Rigid 1112	1	\$94.55
1n	101510	Seat Upholstery Rigid 1113	1	\$94.55
1o	101511	Seat Upholstery Rigid 1114	1	\$94.55
1p	101512	Seat Upholstery Rigid 1115	1	\$94.55
1q	101513	Seat Upholstery Rigid 1116	1	\$94.55
1r	101567	Seat Upholstery Rigid 1117	1	\$94.55
1s	101568	Seat Upholstery Rigid 1118	1	\$94.55
1t	102529	Seat Upholstery Rigid 1119	1	\$94.55
1u	101170	Seat Upholstery Rigid 1210	1	\$94.55
1v	101228	Seat Upholstery Rigid 1211	1	\$94.55
1w	101514	Seat Upholstery Rigid 1212	1	\$94.55
1x	101515	Seat Upholstery Rigid 1213	1	\$94.55
1y	101516	Seat Upholstery Rigid 1214	1	\$94.55
1z	101517	Seat Upholstery Rigid 1215	1	\$94.55
1aa	101518	Seat Upholstery Rigid 1216	1	\$94.55
1ab	101519	Seat Upholstery Rigid 1217	1	\$94.55
1ac	101520	Seat Upholstery Rigid 1218	1	\$94.55
1ad	101569	Seat Upholstery Rigid 1219	1	\$94.55

Item Number	Part Number	Description	Quantity	Retail
1ae	101552	Seat Upholstery Rigid 1310	1	\$94.55
1af	101553	Seat Upholstery Rigid 1311	1	\$94.55
1ag	101521	Seat Upholstery Rigid 1312	1	\$94.55
1ah	101522	Seat Upholstery Rigid 1313	1	\$94.55
1ai	101523	Seat Upholstery Rigid 1314	1	\$94.55
1aj	101524	Seat Upholstery Rigid 1315	1	\$94.55
1ak	101525	Seat Upholstery Rigid 1316	1	\$94.55
1al	101526	Seat Upholstery Rigid 1317	1	\$94.55
1am	101527	Seat Upholstery Rigid 1318	1	\$94.55
1an	101571	Seat Upholstery Rigid 1319	1	\$94.55
1ao	101572	Seat Upholstery Rigid 1320	1	\$94.55
1ap	101570	Seat Upholstery Rigid 1410	1	\$94.55
1aq	101077	Seat Upholstery Rigid 1411	1	\$94.55
1ar	101078	Seat Upholstery Rigid 1412	1	\$94.55
1as	101079	Seat Upholstery Rigid 1413	1	\$94.55
1at	101000	Seat Upholstery Rigid/LFT 1414	1	\$94.55
1au	101005	Seat Upholstery Rigid 1415	1	\$94.55
1av	101012	Seat Upholstery Rigid/LFT 1416	1	\$94.55
1aw	101019	Seat Upholstery Rigid 1417	1	\$94.55
1ax	101026	Seat Upholstery Rigid/LFT 1418	1	\$94.55
1ay	101033	Seat Upholstery Rigid 1419	1	\$94.55
1az	101074	Seat Upholstery Rigid/LFT 1420	1	\$94.55
1ba	101573	Seat Upholstery Rigid 1510	1	\$94.55
1bb	101574	Seat Upholstery Rigid 1511	1	\$94.55
1bc	101575	Seat Upholstery Rigid 1512	1	\$94.55
1bd	101080	Seat Upholstery Rigid 1513	1	\$94.55
1be	101001	Seat Upholstery Rigid/LFT 1514	1	\$94.55
1bf	101006	Seat Upholstery Rigid 1515	1	\$94.55

Item Number	Part Number	Description	Quantity	Retail
1bg	101013	Seat Upholstery Rigid/LFT 1516	1	\$94.55
1bh	101020	Seat Upholstery Rigid 1517	1	\$94.55
1bi	101027	Seat Upholstery Rigid/LFT 1518	1	\$94.55
1bj	101034	Seat Upholstery Rigid 1519	1	\$94.55
1bk	101040	Seat Upholstery Rigid/LFT 1520	1	\$94.55
1bl	101576	Seat Upholstery Rigid 1610	1	\$94.55
1bm	101577	Seat Upholstery Rigid 1611	1	\$94.55
1bn	101578	Seat Upholstery Rigid 1612	1	\$94.55
1bo	101579	Seat Upholstery Rigid 1613	1	\$94.55
1bp	101002	Seat Upholstery Rigid/LFT 1614	1	\$94.55
1bq	101007	Seat Upholstery Rigid 1615	1	\$94.55
1br	101014	Seat Upholstery Rigid/LFT 1616	1	\$94.55
1bs	101021	Seat Upholstery Rigid 1617	1	\$94.55
1bt	101028	Seat Upholstery Rigid/LFT 1618	1	\$94.55
1bu	101035	Seat Upholstery Rigid 1619	1	\$94.55
1bv	101041	Seat Upholstery Rigid/LFT 1620	1	\$94.55
1bw	102454	Seat Upholstery Rigid 1712	1	\$94.55
1bx	101580	Seat Upholstery Rigid 1713	1	\$94.55
1by	101003	Seat Upholstery Rigid/LFT 1714	1	\$94.55
1bz	101008	Seat Upholstery Rigid 1715	1	\$94.55
1ca	101015	Seat Upholstery Rigid/LFT 1716	1	\$94.55
1cb	101022	Seat Upholstery Rigid 1717	1	\$94.55
1cc	101029	Seat Upholstery Rigid/LFT 1718	1	\$94.55
1cd	101036	Seat Upholstery Rigid 1719	1	\$94.55
1ce	101042	Seat Upholstery Rigid/LFT 1720	1	\$94.55
1cf	101004	Seat Upholstery Rigid/LFT 1814	1	\$94.55
1cg	101009	Seat Upholstery Rigid 1815	1	\$94.55
1ch	101016	Seat Upholstery Rigid/LFT 1816	1	\$94.55

Item Number	Part Number	Description	Quantity	Retail
1ci	101023	Seat Upholstery Rigid 1817	1	\$94.55
1cj	101030	Seat Upholstery Rigid/LFT 1818	1	\$94.55
1ck	101037	Seat Upholstery Rigid 1819	1	\$94.55
1cl	101043	Seat Upholstery Rigid/LFT 1820	1	\$94.55
2a	000548	Seat Rail 10" Deep RIGID	2	\$12.32
2b	000549	Seat Rail 11" Deep RIGID	2	\$12.32
2c	000550	Seat Rail 12" Deep RIGID	2	\$12.32
2d	000551	Seat Rail 13" Deep RIGID	2	\$12.64
2e	000241	Seat Rail 14" Deep Rigid/LFT	2	\$15.68
2f	000242	Seat Rail 15" Deep RIGID	2	\$15.68
2g	000243	Seat Rail 16" Deep Rigid/LFT	2	\$15.68
2h	000244	Seat Rail 17" Deep RIGID	2	\$15.68
2i	000245	Seat Rail 18" Deep Rigid/LFT	2	\$15.68
2j	000246	Seat Rail 19" Deep RIGID	2	\$15.68
2k	000247	Seat Rail 20" Deep Rigid/LFT	2	\$15.68
3	100734	M5 x 16 PHMS ZC	14	\$1.58