



*Lightweight Design...Battery-Powered Option...Custom Designed*

## Intelect Transport

2782K – Intelect Transport Ultrasound

## Standard Accessories

27383 – 5 cm<sup>2</sup> Sound Head Applicator

4248 – Conductor™ Ultrasound Gel -  
8.5 oz (250 ml) Bottle (1x)

## Optional Accessories

27478 – Intelect Transport Nickel  
Metal Hydride Battery Pack

27381 – 1 cm<sup>2</sup> Sound Head Applicator

27382 – 2 cm<sup>2</sup> Sound Head Applicator

27384 – 10 cm<sup>2</sup> Sound Head Applicator

27467 – Intelect Transport Carry Bag

2911K – Battery Pack and Carry Bag

2780ASY – Therapy Cart

2891K – Therapy Cart with  
Cart Adapter (Gray)

2884 – Cart Adapter (Gray)  
Allows the Intelect Transport  
unit to attach securely to the  
Therapy Cart.

# Intelect® TRANSPORT ULTRASOUND



## Technical Specifications

Mains Power: 100-240 VAC-50/60 Hz

Weight: 5 lb. (2.3 kg)

Shipping Weight: 16 lb. (7.3 kg)

Dimensions: 12.9"L x 11.3"W x 6.4"H  
(32.8 cm x 28.8 cm x 16.3 cm)

Electrical Safety Class: Class 1, Type B

Safety Tests: UL/IEC/EN 60601-1,  
IEC/EN 60601-1-2, IEC 60601-2-5

- Fully functional 1 and 3.3 MHz frequencies
- Ten user-defined memory positions for User Protocols
- Pulsed and continuous therapy operation (10%, 20%, 50% and 100%)
- Ergonomically designed sound head applicators in variety of sizes (1 cm<sup>2</sup>, 2 cm<sup>2</sup>, 5 cm<sup>2</sup> and 10 cm<sup>2</sup>)
- Soundheads are watertight for use in underwater therapy
- All sound head applicators are interchangeable due to Electronic Signature™ residing on all transducers
- Head warming feature on all sound head applicators
- Unit uniquely designed for tabletop, wall mount, therapy cart or mobile use
- 110V or battery powered option
- Optional Therapy Cart integrates with the unit for maximizing storage of clinical supplies
- Two year warranty



Optional Battery Pack



Applicators



Cart Adapter