


Discontinued Spark Height Adjustable T-Arm



Short assemblies are used with 14-17 Seat Depths. Long assemblies are used with 18"-19" Seat Depths.

Item Number	Part Number	Description	Quantity
1a, 3-5, 6a, 7-15, 20	100238	Height Adj T Armrest Short Desk Right <i>Replaced with 003463</i>	1
1a, 3-5, 6a, 7-15, 20	100242	Height Adj T Armrest Short Desk Left <i>Replaced with 003463</i>	1
1b, 3-5, 6a, 7-15, 20	100239	Height Adj T Armrest Short Full Right <i>Replaced with 003464</i>	1
1b, 3-5, 6a, 7-15, 20	100243	Height Adj T Armrest Short Full Left <i>Replaced with 003464</i>	1
1a, 3-5, 6b, 7-15, 20	100240	Height Adj T Armrest Long Desk Right <i>Replaced with 003465</i>	1
1a, 3-5, 6b, 7-15, 20	100244	Height Adj T Armrest Long Desk Left <i>Replaced with 003465</i>	1
1b, 3-5, 6b, 7-15, 20	100241	Height Adj T Armrest Long Full Right <i>Replaced with 003466</i>	1
1b, 3-5, 6b, 7-15, 20	100245	Height Adj T Armrest Long Full Left <i>Replaced with 003466</i>	1
1a, 2b, 3-5, 6a, 7-15, 19a, 20	501682	Height Adj T-Arm Assm Short Desk w/ Std Pad Left <i>Replaced with 505646</i>	1
1a, 2b, 3-5, 6a, 7-15, 19a, 20	501683	Height Adj T-Arm Assm Short Desk w/ Std Pad Right <i>Replaced with 505646</i>	1
1a, 2b, 3-5, 6a, 7-15, 19a, 20	503051	Height Adj T Armrest, 14"-17" Depths, Desk w/ Std Pad PAIR CAT <i>Replaced with 505647</i>	1
1b, 2a, 3-5, 6a, 7-15, 19a, 20	501684	Height Adj T-Arm Assm Short Full w/ Std Pad Left <i>Replaced with 505644</i>	1
1b, 2a, 3-5, 6a, 7-15, 19a, 20	501685	Height Adj T-Arm Assm Short Full w/ Std Pad Right <i>Replaced with 505644</i>	1
1b, 2a, 3-5, 6a, 7-15, 19a, 20	503053	Height Adj T Armrest, 14"-17" Depths, Full w/ Std Pad PAIR CAT <i>Replaced with 505645</i>	1
1a, 2b, 3-5, 6b, 7-15, 19a, 20	501686	Height Adj T-Arm Assm Long Desk w/ Std Pad Left <i>Replaced with 505493</i>	1
1a, 2b, 3-5, 6b, 7-15, 19a, 20	501687	Height Adj T-Arm Assm Long Desk w/ Std Pad Right <i>Replaced with 505493</i>	1
1a, 2b, 3-5, 6b, 7-15, 19a, 20	503052	Height Adj T Armrest, 18"-20" Depths, Desk w/ Std Pad PAIR CAT <i>Replaced with 505494</i>	1

Item Number	Part Number	Description	Quantity
1b, 2a, 3-5, 6b, 7-15, 19a, 20	501688	Height Adj T-Arm Assm Long Full w/ Std Pad Left <i>Replaced with 505491</i>	1
1b, 2a, 3-5, 6b, 7-15, 19a, 20	501689	Height Adj T-Arm Assm Long Full w/ Std Pad Right <i>Replaced with 505491</i>	1
1b, 2a, 3-5, 6b, 7-15, 19a, 20	503054	Height Adj T Armrest, 18"-20" Depths, Full w/ Std Pad PAIR CAT <i>Replaced with 505492</i>	1
1a, 2d, 3-5, 6a, 7-15, 19b, 20	501690	Height Adj T-Arm Assm Short Desk w/ Waterfall Pad Left <i>Replaced with 505650</i>	1
1a, 2d, 3-5, 6a, 7-15, 19b, 20	501691	Height Adj T-Arm Assm Short Desk w/ Waterfall Pad Right <i>Replaced with 505650</i>	1
1a, 2d, 3-5, 6a, 7-15, 19b, 20	503059	Height Adj T Armrest, 14"-17" Depths, Desk w/ Waterfall Pad PAIR CAT <i>Replaced with 505651</i>	1
1b, 2c, 3-5, 6a, 7-15, 19b, 20	501692	Height Adj T-Arm Assm Short Full w/ Waterfall Pad Left <i>Replaced with 505648</i>	1
1b, 2c, 3-5, 6a, 7-15, 19b, 20	501693	Height Adj T-Arm Assm Short Full w/ Waterfall Pad Right <i>Replaced with 505648</i>	1
1b, 2c, 3-5, 6a, 7-15, 19b, 20	503061	Height Adj T Armrest, 14"-17" Depths, Full w/ Waterfall Pad PAIR CAT <i>Replaced with 505649</i>	1
1a, 2d, 3-5, 6b, 7-15, 19b, 20	501694	Height Adj T-Arm Assm Long Desk w/ Waterfall Pad Left <i>Replaced with 505497</i>	1
1a, 2d, 3-5, 6b, 7-15, 19b, 20	501695	Height Adj T-Arm Assm Long Desk w/ Waterfall Pad Right <i>Replaced with 505497</i>	1
1a, 2d, 3-5, 6b, 7-15, 19b, 20	503060	Height Adj T Armrest, 18"-20" Depths, Desk w/ Waterfall Pad PAIR CAT <i>Replaced with 505498</i>	1
1b, 2c, 3-5, 6b, 7-15, 19b, 20	501696	Height Adj T-Arm Assm Long Full w/ Waterfall Pad Left <i>Replaced with 505495</i>	1
1b, 2c, 3-5, 6b, 7-15, 19b, 20	501697	Height Adj T-Arm Assm Long Full w/ Waterfall Pad Right <i>Replaced with 505495</i>	1
1b, 2c, 3-5, 6b, 7-15, 19b, 20	503062	Height Adj T Armrest, 18"-20" Depths, Full w/ Waterfall Pad PAIR CAT <i>Replaced with 505496</i>	1
1a, 2b, 3-5, 6a, 7-18, 19a, 20-22	505277	Height Adj T-Arm Assembly Short Desk w/ Std Pad and 1" Receiver Left Spark <i>Replaced with 506059</i>	1
1a, 2b, 3-5, 6a, 7-18, 19a, 20-22	505278	Height Adj T-Arm Assembly Short Desk w/ Std Pad and 1" Receiver Right Spark <i>Replaced with 506059</i>	1

Item Number	Part Number	Description	Quantity
1a, 2b, 3-5, 6a, 7-18, 19a, 20-22	505279	Height Adjustable T-Arm w/ Brackets and Hardware, 14"-17" Depths, Desk w/ Std Pad PAIR Spark <i>Replaced with 506060</i>	1
1b, 2a, 3-5, 6a, 7-18, 19a, 20-22	505280	Height Adj T-Arm Assembly Short Full w/ Std Pad and 1" Receiver Left Spark <i>Replaced with 506057</i>	1
1b, 2a, 3-5, 6a, 7-18, 19a, 20-22	505281	Height Adj T-Arm Assembly Short Full w/ Std Pad and 1" Receiver Right Spark <i>Replaced with 506057</i>	1
1b, 2a, 3-5, 6a, 7-18, 19a, 20-22	505282	Height Adj T-Arm w/ Brackets and Hardware, 14"-17" Depths, Full w/ Std Pad PAIR Spark <i>Replaced with 506058</i>	1
1a, 2b, 3-5, 6b, 7-18, 19a, 20-22	505283	Height Adj T-Arm Assembly Long Desk w/ Std Pad and 1" Receiver Left Spark <i>Replaced with 505511</i>	1
1a, 2b, 3-5, 6b, 7-18, 19a, 20-22	505284	Height Adj T-Arm Assembly Long Desk w/ Std Pad and 1" Receiver Right Spark <i>Replaced with 505511</i>	1
1a, 2b, 3-5, 6b, 7-18, 19a, 20-22	505285	Height Adj T-Arm w/ Brackets and Hardware, 18"-19" Depths, Desk w/ Std Pad PAIR Spark <i>Replaced with 505512</i>	1
1b, 2a, 3-5, 6b, 7-18, 19a, 20-22	505286	Height Adj T-Arm Assembly Long Full w/ Std Pad and 1" Receiver Left Spark <i>Replaced with 505509</i>	1
1b, 2a, 3-5, 6b, 7-18, 19a, 20-22	505287	Height Adj T-Arm Assembly Long Full w/ Std Pad and 1" Receiver Right Spark <i>Replaced with 505509</i>	1
1b, 2a, 3-5, 6b, 7-18, 19a, 20-22	505288	Height Adj T-Arm w/ Brackets and Hardware, 18"-19" Depths, Full w/ Std Pad PAIR Spark <i>Replaced with 505510</i>	1
1a, 2d, 3-5, 6a, 7-18, 19b, 20-22	505289	Height Adj T-Arm Assembly Short Desk w/ Waterfall Pad and 1" Receiver Left Spark <i>Replaced with 506063</i>	1
1a, 2d, 3-5, 6a, 7-18, 19b, 20-22	505290	Height Adj T-Arm Assembly Short Desk w/ Waterfall Pad and 1" Receiver Right Spark <i>Replaced with 506063</i>	1
1a, 2d, 3-5, 6a, 7-18, 19b, 20-22	505291	Height Adj T-Arm w/ Brackets and Hardware, 14"-17" Depths, Desk w/ Waterfall Pad PAIR Spark <i>Replaced with 506064</i>	1
1b, 2c, 3-5, 6a, 7-18, 19b, 20-22	505292	Height Adj T-Arm Assembly Short Full w/ Waterfall Pad and 1" Receiver Left Spark <i>Replaced with 506061</i>	1
1b, 2c, 3-5, 6a, 7-18, 19b, 20-22	505293	Height Adj T-Arm Assembly Short Full w/ Waterfall Pad and 1" Receiver Right Spark	1

Item Number	Part Number	Description	Quantity
		<i>Replaced with 506061</i>	
1b, 2c, 3-5, 6a, 7-18, 19b, 20-22	505294	Height Adj T-Arm w/ Brackets and Hardware, 14"-17" Depths, Full w/ Waterfall Pad PAIR Spark <i>Replaced with 506062</i>	1
1a, 2d, 3-5, 6b, 7-18, 19b, 20-22	505295	Height Adj T-Arm Assembly Long Desk w/ Waterfall Pad and 1" Receiver Left Spark <i>Replaced with 505515</i>	1
1a, 2d, 3-5, 6b, 7-18, 19b, 20-22	505296	Height Adj T-Arm Assembly Long Desk w/ Waterfall Pad and 1" Receiver Right Spark <i>Replaced with 505515</i>	1
1a, 2d, 3-5, 6b, 7-18, 19b, 20-22	505297	Height Adj T-Arm w/ Brackets and Hardware, 18"-19" Depths, Desk w/ Waterfall Pad PAIR Spark <i>Replaced with 505516</i>	1
1b, 2c, 3-5, 6b, 7-18, 19b, 20-22	505298	Height Adj T-Arm Assembly Long Full w/ Waterfall Pad and 1" Receiver Left Spark <i>Replaced with 505513</i>	1
1b, 2c, 3-5, 6b, 7-18, 19b, 20-22	505299	Height Adj T-Arm Assembly Long Full w/ Waterfall Pad and 1" Receiver Right Spark <i>Replaced with 505513</i>	1
1b, 2c, 3-5, 6b, 7-18, 19b, 20-22	505300	Height Adj T-Arm w/ Brackets and Hardware, 18"-19" Depths, Full w/ Waterfall Pad PAIR Spark <i>Replaced with 505514</i>	1
16, 17, 18, 21, 22	505276	Height Adjustable T-Arm Receiver Assembly 1" Spark	1
1a, 20	000091	Height Adj T Armrest Upper Desk	1
1b, 20	000101	Height Adj T Armrest Upper Full	1
2a	100500	Arm Pad Full	1
2b	100501	Arm Pad Desk	1
3	000090	Height Adjustable T Armrest Shim	2
4	000089	Height Adjustable T Armrest Post	1
6a	000097	Height Adj T Armrest Side Guard Short <i>Used with 14" - 17" Seat Depths</i>	1
6b	000120	Height Adj T Armrest Side Guard Long <i>Used with 18" - 19" Seat Depths</i>	1
7	000095	Height Adjustable T Armrest Transfer Tube	1
9	200257	Dowel Pin 1/8 x 1/2	1
10	000134	Height Adjustable Armrest Lever Spring	1
11	000099	Height Adjustable Armrest Lever Pin	1

Item Number	Part Number	Description	Quantity
12	000096	Height Adjustable T Armrest Pin Lock Housing	1
8-12	500541	Arm Height Adjustable T Upper Lever Assembly	1
14	000135	Height Adjustable T Armrest Release Lever Spring	1
15	000152	Height Adjustable T Armrest Release Spring Button	1
16	002723	Height Adjustable T-Arm Receiver Spark	1
17	000088	Height Adjustable T Armrest Bracket	1
18	100866	M6 x 20 CLS12.9 DIN7991 FSHCS Blk Zc w/Patch	2
19a	100696	M5 x 30mm PPH Capscrew <i>Used with Standard Arm Pads</i>	2
19b	100654	M5 x 40 DIN7985 PH Pan HD M/S ZC <i>Used with Waterfall Arm Pads</i>	2
20	100803	Tube End Plug 7/8"	2
21	100667	M6 x 12 BTN SHCS BLKZC w/Patch	2
22	000065	M6 x 23.68 Threaded Barrel	2
23a	102452	Waterfall Arm Pad Full	1
23b	102521	Waterfall Arm Pad Desk	1