

Bill To (Altimate Medical Supplier Only)	Ship To
Purchase Order #: _____ Quote #: _____	Name: _____
Account #: _____ Ordered By: _____	Address: _____
Supplier Company: _____	Address: _____
Address: _____ City, State, Zip: _____	City, State: _____
Phone: _____ Fax: _____	Zip: _____
Email: _____ Payer Source: _____	Phone: _____
Name of ATP (Required): _____	Tag For: _____

Bantam Extra Small User Specifications	
Approximate Height Range	28"-40"
Weight Limit	50 lbs.
Seat Depth Range (seat to back surface)	7"-12"
Seat to Foot Plate Range (from seat pivot)	4.5"-16.5"

Base-must select one

The Base includes: Two 5" Front Wheels, Height-Adjustable Foot Plates, Planar Seat and Back, Chest Strap, Black Upholstery, Two 3" Rear Locking Swivel Casters.

Default items include: Frame Color, Knee Pads, Knee Pad Size, Foot Plates, Foot Plate Size and Back Support.

Note: The default items will be added if no selection is made.

<input type="checkbox"/> PT50001-1 EasyStand® Bantam Extra Small Supine Ready Base \$2,173 -Supine CANNOT be ordered off-unit	
<input type="checkbox"/> PT50234 Supine (head support required, shadow tray is recommended)*	\$450
Lift Mechanisms Accessories	
<input type="checkbox"/> PA5614 Gas Spring Lift Lockout	\$108
<input type="checkbox"/> PT50236 Dual Control (both foot and hand controls, push handle PT50280 required)	\$108

<input type="checkbox"/> PT50001-1 EasyStand® Bantam Extra Small Base \$2,173 -Supine not available with Pow'r Up Lift or Hydraulic Actuator	
Lift Mechanisms-must select one	
<input type="checkbox"/> PT50022 Pow'r Up Lift (battery powered, not available with hydraulic actuator)	\$1,222
<input type="checkbox"/> PT50020 Hydraulic Actuator with Removable Handle (not available with pow'r up lift)	\$730
Lift Mechanisms Accessories	
<input type="checkbox"/> PNG50044 Quad Grip Handle Extension (hydraulic actuator required)	\$153
<input type="checkbox"/> PNG50070 T-Style Handle Extension (hydraulic actuator required)	\$153

Frame Color-must select one

- Green default N/C
 - Purple N/C
 - White N/C
 - Orange N/C
- Go to easystand.com to see specific frame colors*

Seat Adjustment-CANNOT be ordered off-unit

- PNG50499 Easy-Adjust Seat Depth (roller bearing allows for smoother seat depth adjustment) \$108

Front Frame*

- PT50116 Front Frame (5" Front Swivel Casters, replaces front wheels) \$179

Mobile

- PT50177 Mobile (not available with supine or front swivel casters) \$1,373

Mobile Tray

- PT50246 Mobile Tray \$479
(depth and height adjustable, not available w/form to fit backs, lateral supports or elbow stop)

Shadow Trays*

- PT50204 Black Molded Tray-11"Lx19"W (depth & height adjustable) \$473
- PA5600 Black Molded Swing-Away Tray-11"Lx19"W (depth, height & angle adjustable) \$709
- PA5602 Clear Swing-Away Tray-11"Lx19"W (depth, height & angle adjustable) \$844
- PT50074 Black Molded Tray-21"Lx24"W (depth, height & angle adjustable) \$592
- PT50076 Clear Tray-21"Lx24"W (depth, height & angle adjustable) \$861

Tray Accessories-tray required*

- PNG50042 Hand Grips (for clear tray only, pair) \$153
- PT50156 Elbow Pad-9.25" Cutout (21"Lx24"W tray only) \$144
- PT50160 Elbow Pad-11" Cutout (21"Lx24"W tray only) \$144
- PA5604 Forearm Wings-Black Molded Tray (11"Lx19"W tray only, pair) \$144
- PA5606 Forearm Wings-Clear Tray (11"Lx19"W tray only, pair) \$144
- PA3038 Hygienic Padded Cover (11"Lx19"W tray only, must select pattern) \$199
- PT30093 Hygienic Padded Cover (21"Lx24"W tray only, must select pattern) \$199

* not available with mobile

Knee Pads-must select one

<input type="checkbox"/> PT50262-1 Slide-In Knee Pads (pair) default	\$328
<input type="checkbox"/> PT50264-1 Swing-Away Knee Pads (pair)	\$543

Knee Pad Size-must select one, inside width

<input type="checkbox"/> PT50266 2.5"	N/C
<input type="checkbox"/> PT50268 3.25" default	N/C
<input type="checkbox"/> PT50270 4.25"	N/C

Foot Plates

<input type="checkbox"/> PT50252 Multi-Adjustable Foot Plates (replaces height-adjustable foot plates, pair)	\$216
--	-------

Foot Holder Size-must select one

<input type="checkbox"/> PT50254 6"Lx3"W	N/C
<input type="checkbox"/> PT50256 7.75"Lx3.25"W default	N/C
<input type="checkbox"/> PT50258 9.75"Lx4"W	N/C

Foot Plate Accessories

<input type="checkbox"/> PT50080 Foot Straps-10"L (length over top of foot, pair)	\$46
<input type="checkbox"/> PY5636 Foot Straps-15"L (length over top of foot, pair)	\$46
<input type="checkbox"/> PT50082 Secure Foot Straps-10"L (length over top of foot, two pair)	\$88
<input type="checkbox"/> PY5638 Secure Foot Straps-15"L (length over top of foot, two pair)	\$88
<input type="checkbox"/> PT50278 Calf Pad	\$142

Seats

<input type="checkbox"/> PT50044 Form to Fit Seat (replaces planar seat)	\$354
--	-------

Seat Accessories

<input type="checkbox"/> PT30082 Hygienic Cover for Planar Seat (must select pattern)	\$80
<input type="checkbox"/> PT30083 Hygienic Cover for Form to Fit Seat (must select pattern)	\$80
<input type="checkbox"/> PT50098 Hip Supports (6"-9" range, pair, not available with form to fit seat)	\$363
<input type="checkbox"/> PA3039 Hygienic Cover for Hip Supports (must select pattern, pair)	\$90

Backs-height from seat, must select one

<input type="checkbox"/> PT50050 Planar Back-11"-13"H Range default	N/C
<input type="checkbox"/> PT50238 Planar Back-15"-17"H Range	N/C
<input type="checkbox"/> PT50052 Form to Fit Back-11"-13"H Range (replaces planar back)	\$354
<input type="checkbox"/> PT50240 Form to Fit Back-15"-17"H Range (replaces planar back)	\$354

Back Accessories

<input type="checkbox"/> PT30078 Hygienic Cover for Planar Back-11"-13"H Range (must select pattern)	\$80
<input type="checkbox"/> PT30089 Hygienic Cover for Planar Back-15"-17"H Range (must select pattern)	\$80
<input type="checkbox"/> PT30079 Hygienic Cover for PT50052 (must select pattern)	\$80
<input type="checkbox"/> PT30090 Hygienic Cover for PT50240 (must select pattern)	\$80

Notes:

Positioning and Upper Body Support

<input type="checkbox"/> PT50008 Lateral Supports-8"-12" Range Curved Pad (for planar back only, pair)*	\$363
<input type="checkbox"/> PT50274 Lateral Supports-8"-12" Range Flat Pad (for planar back only, pair)*	\$363
<input type="checkbox"/> PT50248 Lateral Supports-5"-9" Range Curved Pad (for planar back only, pair)*	\$363
<input type="checkbox"/> PT50276 Lateral Supports-5"-9" Range Flat Pad (for planar back only, pair)*	\$363
<input type="checkbox"/> PB3049 Hygienic Cover for Lateral Supports (must select pattern, pair)*	\$102
<input type="checkbox"/> PT50216 Lateral Supports with Elbow Stop-8"-12" Range (for planar back only, pair)*	\$371
<input type="checkbox"/> PT50218 Lateral Supports with Elbow Stop and Arm Rest-8"-12" Range* (for planar back only, pair)	\$410
<input type="checkbox"/> PT50220 Elbow Stop with Arm Rest-8"-12" Range (for planar and form to fit back, pair)*	\$371
<input type="checkbox"/> PT30066 Positioning Belt (fits hip circumference up to 23", velcro with D-ring)	\$57
<input type="checkbox"/> PT30094 Positioning Belt (fits hip circumference up to 33", velcro with D-ring)	\$57
<input type="checkbox"/> PT50086 X-Style Chest Vest-9.5"Lx9"W	\$277
<input type="checkbox"/> PT50088 X-Style Chest Vest-11"Lx9.5"W	\$277
<input type="checkbox"/> PT50302 Non-Stretch X-Style Chest Vest-9.5"Lx9"W	\$323
<input type="checkbox"/> PNG30318 Non-Stretch X-Style Chest Vest-11"Lx9.5"W	\$323
<input type="checkbox"/> PT50128 High Mount Chest Vest Bracket	\$78

Head Support-supine requires a head support

<input type="checkbox"/> PT50062 Head Support-5"Hx8"W	\$363
<input type="checkbox"/> PT50064 Head Support-6"Hx10"W	\$363
<input type="checkbox"/> PT50284 Head Support-Form to Fit-5"Hx10"W	\$363

Head Support Accessories

<input type="checkbox"/> PA3035 Hygienic Cover for Head Support-5"Hx8"W (must select pattern)	\$57
<input type="checkbox"/> PA3036 Hygienic Cover for Head Support-6"Hx10"W (must select pattern)	\$57
<input type="checkbox"/> PA3037 Hygienic Cover for Head Support-Form to Fit (must select pattern)	\$57

Other Accessories

<input type="checkbox"/> PT50280 Push Handle (required for dual control)	\$153
<input type="checkbox"/> P80809 Tool Pouch	\$25
<input type="checkbox"/> P80246 Angle Locator	\$36

Hygienic Cover Patterns (must select one if hygienic covers were ordered)

<input type="checkbox"/> Rainbow Unicorn	<i>Go to easystand.com to see specific fabric patterns</i>
<input type="checkbox"/> Sea Creatures	
<input type="checkbox"/> Fancy Pants Elephants	
<input type="checkbox"/> Space Race	
<input type="checkbox"/> All Smiles Animals	

* not available with mobile

AltimateMedical

Ph: 800.342.8968 or 507.697.6393
 Fax: 877.342.8968 or 507.697.6900

info@easystand.com easystand.com

Email form to webquotes@easystand.com