

## Catalyst 5 Wheels - Superlight Spoke

---



Not available with Heavy Duty option. <p>

If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles may be needed. <p>

Please note tab mounted handrims must be installed prior to the tire being installed.

<b>Item Number</b>	<b>Part Number</b>	<b>Description</b>	<b>Quantity</b>	<b>Retail</b>
1a	101159	24" Superlight Spoke Wheel	1	\$191.46
1b	101158	25" Superlight Spoke Wheel	1	\$191.46
2a	102827	Replacement Spoke Superlight 24"	1	\$6.70
2b	102828	Replacement Spoke Superlight 25"	1	\$6.97