



Account No. _____ Drop Ship: _____
 Date: _____ Name _____
 P.O. Number: _____ Address _____
 Buyer: _____ City _____
 Marked For: _____ State _____ Zip _____
 E Mail: _____ Tel. _____ Fax _____




A (Axilla to Floor Height) _____
 W (Wrist to Floor) _____
 G (Groin to Floor) _____
 H (Hip Width) _____



Features included in standard price







- Pony**
 Rear drag wheel with directional stabilizer
 Leg & Foot Separator
 Height adjustable seat with spring
 Height adjustable removable back support
 Height adjustable chest support
 5 wheels (1 rear & 2 front with locks)
 Height and depth adjustable handle bar







- Bronco**
 4 Locking casters with directional stabilizers
 Gas spring assisted torso angle adjustment
 Leg & Foot Separator
 Height adjustable seat
 Height adjustable removable back support
 Height adjustable chest support
 Height and depth adjustable handle bar

Pony	Description	HCPCS: E8002	Size	Part No.	Price
	Pony, Yellow		size 0	86400	2,325 <input type="checkbox"/>
	Pony, Purple		size 1	86401	2,325 <input type="checkbox"/>
	Pony, Lilac		size 2	86402	2,325 <input type="checkbox"/>
Bronco					
	Bronco, Purple		size 1	86405-F	2,750 <input type="checkbox"/>
	Bronco, Lilac		size 2	86406-F	2,750 <input type="checkbox"/>





Accessories

	Hip & thigh supports, set	6" w x 3" h	size 0	86408	330	<input type="checkbox"/>
	Hip & thigh supports, set	6" w x 4" h	size 1	86410	330	<input type="checkbox"/>
	Hip & thigh supports, set	6" w x 4" h	size 2	86412	330	<input type="checkbox"/>
	Hip & thigh supports, set	6" w x 4" h	size 3	86413	330	<input type="checkbox"/>
	Minimum/Maximum Width	Size 0 6" x 8 3/4"	Size 1 6" x 8 3/4"	Size 2 6 1/2" x 10 3/4"	Size 3 10" x 14 1/2"	
	Abdominal cushion		size 1	86420	195	<input type="checkbox"/>
	Abdominal cushion		size 2	86422	195	<input type="checkbox"/>
	Abdominal cushion		size 3	86423	195	<input type="checkbox"/>
	Padded Armrest		size 0	86427	485	<input type="checkbox"/>
	Padded Armrest		size 1	86428	485	<input type="checkbox"/>
	Padded Armrest		size 2	86429	485	<input type="checkbox"/>
	Padded Armrest		size 3	86433	485	<input type="checkbox"/>
	Vertical hand grips - variable mount		all sizes	36-0432	71	<input type="checkbox"/>


Select a Head support, Swing-away Bracket, and Headrest Bracket from below

	Pad Style					
	Head support no. 1, Black Polyurethane (cover optional)	6171550		205	<input type="checkbox"/>	
	Head support no. 2, Black Polyurethane (cover optional)	6172550		205	<input type="checkbox"/>	
	Head support no. 3, Black Polyurethane (cover optional)	6173550		205	<input type="checkbox"/>	
	Head support no. 4, Black Polyurethane (cover optional)	6174550		205	<input type="checkbox"/>	
	Anatomic head support*, Small - A = 3", B = 8"	5099001		255	<input type="checkbox"/>	
	Anatomic head support*, Medium - A = 3", B = 10 1/4"	5099002		255	<input type="checkbox"/>	
Anatomic head support*, Large - A = 4 1/4", B = 13"	5099003		255	<input type="checkbox"/>		
	<i>*Includes black Ducutex cover</i>					
	Headband System (For Anatomic Head Support only)	99500		230	<input type="checkbox"/>	
	Extra headband strap (For Anatomic Head Support only)	99400-5		56	<input type="checkbox"/>	
	Swing-away Mounting Bracket - Height, Depth, Angle & Lateral Adjustable					
	Turtle bar system, short	98619-25		230	<input type="checkbox"/>	
	Turtle bar system, long	98619-26		230	<input type="checkbox"/>	
	Pony/Bronco Headrest bracket (required)	all sizes	861424-3	258	<input type="checkbox"/>	

Accessories and Spare parts for Pony and Bronco

	Knee pad, Pony	one size	86450	220	<input type="checkbox"/>
	Knee pad, Bronco	small	86461	170	<input type="checkbox"/>
	Knee pad, Bronco	large	86462	170	<input type="checkbox"/>
	Push handle, Pony	one size	86416	255	<input type="checkbox"/>
	Push handle, Bronco	one size	86417	235	<input type="checkbox"/>
	Directional stabilizer for front wheel	size 0-2	86448	370	<input type="checkbox"/>
	Directional stabilizer for front wheel	size 3	86449	370	<input type="checkbox"/>

Measurements



Pony	Size 0 inch	Size 1 inch	Size 2 inch
Length (L)	26¼"	28"	40½"
Width (W)	23"	27"	30"
Height (H)	20¼-27"	23½-32"	29¼-40½"
Seat height (SH)	9¼-14¾"	14¾-19½"	18¾-25"
Angle range (A)	32°	32°	32°
Weight	15.5 lb	22 lb	30 lb
Maximum user weight	66 lb	88 lb	110 lb
Approximate user age	0 to 3	3 to 6	6 to 9
Bronco	Size 1 inch	Size 2 inch	
Length (L)	41¾"	48¾"	
Width (W)	25¾"	28½"	
Trunk pad height (H)	22½-31½"	27¼-40¼"	
Seat height (SH)	12¾-20"	16¾-23¾"	
Tilt angle (A)	17°	17°	
Weight	30 lb	34 lb	
Maximum user weight	88 lb	110 lb	
Approximate user age	3 to 6	6 to 9	