


Flip Attendant Foot Lock



This design effective starting with serial number FLP007536. If retrofitting/replacing, for a full assembly please select a wheel size, axle tube and the Attendant Foot Lock Kit.

| Item Number | Part Number | Description | Quantity |
|--------------------|-------------|--------------------------------------|----------|
| 3-5, 6a, 7a, 13-17 | 505470 | Attendant Foot Lock Kit 10 Wide FLIP | 1 |
| 3-5, 6b, 7b, 13-17 | 505471 | Attendant Foot Lock Kit 11 Wide FLIP | 1 |
| 3-5, 6c, 7c, 13-17 | 505472 | Attendant Foot Lock Kit 12 Wide FLIP | 1 |
| 3-5, 6d, 7d, 13-17 | 505473 | Attendant Foot Lock Kit 13 Wide FLIP | 1 |
| 3-5, 6e, 7e, 13-17 | 505474 | Attendant Foot Lock Kit 14 Wide FLIP | 1 |
| 3-5, 6f, 7f, 13-17 | 505475 | Attendant Foot Lock Kit 15 Wide FLIP | 1 |
| 3-5, 6g, 7g, 13-17 | 505476 | Attendant Foot Lock Kit 16 Wide FLIP | 1 |

| Item Number | Part Number | Description | Quantity |
|--------------------------------------|-------------|---|----------|
| 3-5, 6h, 7h, 13-17 | 505477 | Attendant Foot Lock Kit 17 Wide FLIP <i>Used with axle receiver 100756 (Discontinued 3/3/2023)</i> | 1 |
| 3-5, 6i, 7i, 13-17 | 505478 | Attendant Foot Lock Kit 18 Wide FLIP <i>Used with axle receiver 114480</i> | 1 |
| 4, 6a, 7a, 8a, 13, 18a, 18d, 19a, 20 | 505479 | Strut Kit 10" Wide Flip w/ Attendant Foot Lock | 1 |
| 4, 6b, 7b, 8b, 13, 18b, 18e, 19b, 20 | 505480 | Strut Kit 11" Wide Flip w/ Attendant Foot Lock | 1 |
| 4, 6c, 7c, 8c, 13, 18c, 18f, 19c, 20 | 505481 | Strut Kit 12" Wide Flip w/ Attendant Foot Lock | 1 |
| 4, 6d, 7d, 8d, 13, 18d, 18g, 19d, 20 | 505482 | Strut Kit 13" Wide Flip w/ Attendant Foot Lock | 1 |
| 4, 6e, 7e, 8e, 13, 18e, 18h, 19e, 20 | 505483 | Strut Kit 14" Wide Flip w/ Attendant Foot Lock | 1 |
| 4, 6f, 7f, 8f, 13, 18f, 18i, 19f, 20 | 505484 | Strut Kit 15" Wide Flip w/ Attendant Foot Lock | 1 |
| 4, 6g, 7g, 8g, 13, 18g, 18j, 19g, 20 | 505485 | Strut Kit 16" Wide Flip w/ Attendant Foot Lock | 1 |
| 4, 6h, 7h, 8h, 13, 18h, 18k, 19h, 20 | 505486 | Strut Kit 17" Wide Flip w/ Attendant Foot Lock | 1 |
| 4, 6i, 7i, 8i, 13, 18i, 18l, 19i, 20 | 505487 | Strut Kit 18" Wide Flip w/ Attendant Foot Lock | 1 |
| 1a | 502853 | Flip Attendant Foot Lock Wheel, 12" Pneumatic EACH | 1 |
| 1b | 502854 | Flip Attendant Foot Lock Wheel, 12" Airless EACH | 1 |
| 1c | 502855 | Flip Attendant Foot Lock Wheel, 12" Full Poly EACH | 1 |
| 1d | 502856 | Flip Attendant Foot Lock Wheel, 16" Pneumatic EACH | 1 |
| 1e | 502857 | Flip Attendant Foot Lock Wheel, 16" Airless EACH | 1 |
| 1f | 502858 | Flip Attendant Foot Lock Wheel, 16" Full Poly EACH | 1 |
| 1g | 113604 | Flip Attendant Foot Lock Wheel, 16" Low Poly EACH | 1 |
| 2a | 101108 | Axle Quick Release w/ Button 1/2" x 99MM | 2 |
| 2b | 101191 | Axle Quick Release w/ Button 1/2" x 110MM | 2 |
| 3 | 001928 | Lever, Hub Lock | 2 |
| 5 | 001929 | Spring, Hub Lock | 2 |
| 4, 6a, 13 | 505461 | Locking Rod Hub Lock 10 Wide Kit | 1 |

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|---|----------|
| 4, 6b, 13 | 505462 | Locking Rod Hub Lock 11 Wide Kit | 1 |
| 4, 6c, 13 | 505463 | Locking Rod Hub Lock 12 Wide Kit | 1 |
| 4, 6d, 13 | 505464 | Locking Rod Hub Lock 13 Wide Kit | 1 |
| 4, 6e, 13 | 505465 | Locking Rod Hub Lock 14 Wide Kit | 1 |
| 4, 6f, 13 | 505466 | Locking Rod Hub Lock 15 Wide Kit | 1 |
| 4, 6g, 13 | 505467 | Locking Rod Hub Lock 16 Wide Kit | 1 |
| 4, 6h, 13 | 505468 | Locking Rod Hub Lock 17 Wide Kit | 1 |
| 4, 6i, 13 | 505469 | Locking Rod Hub Lock 18 Wide Kit | 1 |
| 7a | 001939 | Link Rod Hub Lock 10 Wide | 1 |
| 7b | 001940 | Link Rod Hub Lock 11 Wide | 1 |
| 7c | 001941 | Link Rod Hub Lock 12 Wide <i>Used with the reclining backrest.</i> | 1 |
| 7d | 001942 | Link Rod Hub Lock 13 Wide <i>Replaced part number 100756</i> | 1 |
| 7e | 001943 | Link Rod Hub Lock 14 Wide | 1 |
| 7f | 001944 | Link Rod Hub Lock 15 Wide | 1 |
| 7g | 001945 | Link Rod Hub Lock 16 Wide | 1 |
| 7h | 001946 | Link Rod Hub Lock 17 Wide | 1 |
| 7i | 001947 | Link Rod Hub Lock 18 Wide | 1 |
| 8a, 20 | 115883 | Axle Tube with Wear Strip, 10 Wide | 1 |
| 8b, 20 | 115884 | Axle Tube with Wear Strip, 11 Wide | 1 |
| 8c, 20 | 115885 | Axle Tube with Wear Strip, 12 Wide | 1 |
| 8d, 20 | 115886 | Axle Tube with Wear Strip, 13 Wide | 1 |
| 8e, 20 | 115887 | Axle Tube with Wear Strip, 14 Wide | 1 |
| 8f, 20 | 115888 | Axle Tube with Wear Strip, 15 Wide | 1 |
| 8g, 20 | 115889 | Axle Tube with Wear Strip, 16 Wide | 1 |
| 8h, 20 | 115890 | Axle Tube with Wear Strip, 17 Wide | 1 |
| 8i, 20 | 115891 | Axle Tube with Wear Strip, 18 Wide | 1 |
| 9a | 001838 | Axle Bar, Standard | 2 |

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|--|----------|
| 9b | 001839 | Axle Bar, Self Propel | 2 |
| 9c | 003838 | Axle Bar, Extended | 1 |
| 10 | 114480 | Axle Receiver, 2.5" Length | 2 |
| 11 | 100514 | Axle Receiver Nut M18 x D12 x 8mm | 2 |
| 12 | 100583 | Axle Plate Washer 18mm | 2 |
| 13 | 103503 | 3/8" External Self-Locking Shaft Ring | 2 |
| 14 | 100670 | M6 x 45 BTN SHCS BLK ZC | 4 |
| 15 | 100669 | M5 x 10 BTN SHCS BLKZC w/Patch | 2 |
| 16 | 100658 | M6 Nylock Nut Blk Zn | 4 |
| 17 | 100746 | M6 Flat Washer Blk Zn | 2 |
| 18a | 002719 | Frame Growth Strut 6.5" <i>Used on 10" wide chairs</i> | 1 |
| 18b | 002629 | Frame Growth Strut 7.5" <i>Used on 11" wide chairs</i> | 1 |
| 18c | 002630 | Frame Growth Strut 8.5" <i>Used on 12" wide chairs</i> | 1 |
| 18d | 002631 | Frame Growth Strut 9.5" <i>Used on 10" and 13" wide chairs</i> | 1 |
| 18e | 002632 | Frame Growth Strut 10.5" <i>Used on 11" and 14" wide chairs</i> | 1 |
| 18f | 002633 | Frame Growth Strut 11.5" <i>Used on 12" and 15" wide chairs</i> | 1 |
| 18g | 002634 | Frame Growth Strut 12.5" <i>Used on 13" and 16" wide chairs</i> | 1 |
| 18h | 002635 | Frame Growth Strut 13.5" <i>Used on 14" and 17" wide chairs</i> | 1 |
| 18i | 002636 | Frame Growth Strut 14.5" <i>Used on 15" and 18" wide chairs</i> | 1 |
| 18j | 002637 | Frame Growth Strut 15.5" <i>Used on 16" wide chairs</i> | 1 |
| 18k | 002638 | Frame Growth Strut 16.5" <i>Used on 17" wide chairs</i> | 1 |
| 18l | 002639 | Frame Growth Strut 17.5" | 1 |

| Item Number | Part Number | Description | Quantity |
|--------------------|--------------------|--------------------------------|-----------------|
| | | <i>Used on 18" wide chairs</i> | |
| 19a | 001879 | Strut Tube Seat 10 Wide | 1 |
| 19b | 001880 | Strut Tube Seat 11 Wide | 1 |
| 19c | 001881 | Strut Tube Seat 12 Wide | 1 |
| 19d | 001882 | Strut Tube Seat 13 Wide | 1 |
| 19e | 001883 | Strut Tube Seat 14 Wide | 1 |
| 19f | 001884 | Strut Tube Seat 15 Wide | 1 |
| 19g | 001885 | Strut Tube Seat 16 Wide | 1 |
| 19h | 001886 | Strut Tube Seat 17 Wide | 1 |
| 19i | 001887 | Strut Tube Seat 18 Wide | 1 |