




Focus CR Adjustable Height Rigidizer Bar - Growth



Backpost	
Part No.	Width
000185	14.3"
000186	15.3"
000187	16.3"
000188	17.3"
000189	18.3"
000190	19.3"
000191	20.3"
000529	21.3"
000530	22.3"

Has your width changed? If yes, please select accordingly.

Item Number	Part Number	Description	Quantity
1a	000185	Adjustable Height Rigidizer Bar 14" Wide	1
1b	000186	Adjustable Height Rigidizer Bar 15" Wide	1

Item Number	Part Number	Description	Quantity
1c	000187	Adjustable Height Rigidizer Bar 16" Wide	1
1d	000188	Adjustable Height Rigidizer Bar 17" Wide	1
1e	000189	Adjustable Height Rigidizer Bar 18" Wide	1
1f	000190	Adjustable Height Rigidizer Bar 19" Wide	1
1g	000191	Adjustable Height Rigidizer Bar 20" Wide	1
1h	000529	Adjustable Height Rigidizer Bar 21"	1
1i	000530	Adjustable Height Rigidizer Bar 22"	1
2	000987	Threaded Fastener	4
3	100669	M5 x 10 BTN SHCS BLKZC w/Patch	4
4	100797	M6 x 40 BHCS BLK ZC	4
5	100558	M6 Nyl Ins Domed Cap L/N Zc	4