


Catalyst Rear Anti-Tip



If changing to or from a reclining configuration, new sized anti-tips may be required. If retrofitting, Anti-Tip Receivers are needed.

Item Number	Part Number	Description	Quantity
1, 2	501509	Anti Tip Extra Long <i>Commonly used with Reclining Backrest Option</i>	1
1, 3	501508	Anti Tip Long <i>Commonly used with 25/26" Wheels on Catalyst. Used as needed with high seat heights.</i>	1
1, 4	501507	Anti Tip Medium <i>Commonly used with 22/24" Wheels on Catalyst</i>	1
1, 5	501506	Anti Tip Short <i>Commonly used with 18/20" Wheels on Catalyst</i>	1
1, 6	503131	Anti Tip Extra Short <i>Used as needed with low seat heights</i>	1

Item Number	Part Number	Description	Quantity
1	100533	Detent Button	1
2	200524	Anti Tip Extra Long 24-26 RG Cat w/Reclining Back	1
3	100228	Anti Tip Long 25/26CAT 22ROGUE 22-26RG ALX/LW XP 26ETHOS	1
4	100226	Anti Tip Medium 22/24CAT 18/20RIGID 22/24/25ETHOS	1
5	200522	Anti Tip Short 18/20CAT	1
6	102852	Anti Tip Extra Short Assembly, Super Low	1