


# Rogue ALX Fixed Height Backrest Frame with Non-Adjustable Rigidizer Bar (Formerly Tsunami)



Backrest 10"	10.13"
Backrest 11"	11.13"
Backrest 12"	12.13"
Backrest 13"	13.13"
Backrest 14"	14.13"
Backrest 15"	15.13"
Backrest 16"	16.13"
Backrest 17"	17.13"
Backrest 18"	18.13"
Backrest 19"	19.13"
Backrest 20"	20.13"
Backrest 21"	21.13"
Backrest 22"	22.13"

Item Number	Part Number	Description	Quantity
1a	108953	Fixed Height Rigid Backrest 7.5", 10" Wide <i>Used with 9" back height with omit push handles, 13" back height with fold down push handles or 10" back height with standard or low profile push handles</i>	1
1b	108952	Fixed Height Rigid Backrest 8.5", 10" Wide <i>Used with 10" back height with omit push handles, 14" back height with fold down push handles or 11" back height with standard or low profile push handles</i>	1
1c	108951	Fixed Height Rigid Backrest 9.5", 10" Wide <i>Used with 11" back height with omit push handles, 15" back height with fold down push handles or 12" back height with standard or low profile push handles</i>	1
1d	108950	Fixed Height Rigid Backrest 10.5", 10" Wide <i>Used with 12" back height with omit push handles, 16" back height with fold down push handles or 13" back height with standard or low profile push handles</i>	1
1e	108949	Fixed Height Rigid Backrest 11.5", 10" Wide <i>Used with 13" back height with omit push handles, 17" back height with fold down push handles or 14" back height with standard or low profile push handles</i>	1
1f	108948	Fixed Height Rigid Backrest 12.5", 10" Wide <i>Used with 14" back height with omit push handles, 18" back height with fold down push handles or 15" back height with standard or low profile push handles</i>	1
1g	108947	Fixed Height Rigid Backrest 13.5", 10" Wide <i>Used with 15" back height with omit push handles or 16" back height with standard or low profile push handles</i>	1
1h	108946	Fixed Height Rigid Backrest 14.5", 10" Wide <i>Used with 16" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1i	108945	Fixed Height Rigid Backrest 15.5", 10" Wide <i>Used with 17" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1j	108944	Fixed Height Rigid Backrest 16.5", 10" Wide <i>Used with 18" back height with omit push handles</i>	1
1k	108964	Fixed Height Rigid Backrest 7.5", 11" Wide <i>Used with 9" back height with omit push handles, 13" back height with fold down push handles or 10" back height with standard or low profile push handles</i>	1
1l	108963	Fixed Height Rigid Backrest 8.5", 11" Wide <i>Used with 10" back height with omit push handles, 14" back height with fold down push handles or 11" back height with standard or low profile push handles</i>	1
1m	108962	Fixed Height Rigid Backrest 9.5", 11" Wide <i>Used with 11" back height with omit push handles, 15" back height with fold down push handles or 12" back height with standard or low profile push handles</i>	1

Item Number	Part Number	Description	Quantity
1n	108961	Fixed Height Rigid Backrest 10.5", 11" Wide <i>Used with 12" back height with omit push handles, 16" back height with fold down push handles or 13" back height with standard or low profile push handles</i>	1
1o	108960	Fixed Height Rigid Backrest 11.5", 11" Wide <i>Used with 13" back height with omit push handles, 17" back height with fold down push handles or 14" back height with standard or low profile push handles</i>	1
1p	108959	Fixed Height Rigid Backrest 12.5", 11" Wide <i>Used with 14" back height with omit push handles, 18" back height with fold down push handles or 15" back height with standard or low profile push handles</i>	1
1q	108958	Fixed Height Rigid Backrest 13.5", 11" Wide <i>Used with 15" back height with omit push handles or 16" back height with standard or low profile push handles</i>	1
1r	108957	Fixed Height Rigid Backrest 14.5", 11" Wide <i>Used with 16" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1s	108956	Fixed Height Rigid Backrest 15.5", 11" Wide <i>Used with 17" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1t	108955	Fixed Height Rigid Backrest 16.5", 11" Wide <i>Used with 18" back height with omit push handles</i>	1
1u	108975	Fixed Height Rigid Backrest 7.5", 12" Wide <i>Used with 9" back height with omit push handles, 13" back height with fold down push handles or 10" back height with standard or low profile push handles</i>	1
1v	108974	Fixed Height Rigid Backrest 8.5", 12" Wide <i>Used with 10" back height with omit push handles, 14" back height with fold down push handles or 11" back height with standard or low profile push handles</i>	1
1w	108973	Fixed Height Rigid Backrest 9.5", 12" Wide <i>Used with 11" back height with omit push handles, 15" back height with fold down push handles or 12" back height with standard or low profile push handles</i>	1
1x	108972	Fixed Height Rigid Backrest 10.5", 12" Wide <i>Used with 12" back height with omit push handles, 16" back height with fold down push handles or 13" back height with standard or low profile push handles</i>	1
1y	108971	Fixed Height Rigid Backrest 11.5", 12" Wide <i>Used with 13" back height with omit push handles, 17" back height with fold down push handles or 14" back height with standard or low profile push handles</i>	1
1z	108970	Fixed Height Rigid Backrest 12.5", 12" Wide <i>Used with 14" back height with omit push handles, 18" back height with fold down push handles or 15" back height with standard or low profile push handles</i>	1

Item Number	Part Number	Description	Quantity
1aa	108969	Fixed Height Rigid Backrest 13.5", 12" Wide <i>Used with 15" back height with omit push handles or 16" back height with standard or low profile push handles</i>	1
1ab	108968	Fixed Height Rigid Backrest 14.5", 12" Wide <i>Used with 16" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1ac	108967	Fixed Height Rigid Backrest 15.5", 12" Wide <i>Used with 17" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1ad	108966	Fixed Height Rigid Backrest 16.5", 12" Wide <i>Used with 18" back height with omit push handles</i>	1
1ae	108986	Fixed Height Rigid Backrest 7.5", 13" Wide <i>Used with 9" back height with omit push handles, 13" back height with fold down push handles or 10" back height with standard or low profile push handles</i>	1
1af	108985	Fixed Height Rigid Backrest 8.5", 13" Wide <i>Used with 10" back height with omit push handles, 14" back height with fold down push handles or 11" back height with standard or low profile push handles</i>	1
1ag	108984	Fixed Height Rigid Backrest 9.5", 13" Wide <i>Used with 11" back height with omit push handles, 15" back height with fold down push handles or 12" back height with standard or low profile push handles</i>	1
1ah	108983	Fixed Height Rigid Backrest 10.5", 13" Wide <i>Used with 12" back height with omit push handles, 16" back height with fold down push handles or 13" back height with standard or low profile push handles</i>	1
1ai	108982	Fixed Height Rigid Backrest 11.5", 13" Wide <i>Used with 13" back height with omit push handles, 17" back height with fold down push handles or 14" back height with standard or low profile push handles</i>	1
1aj	108981	Fixed Height Rigid Backrest 12.5", 13" Wide <i>Used with 14" back height with omit push handles, 18" back height with fold down push handles or 15" back height with standard or low profile push handles</i>	1
1ak	108980	Fixed Height Rigid Backrest 13.5", 13" Wide <i>Used with 15" back height with omit push handles or 16" back height with standard or low profile push handles</i>	1
1al	108979	Fixed Height Rigid Backrest 14.5", 13" Wide <i>Used with 16" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1am	108978	Fixed Height Rigid Backrest 15.5", 13" Wide <i>Used with 17" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1an	108977	Fixed Height Rigid Backrest 16.5", 13" Wide <i>Used with 18" back height with omit push handles</i>	1

Item Number	Part Number	Description	Quantity
1ao	108997	Fixed Height Rigid Backrest 7.5", 14" Wide <i>Used with 9" back height with omit push handles, 13" back height with fold down push handles or 10" back height with standard or low profile push handles</i>	1
1ap	108996	Fixed Height Rigid Backrest 8.5", 14" Wide <i>Used with 10" back height with omit push handles, 14" back height with fold down push handles or 11" back height with standard or low profile push handles</i>	1
1aq	108995	Fixed Height Rigid Backrest 9.5", 14" Wide <i>Used with 11" back height with omit push handles, 15" back height with fold down push handles or 12" back height with standard or low profile push handles</i>	1
1ar	108994	Fixed Height Rigid Backrest 10.5", 14" Wide <i>Used with 12" back height with omit push handles, 16" back height with fold down push handles or 13" back height with standard or low profile push handles</i>	1
1as	108993	Fixed Height Rigid Backrest 11.5", 14" Wide <i>Used with 13" back height with omit push handles, 17" back height with fold down push handles or 14" back height with standard or low profile push handles</i>	1
1at	108992	Fixed Height Rigid Backrest 12.5", 14" Wide <i>Used with 14" back height with omit push handles, 18" back height with fold down push handles or 15" back height with standard or low profile push handles</i>	1
1au	108991	Fixed Height Rigid Backrest 13.5", 14" Wide <i>Used with 15" back height with omit push handles or 16" back height with standard or low profile push handles</i>	1
1av	108990	Fixed Height Rigid Backrest 14.5", 14" Wide <i>Used with 16" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1aw	108989	Fixed Height Rigid Backrest 15.5", 14" Wide <i>Used with 17" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1ax	108988	Fixed Height Rigid Backrest 16.5", 14" Wide <i>Used with 18" back height with omit push handles</i>	1
1ay	109008	Fixed Height Rigid Backrest 7.5", 15" Wide <i>Used with 9" back height with omit push handles, 13" back height with fold down push handles or 10" back height with standard or low profile push handles</i>	1
1az	109007	Fixed Height Rigid Backrest 8.5", 15" Wide <i>Used with 10" back height with omit push handles, 14" back height with fold down push handles or 11" back height with standard or low profile push handles</i>	1
1ba	109006	Fixed Height Rigid Backrest 9.5", 15" Wide <i>Used with 11" back height with omit push handles, 15" back height with fold down push handles or 12" back height with standard or low profile push handles</i>	1

Item Number	Part Number	Description	Quantity
1bb	109005	Fixed Height Rigid Backrest 10.5", 15" Wide <i>Used with 12" back height with omit push handles, 16" back height with fold down push handles or 13" back height with standard or low profile push handles</i>	1
1bc	109004	Fixed Height Rigid Backrest 11.5", 15" Wide <i>Used with 13" back height with omit push handles, 17" back height with fold down push handles or 14" back height with standard or low profile push handles</i>	1
1bd	109003	Fixed Height Rigid Backrest 12.5", 15" Wide <i>Used with 14" back height with omit push handles, 18" back height with fold down push handles or 15" back height with standard or low profile push handles</i>	1
1be	109002	Fixed Height Rigid Backrest 13.5", 15" Wide <i>Used with 15" back height with omit push handles or 16" back height with standard or low profile push handles</i>	1
1bf	109001	Fixed Height Rigid Backrest 14.5", 15" Wide <i>Used with 16" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1bg	109000	Fixed Height Rigid Backrest 15.5", 15" Wide <i>Used with 17" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1bh	108999	Fixed Height Rigid Backrest 16.5", 15" Wide <i>Used with 18" back height with omit push handles</i>	1
1bi	109019	Fixed Height Rigid Backrest 7.5", 16" Wide <i>Used with 9" back height with omit push handles, 13" back height with fold down push handles or 10" back height with standard or low profile push handles</i>	1
1bj	109018	Fixed Height Rigid Backrest 8.5", 16" Wide <i>Used with 10" back height with omit push handles, 14" back height with fold down push handles or 11" back height with standard or low profile push handles</i>	1
1bk	109017	Fixed Height Rigid Backrest 9.5", 16" Wide <i>Used with 11" back height with omit push handles, 15" back height with fold down push handles or 12" back height with standard or low profile push handles</i>	1
1bl	109016	Fixed Height Rigid Backrest 10.5", 16" Wide <i>Used with 12" back height with omit push handles, 16" back height with fold down push handles or 13" back height with standard or low profile push handles</i>	1
1bm	109015	Fixed Height Rigid Backrest 11.5", 16" Wide <i>Used with 13" back height with omit push handles, 17" back height with fold down push handles or 14" back height with standard or low profile push handles</i>	1
1bn	109014	Fixed Height Rigid Backrest 12.5", 16" Wide <i>Used with 14" back height with omit push handles, 18" back height with fold down push handles or 15" back height with standard or low profile push handles</i>	1

Item Number	Part Number	Description	Quantity
1bo	109013	Fixed Height Rigid Backrest 13.5", 16" Wide <i>Used with 15" back height with omit push handles or 16" back height with standard or low profile push handles</i>	1
1bp	109012	Fixed Height Rigid Backrest 14.5", 16" Wide <i>Used with 16" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1bq	109011	Fixed Height Rigid Backrest 15.5", 16" Wide <i>Used with 17" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1br	109010	Fixed Height Rigid Backrest 16.5", 16" Wide <i>Used with 18" back height with omit push handles</i>	1
1bs	109030	Fixed Height Rigid Backrest 7.5", 17" Wide <i>Used with 9" back height with omit push handles, 13" back height with fold down push handles or 10" back height with standard or low profile push handles</i>	1
1bt	109029	Fixed Height Rigid Backrest 8.5", 17" Wide <i>Used with 10" back height with omit push handles, 14" back height with fold down push handles or 11" back height with standard or low profile push handles</i>	1
1bu	109028	Fixed Height Rigid Backrest 9.5", 17" Wide <i>Used with 11" back height with omit push handles, 15" back height with fold down push handles or 12" back height with standard or low profile push handles</i>	1
1bv	109027	Fixed Height Rigid Backrest 10.5", 17" Wide <i>Used with 12" back height with omit push handles, 16" back height with fold down push handles or 13" back height with standard or low profile push handles</i>	1
1bw	109026	Fixed Height Rigid Backrest 11.5", 17" Wide <i>Used with 13" back height with omit push handles, 17" back height with fold down push handles or 14" back height with standard or low profile push handles</i>	1
1bx	109025	Fixed Height Rigid Backrest 12.5", 17" Wide <i>Used with 14" back height with omit push handles, 18" back height with fold down push handles or 15" back height with standard or low profile push handles</i>	1
1by	109024	Fixed Height Rigid Backrest 13.5", 17" Wide <i>Used with 15" back height with omit push handles or 16" back height with standard or low profile push handles</i>	1
1bz	109023	Fixed Height Rigid Backrest 14.5", 17" Wide <i>Used with 16" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1ca	109022	Fixed Height Rigid Backrest 15.5", 17" Wide <i>Used with 17" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1cb	109021	Fixed Height Rigid Backrest 16.5", 17" Wide <i>Used with 18" back height with omit push handles</i>	1

Item Number	Part Number	Description	Quantity
1cc	109041	Fixed Height Rigid Backrest 7.5", 18" Wide <i>Used with 9" back height with omit push handles, 13" back height with fold down push handles or 10" back height with standard or low profile push handles</i>	1
1cd	109040	Fixed Height Rigid Backrest 8.5", 18" Wide <i>Used with 10" back height with omit push handles, 14" back height with fold down push handles or 11" back height with standard or low profile push handles</i>	1
1ce	109039	Fixed Height Rigid Backrest 9.5", 18" Wide <i>Used with 11" back height with omit push handles, 15" back height with fold down push handles or 12" back height with standard or low profile push handles</i>	1
1cf	109038	Fixed Height Rigid Backrest 10.5", 18" Wide <i>Used with 12" back height with omit push handles, 16" back height with fold down push handles or 13" back height with standard or low profile push handles</i>	1
1cg	109037	Fixed Height Rigid Backrest 11.5", 18" Wide <i>Used with 13" back height with omit push handles, 17" back height with fold down push handles or 14" back height with standard or low profile push handles</i>	1
1ch	109036	Fixed Height Rigid Backrest 12.5", 18" Wide <i>Used with 14" back height with omit push handles, 18" back height with fold down push handles or 15" back height with standard or low profile push handles</i>	1
1ci	109035	Fixed Height Rigid Backrest 13.5", 18" Wide <i>Used with 15" back height with omit push handles or 16" back height with standard or low profile push handles</i>	1
1cj	109034	Fixed Height Rigid Backrest 14.5", 18" Wide <i>Used with 16" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1ck	109033	Fixed Height Rigid Backrest 15.5", 18" Wide <i>Used with 17" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1cl	109032	Fixed Height Rigid Backrest 16.5", 18" Wide <i>Used with 18" back height with omit push handles</i>	1
1cm	109052	Fixed Height Rigid Backrest 7.5", 19" Wide <i>Used with 9" back height with omit push handles, 13" back height with fold down push handles or 10" back height with standard or low profile push handles</i>	1
1cn	109051	Fixed Height Rigid Backrest 8.5", 19" Wide <i>Used with 10" back height with omit push handles, 14" back height with fold down push handles or 11" back height with standard or low profile push handles</i>	1
1co	109050	Fixed Height Rigid Backrest 9.5", 19" Wide <i>Used with 11" back height with omit push handles, 15" back height with fold down push handles or 12" back height with standard or low profile push handles</i>	1



Item Number	Part Number	Description	Quantity
1cp	109049	Fixed Height Rigid Backrest 10.5", 19" Wide <i>Used with 12" back height with omit push handles, 16" back height with fold down push handles or 13" back height with standard or low profile push handles</i>	1
1cq	109048	Fixed Height Rigid Backrest 11.5", 19" Wide <i>Used with 13" back height with omit push handles, 17" back height with fold down push handles or 14" back height with standard or low profile push handles</i>	1
1cr	109047	Fixed Height Rigid Backrest 12.5", 19" Wide <i>Used with 14" back height with omit push handles, 18" back height with fold down push handles or 15" back height with standard or low profile push handles</i>	1
1cs	109046	Fixed Height Rigid Backrest 13.5", 19" Wide <i>Used with 15" back height with omit push handles or 16" back height with standard or low profile push handles</i>	1
1ct	109045	Fixed Height Rigid Backrest 14.5", 19" Wide <i>Used with 16" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1cu	109044	Fixed Height Rigid Backrest 15.5", 19" Wide <i>Used with 17" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1cv	109043	Fixed Height Rigid Backrest 16.5", 19" Wide <i>Used with 18" back height with omit push handles</i>	1
1cw	109063	Fixed Height Rigid Backrest 7.5", 20" Wide <i>Used with 9" back height with omit push handles, 13" back height with fold down push handles or 10" back height with standard or low profile push handles</i>	1
1cx	109062	Fixed Height Rigid Backrest 8.5", 20" Wide <i>Used with 10" back height with omit push handles, 14" back height with fold down push handles or 11" back height with standard or low profile push handles</i>	1
1cy	109061	Fixed Height Rigid Backrest 9.5", 20" Wide <i>Used with 11" back height with omit push handles, 15" back height with fold down push handles or 12" back height with standard or low profile push handles</i>	1
1cz	109060	Fixed Height Rigid Backrest 10.5", 20" Wide <i>Used with 12" back height with omit push handles, 16" back height with fold down push handles or 13" back height with standard or low profile push handles</i>	1
1da	109059	Fixed Height Rigid Backrest 11.5", 20" Wide <i>Used with 13" back height with omit push handles, 17" back height with fold down push handles or 14" back height with standard or low profile push handles</i>	1
1db	109058	Fixed Height Rigid Backrest 12.5", 20" Wide <i>Used with 14" back height with omit push handles, 18" back height with fold down push handles or 15" back height with standard or low profile push handles</i>	1

Item Number	Part Number	Description	Quantity
1dc	109057	Fixed Height Rigid Backrest 13.5", 20" Wide <i>Used with 15" back height with omit push handles or 16" back height with standard or low profile push handles</i>	1
1dd	109056	Fixed Height Rigid Backrest 14.5", 20" Wide <i>Used with 16" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1de	109055	Fixed Height Rigid Backrest 15.5", 20" Wide <i>Used with 17" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1df	109054	Fixed Height Rigid Backrest 16.5", 20" Wide <i>Used with 18" back height with omit push handles</i>	1
1dg	109074	Fixed Height Rigid Backrest 7.5", 21" Wide <i>Used with 9" back height with omit push handles, 13" back height with fold down push handles or 10" back height with standard or low profile push handles</i>	1
1dh	109073	Fixed Height Rigid Backrest 8.5", 21" Wide <i>Used with 10" back height with omit push handles, 14" back height with fold down push handles or 11" back height with standard or low profile push handles</i>	1
1di	109072	Fixed Height Rigid Backrest 9.5", 21" Wide <i>Used with 11" back height with omit push handles, 15" back height with fold down push handles or 12" back height with standard or low profile push handles</i>	1
1dj	109071	Fixed Height Rigid Backrest 10.5", 21" Wide <i>Used with 12" back height with omit push handles, 16" back height with fold down push handles or 13" back height with standard or low profile push handles</i>	1
1dk	109070	Fixed Height Rigid Backrest 11.5", 21" Wide <i>Used with 13" back height with omit push handles, 17" back height with fold down push handles or 14" back height with standard or low profile push handles</i>	1
1dl	109069	Fixed Height Rigid Backrest 12.5", 21" Wide <i>Used with 14" back height with omit push handles, 18" back height with fold down push handles or 15" back height with standard or low profile push handles</i>	1
1dm	109068	Fixed Height Rigid Backrest 13.5", 21" Wide <i>Used with 15" back height with omit push handles or 16" back height with standard or low profile push handles</i>	1
1dn	109067	Fixed Height Rigid Backrest 14.5", 21" Wide <i>Used with 16" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1do	109066	Fixed Height Rigid Backrest 15.5", 21" Wide <i>Used with 17" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1dp	109065	Fixed Height Rigid Backrest 16.5", 21" Wide <i>Used with 18" back height with omit push handles</i>	1

Item Number	Part Number	Description	Quantity
1dq	109085	Fixed Height Rigid Backrest 7.5", 22" Wide <i>Used with 9" back height with omit push handles, 13" back height with fold down push handles or 10" back height with standard or low profile push handles</i>	1
1dr	109084	Fixed Height Rigid Backrest 8.5", 22" Wide <i>Used with 10" back height with omit push handles, 14" back height with fold down push handles or 11" back height with standard or low profile push handles</i>	1
1ds	109083	Fixed Height Rigid Backrest 9.5", 22" Wide <i>Used with 11" back height with omit push handles, 15" back height with fold down push handles or 12" back height with standard or low profile push handles</i>	1
1dt	109082	Fixed Height Rigid Backrest 10.5", 22" Wide <i>Used with 12" back height with omit push handles, 16" back height with fold down push handles or 13" back height with standard or low profile push handles</i>	1
1du	109081	Fixed Height Rigid Backrest 11.5", 22" Wide <i>Used with 13" back height with omit push handles, 17" back height with fold down push handles or 14" back height with standard or low profile push handles</i>	1
1dv	109080	Fixed Height Rigid Backrest 12.5", 22" Wide <i>Used with 14" back height with omit push handles, 18" back height with fold down push handles or 15" back height with standard or low profile push handles</i>	1
1dw	109079	Fixed Height Rigid Backrest 13.5", 22" Wide <i>Used with 15" back height with omit push handles or 16" back height with standard or low profile push handles</i>	1
1dx	109078	Fixed Height Rigid Backrest 14.5", 22" Wide <i>Used with 16" back height with omit push handles or 17" back height with standard or low profile push handles</i>	1
1dy	109077	Fixed Height Rigid Backrest 15.5", 22" Wide <i>Used with 17" back height with omit push handles or 18" back height with standard or low profile push handles</i>	1
1dz	109076	Fixed Height Rigid Backrest 16.5", 22" Wide <i>Used with 18" back height with omit push handles</i>	1
2	100804	Tube End Plug 1" <i>Used with omit push handles</i>	2
3	000034	Sleeve 3/4" x 1 1/2" <i>Used with standard push handles, low profile push handles or fold down push handles</i>	2