



Raz-SP

Self-Propel

Mobile Shower Commode Chair

The **RAZ-SP** (Self-Propel) is available in 4 widths, 2 lengths and 2 height ranges. The stainless steel frame is designed with unobstructed side as well as front access for the user or attendant.



FEATURES



ALSO AVAILABLE

- Long frame — 3" (8 cm) longer for tall clients
- Low frame - 4" seat-to-floor height
- 16" and 20" (40 and 51 cm) widths available
- Reclining Back

Extensive range of options available

SPECIFICATIONS	Raz-SP 22" wheels	Raz-SP 24" wheels
WEIGHT CAPACITY	400 lbs (181 kg)	400 lbs (181 kg)
FRAME WIDTH (seat and back frame)	Standard: 18" (46 cm)	Standard: 18" (46 cm)
OVERALL WIDTH (with standard armpad)	26.5" (67 cm)	26.5" (67 cm)
OVERALL LENGTH (with footrest)	39" (99 cm)	39" (99 cm)
SEAT WIDTH (commode seat)	18" (46 cm)	18" (46 cm)
SEAT DEPTH (front of back post to front of seat)	16" or 19" (40.5 or 48 cm)	16" or 19" (40.5 or 48 cm)
CLEARANCE HEIGHT UNDER SEAT	16", 17", 18", 20" (41 cm, 43 cm, 46 cm, 51 cm) without pan 15.5", 16.5", 17.5", 19.5" (39 cm, 42 cm, 44 cm, 50 cm) with pan	15.75", 16.75", 17.75", 19.75" (40 cm, 43 cm, 45 cm, 50 cm) without pan brackets 15.25", 16.25", 17.25", 19.25" (39 cm, 41 cm, 44 cm, 49 cm) with pan brackets
ARMREST HEIGHT ABOVE SEAT (17" Adjusta-Back)	7" - 10" in ½" increments (18 - 25 cm)	7" - 10" in ½" increments (18 - 25 cm)
SEAT-TO-FLOOR HEIGHT	20" - 24" in 1" increments (51 - 61 cm)	20" - 24" in 1" increments (51 - 61 cm)
PRODUCT WEIGHT (with footrests)	47 lbs (21.3 kg)	49 lbs (22kg)
FRAME AND HARDWARE MATERIAL	Stainless Steel (Lifetime warranty against corrosion)	Stainless Steel (Lifetime warranty against corrosion)

All measurements taken with standard molded seat and standard arm pad.

Two-year limited warranty covers the frame against manufacturing defects. A one-year limited warranty covers all other components, subject to normal wear and tear exclusions.

Specifications, and options are subject to change without prior notice.

MANUFACTURED BY:

RAZ Design Inc. 135 Railside Road, Toronto, Ontario Canada, M3A 1B7 T: (877) 720-5678 E: info@razdesigninc.com

5685HYDWH