

## Rogue Tires - Sentinel

---



<p>

Utilize the Shox Tire as a comparable alternative if using on a Mag wheel. <p>

Please note tab mounted handdrums must be installed prior to the tire being installed.

Item Number	Part Number	Description	Quantity	Retail
1a	102974	Tire 20" x 1" Sentinel Black (25-451) <i>Used on 20" Spoke and Spinergy Wheels</i>	1	\$52.53
1b	103165	Tire 22" x 1" Sentinel Black (25-501) <i>Used on 22" Spoke and Spinergy Wheels</i>	1	\$52.53
1c	102976	Tire 24" x 1" Sentinel Black (25-540) <i>Used on 24" Spoke and Spinergy Wheels</i>	1	\$68.29
1d	102977	Tire 25" x 1" Sentinel Black (25-559) <i>Used on 25" Spoke and Spinergy Wheels</i>	1	\$68.29
1e	102978	Tire 26" x 1" Sentinel Black (25-590) <i>Used on 26" Spoke and Spinergy Wheels</i>	1	\$68.29