

# 16" Mag Wheel



*If a drum brake is present, see "16" Mag Wheel (Focus Drum Brake)" or "16" Mag Wheel (Flip Drum Brake)" page for correct wheel assembly. If an Attendant Foot Lock is present on a Flip, please see the 12" & 16" Mag Wheel (Attendant Foot Lock) parts manual page. If changing wheel size and the Attendant Foot Lock is present on a Focus, a new linkage may be required; reference the corresponding Focus CR Attendant Foot Lock page. If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles, or axle plate and wheel locks may be needed.*

Item Number	Part Number	Description	Quantity
1a	001159	16" Mag Wheel w/ Pneumatic Assembly Tilt	1
1b	001160	16" Mag Wheel w/ Pneumatic w/ Airless Insert Assembly Tilt	1
2	001161	16" Mag Wheel w/ Full Poly Assembly Tilt	1
3	003453	16" Mag Wheel w/ Low Poly Assembly Tilt	1