

## Best Value and Choice for More Patient Care Environments

The Hoyer® HPL450 series lift is a flexible mobile patient lift that is the perfect crossover product for both nursing and home care. Designed with upscale features, high lifting capacity, and lightweight design allowing the HPL450 to work well in almost all settings. The HPL450 provides additional value with optional electric leg and HML450 hydraulic versions.



### OVERSIZED PUSH HANDLE



Optimally positioned and ergonomically designed for ease of maneuverability.

### FOOT PEDAL LEG ADJUSTMENT



The ergonomically designed foot pedal is perfectly positioned to ensure the user's foot remains level during operation, helping to avoid risk of injury.

### OPTIONAL POWERED BASE



Powered leg adjustment allows the caregiver to open and close the legs without bending, twisting or using force.

### OPTIONAL SCALE



Simple to operate scale has a "zeroing out" feature that allows you to weigh the patient during transfer process, without including the weight of the sling.



**SAFE WORKING LOAD: 450 LB / 204 KG**

AVAILABLE CONFIGURATIONS	
Powered w/Electric Leg	HPL450-EL
Powered w/ Manual Leg	HPL450-ML
Hydraulic	HML450

SPECIFICATION		
Maximum Overall Length	51.1"	1300 mm
Maximum Overall Height	79.8"	2027 mm
Minimum Overall Height	54.3"	1380 mm
Legs Open - External Width	48.4"	1250 mm
Legs Open - Internal Width	44.9"	1140 mm
Legs Closed - External Width	25.0"	660 mm
Legs Closed - Internal Width	22.8"	580 mm
Overall Height of Legs	3.9"	100 mm
Spreader Bar - Maximum Height	73.5"	1867 mm
Spreader Bar - Minimum Height	22.2"	565 mm

WEIGHTS		
Mast, Powered Base and Boom	86.6 lb	39.3 kg
Power Pack	6.1 lb	2.8 kg
Powered Base Total	92.5 lb	42.0 kg
Mast, Manual Base and Boom	86.4 lb	39.2 kg
Manual Base Total	92.0 lb	41.7 kg

Specifications may vary due to manufacturing tolerances.