



600

Available in  
**Black or White**

## Rehab Shower/Commode Chair

- Packages Now Available
- Stainless Steel Frame
- Stainless Steel Frame- New Frame and Seat Sizes  
18"Wx18"D, 18"Wx20"D, 20"Wx18"D and 22"Wx18"D
- Three Mobility Choices:  
Four Casters, 20", or 24" Wheels
- 350 lb. Weight Capacity

All ActiveAid chairs are stainless steel for long lasting and durable use. The large selection of seats allows for customizing for the individual. The 600 has many available options so the chair can be fitted to the user or it could be set up to be used in a multi-use facility. From the least involved client to the most involved, this chair is readily able to fit their needs.

The 600 series shower/commode chairs have features that address the needs of the home based users, long-term care, home care clients, and care givers. This is the go to chair because of all the options and sizes available. From 5" casters to large wheels and many choices to pick from, this series is designed to make those daily activities just a little bit easier.



## FRAMES



**P13675**  
18"Wx18"D  
Stainless Steel  
Frame



**P13676**  
18"Wx20"D  
Stainless Steel  
Frame



**P13677**  
20"Wx18"D  
Stainless Steel  
Frame



**P13678**  
22"Wx18"D  
Stainless Steel  
Frame



**60800-2**  
5" Dual Locking  
Casters  
(pair)



**P60194**  
20" Wheels  
(pair)



**P60551**  
20" Pneumatic  
Street Tread  
Wheels  
(pair)

## AXLES



**P60196**  
24" Wheels  
(pair)



**P60549**  
24" Pneumatic  
Street Tread  
Wheels  
(pair)



**P61059**  
24" Flat Free Street  
Tread Wheels  
(pair)



**P61102**  
20" PVC Handrims  
(replaces standard  
handrims, pair)



**P61104**  
24" PVC Handrims  
(replaces standard  
handrims, pair)



**P13649**  
Bolted Axles with  
Wheel Mounts and  
Tipping Levers  
(pair)



**P13650**  
Quick Release Axles  
with Wheel Mounts  
and Tipping Levers  
(pair)

## WHEEL LOCKS



**P60532**  
Wheel Locks  
(pair)



**P13135**  
Wheel Locks with  
Extended Handle  
(pair)



**60841-2**  
5" Casters  
(pair)



**P13135**  
5" Dual Locking  
Casters  
(pair)



**P60993\***  
18"Wx18"D  
Comfortuff 1.5" waterfall foam  
Front/Rear Open Seat  
A-5.5"xB-8" opening



**P60974**  
18"Wx18"D  
Comfortuff sloped foam  
Front Open Seat  
A-4"xB-6.5" opening



**P60701\***  
18"Wx18"D  
Ensolite® 1" foam  
Front/Rear Open Seat  
A-4.55"xB-7.75" opening

\*Seat works in all 4 positions

## FRONT CASTERS

## SEATS



**P61058**  
18"Wx18"D  
Ensolite® 2" waterfall foam  
Front Open Contoured Seat  
A-5.5"xB-8.25" opening



**P61086**  
18"Wx18"D  
Ensolite® 2" foam  
Solid Seat



**P60675**  
18"Wx18"D  
Ensolite® 1" foam  
Cloverleaf Seat  
A-7.5"xB-9" oval opening



**P60913\***  
18"Wx18"D  
Ensolite® 2" waterfall foam  
Front/Rear Open Seat  
A-5.25"xB-8" opening



**P61121**  
18"Wx18"D  
Ensolite® 1" foam  
Front/Rear Open Seat  
A-4.5"xB-7.5" opening



**PP60903**  
18"Wx20"D  
Ensolite® 1" foam  
Side Open Seat  
A-4.5"xB-7.75" opening



**P60903**  
20"Wx18"D  
Ensolite® 1" foam  
Front/Rear Open Seat  
A-4.5"xB-7.75" opening

## ARMRESTS

## LEG RESTS



**P61038**  
20"Wx18"D  
Comfortuff 1.5" waterfall foam  
Front Open Seat  
A-5"xB-8" opening



**P60905**  
22"Wx18"D  
Ensolite® 1" foam  
Front/Rear Open Seat  
A-4.5"xB-7.75" opening




**P13278**  
18"D Left Hand  
Armrest with Pad



**P13277**  
18"D Right Hand  
Armrest with Pad



**P13663**  
20"D Left Hand  
Armrest with Pad



**P13664**  
20"D Right Hand  
Armrest with Pad



**P13665**  
Left Leg Rest



**P13666**  
Right Leg Rest



**P13582**  
Left Elevating Leg  
Rest



**P13583**  
Right Elevating Leg  
Rest



**PL13581**  
Left Gooseneck  
Leg Rest



**PR13581**  
Right Gooseneck  
Leg Rest




**P13667**  
18"W 1-pc. Foot  
Plate/Leg Rest



**P13668**  
20"W 1-pc. Foot  
Plate/Leg Rest

## FOOT PLATES



**P13669**  
22"W 1-pc. Foot  
Plate/Leg Rest



**PL11099**  
Left Foot Plate



**PR11099**  
Right Foot Plate



**PL11100**  
Left Foot Plate with  
Heel Loop



**PR11100**  
Right Foot Plate  
with Heel Loop



**P61137**  
Left Oversized  
Aluminum  
Foot Plate



**P61138**  
Right Oversized  
Aluminum  
Foot Plate

## BACKS



**PL11101**  
Left Padded  
Foot Plate



**PR11101**  
Right Padded  
Foot Plate



**P13651**  
18"W Sling Back  
with Stroller  
Handles



**P13652**  
18"W Solid Back



**P13658**  
20"W Sling Back  
with Stroller  
Handles



**P13659**  
20"W Solid Back



**P13661**  
22"W Sling Back  
with Stroller  
Handles

## BELTS



**P13713**  
22"W Solid Back



**PW60909**  
Waist Belt  
(48" circumference with  
12" velcro overlap)



**PT60909**  
Trunk Belt  
(48" circumference with  
12" velcro overlap)



**PW11094**  
Waist Belt  
(for 20"W and 22"W frame  
only, 48" circumference  
with 12" velcro overlap)



**PT11094**  
Trunk Belt  
(for 20"W and 22"W frame  
only, 48" circumference  
with 12" velcro overlap)



**P60879**  
Calf Strap



**P11093**  
Bodypoint®  
Trunk Belt



**P11002**  
Bodypoint®  
Trunk Belt  
(pair)



**P11003**  
Bodypoint® Trunk  
Belt & Waist Belt



**P11112**  
Bodypoint®  
XL Trunk Belt



**P11113**  
Bodypoint®  
Belt w/Quick  
Release Knobs



**P61081**  
Bodypoint®  
Calf Strap  
(only available with 18" back)



**P11114**  
Bodypoint®  
Knee Belt

## POSITIONING AIDS



**PL13574**  
Left Hand  
Swing-Away  
Adjustable Lateral



**PR13574**  
Right Hand  
Swing-Away  
Adjustable Lateral



**PL13618**  
Left Hand  
Arm-Mounted Hip  
Guide



**PR13618**  
Right Hand  
Arm-Mounted Hip  
Guide



**P61006**  
Abductor/Deflector  
Large



**P13670**  
Front  
Anti-Tippers

## PANS



**P13167**  
18"Wx18"D  
Pan/Hanger



**P13672**  
18"Wx20"D  
Pan/Hanger




**P13673**  
20"Wx18"D  
Pan/Hanger




**P13674**  
22"Wx18"D  
Pan/Hanger

Seat Opening  
Measurements



Cloverleaf Seat Opening  
Measurements



full option details at  
[activeaid.com](http://activeaid.com)





600 with 20" Wheels

## SPECIFICATIONS

Weight Capacity	350 lbs.
Back Height from Seat	16.5"
Seat Width Between Arms	18.5", 20.5", 22.5"
Seat Width	18", 20", 22"
Seat Depth	18", 20"
Floor to Bottom of Seat Frame	19.5"
Width Between Arm Pads	20.5"
Top of Seat to Top of Armrest Pads	8"
Overall Width	24.75"
Maximum Back Slack	1"
Minimum Overall length with Casters	40"-41"
Maximum Overall Length with Footrest with Wheels	42"-43"
Leg Rest Range	9"-21"
Gooseneck Leg Rest	16.5"-19.5"
Under Seat Width Clearance- Casters	16"
Wheels	14.5"