<u>ziggo</u>

THE LIGHTWEIGHT, FOLDING MANUAL WHEELCHAIR TO TAKE YOU PLACES



AVAILABLE IN 3 SIZES:

12", 14" & 16" (WIDTH)

DESIGNED FOR EASY

TRANSPORT & STORAGE

LET THE ZiggO TAKE YOU PLACES

INDOORS AND OUTDOORS



718.943.6452 www. circlespecialty.com info@circlespecialty.com

THE ZiggO is a lightweight manual wheelchair

designed for kids and teens **on-the-go**. Loved by users and caregivers alike, the <code>ziggO</code> is practical, durable, adjustable and cost effective, allowing for maximum maneuverability and easy transport, be it indoors and outdoors. Meeting the standards of therapists and clinicians, the <code>ziggO</code> provides comfortable support, functionality and positioning options.

Experience the little things in life with great ease.

Colored steel frame

Folds easily for compact storage and transport

Comfortable nylon upholstery with back pocket

Flip-back padded desk length arms

Fold-down backrest on 16" Ziggo model.

Pelvic belt

Quick-release treaded PU rear wheels

Solid PVC front wheels

Push-to-lock brakes

Height adjustable swing-away footrests with pastic footplates

Calf strap for extra support and comfort

Variety of optional accessories to meet your individual needs

3-year limited warranty



	ZG1200	ZG1400	ZG1600
Seat Width	12"	14"	16"
Seat Depth	12"	16"	16"
Back Height	14"	16"	16"
Seat Height Front	18"	18.5"	18.5"
Seat Height Back	17"	17.5"	17.5"
Overall Open Width	18.5"	21"	23"
Overall Length w/Riggings	36"	41"	42"
Weight w/o Riggings	28 lbs	36.5 lbs	37 lbs
Weight Capacity	150 lbs	150 lbs	250 lbs

SPECIFICATIONS

Wheelchair accessories are designed to increase comfort and improve support and mobility. Based on your requirements and circumstances, you can purchase accessories catered to your needs and preferences.

Installation is simple, giving you the ease to remove or add accessories as required.



Height, length and angle adjustable swing-away elevating legrests with padded calf pads.

ZGELR01 for use with ZG1200

ZGELRO2 for use with ZG1400 & ZG1600

Padded trunk harness with adjustable straps provides anterior torso support for improved overall positioning.

ZGTH01 (use with all Ziggo products)





Increases wheelchair stability and prevents wheelchair from tipping over backwards.

ZGAT01 (use with all Ziggo products)

Height, depth and angle adjustable headrest provides posterior head and neck support.





Lightweight foam cushion improves comfort and sitting tolerance.

ZGCS12 for use with ZG1200

ZGCS14 for use with ZG1400

ZGCS16 for use with ZG1600

ZG8100 (use with all Ziggo products)