


# Discontinued 2 Rigid Tension Adjustable Back Upholstery



Quantity Chart			
Back Height	Item 6	Item 5	Item 4
10"-13.5"	1	2	2

Quantity Chart			
Back Height	Item 6	Item 5	Item 4
14"-16.5"	2	2	4

Quantity Chart			
Back Height	Item 6	Item 5	Item 4
17"-20"	3	2	6

Open Back Panels are used when integrated push handles are present. Closed Back Panels are when backposts without push handles are present. Tension Straps (7a-7d) are used with the combination of Fold Down Push Handle and Tension Adjustable Back Upholstery. Tension Straps 6a - 6f and 7a - 7d are based on chair width.

Item Number	Part Number	Description	Quantity
1a, 5a	504747	Backrest Cover with Two Straps 10W Rigid, 9-12H <i>Contains 002485 cover</i>	1
1b, 5a	504180	Backrest Cover with Two Straps 10W Rigid, 11-14H <i>Contains 002486 cover</i>	1
1c, 5a	504181	Backrest Cover with Two Straps 10W Rigid, 13-16H <i>Contains 002487 cover</i>	1
1d, 5a	504182	Backrest Cover with Two Straps 10W Rigid, 15-18H <i>Contains 002488 cover</i>	1
1e, 5b	504748	Backrest Cover with Two Straps 11W Rigid, 9-12H <i>Contains 002489 cover</i>	1
1f, 5b	504183	Backrest Cover with Two Straps 11W Rigid, 11-14H <i>Contains 002490 cover</i>	1

Item Number	Part Number	Description	Quantity
1g, 5b	504184	Backrest Cover with Two Straps 11W Rigid, 13-16H <i>Contains 002491 cover</i>	1
1h, 5b	504185	Backrest Cover with Two Straps 11W Rigid, 15-18H <i>Contains 002492 cover</i>	1
1i, 5c	504749	Backrest Cover with Two Straps 12W Rigid, 9-12H <i>Contains 002493 cover</i>	1
1j, 5c	504186	Backrest Cover with Two Straps 12W Rigid, 11-14H <i>Contains 002494 cover</i>	1
1k, 5c	504187	Backrest Cover with Two Straps 12W Rigid, 13-16H <i>Contains 002495 cover</i>	1
1l, 5c	504188	Backrest Cover with Two Straps 12W Rigid, 15-20H <i>Contains 002496 cover</i>	1
1m, 5d	504750	Backrest Cover with Two Straps 13W Rigid, 9-12H <i>Contains 002497 cover</i>	1
1n, 5d	504189	Backrest Cover with Two Straps 13W Rigid, 11-14H <i>Contains 002498 cover</i>	1
1o, 5d	504190	Backrest Cover with Two Straps 13W Rigid, 13-16H <i>Contains 002499 cover</i>	1
1p, 5d	504191	Backrest Cover with Two Straps 13W Rigid, 15-18H <i>Contains 002500 cover</i>	1
1q, 5d	504192	Backrest Cover with Two Straps 13W Rigid, 17-20H <i>Contains 003043 cover</i>	1
1r, 5e	504751	Backrest Cover with Two Straps 14W Rigid, 9-12H <i>Contains 002501 cover</i>	1
1s, 5e	504193	Backrest Cover with Two Straps 14W Rigid, 11-14H <i>Contains 002502 cover</i>	1
1t, 5e	504194	Backrest Cover with Two Straps 14W Rigid, 13-16H <i>Contains 002503 cover</i>	1
1u, 5e	504195	Backrest Cover with Two Straps 14W Rigid, 15-18H <i>Contains 002504 cover</i>	1
1v, 5e	504196	Backrest Cover with Two Straps 14W Rigid, 17-20H <i>Contains 003044 cover</i>	1
1w, 5f	504752	Backrest Cover with Two Straps 15W Rigid, 9-12H <i>Contains 002505 cover</i>	1
1x, 5f	504197	Backrest Cover with Two Straps 15W Rigid, 11-14H <i>Contains 002506 cover</i>	1

Item Number	Part Number	Description	Quantity
1y, 5f	504198	Backrest Cover with Two Straps 15W Rigid, 13-16H <i>Contains 002507 cover</i>	1
1z, 5f	504199	Backrest Cover with Two Straps 15W Rigid, 15-18H <i>Contains 002508 cover</i>	1
1aa, 5f	504200	Backrest Cover with Two Straps 15W Rigid, 17-20H <i>Contains 003045 cover</i>	1
1ab, 5g	504753	Backrest Cover with Two Straps 16W Rigid, 9-12H <i>Contains 002509 cover</i>	1
1ac, 5g	504201	Backrest Cover with Two Straps 16W Rigid, 11-14H <i>Contains 002510 cover</i>	1
1ad, 5g	504202	Backrest Cover with Two Straps 16W Rigid, 13-16H <i>Contains 002511 cover</i>	1
1ae, 5g	504203	Backrest Cover with Two Straps 16W Rigid, 15-18H <i>Contains 002512 cover</i>	1
1af, 5g	504204	Backrest Cover with Two Straps 16W Rigid, 17-20H <i>Contains 003046 cover</i>	1
1ag, 5h	504754	Backrest Cover with Two Straps 17W Rigid, 9-12H <i>Contains 003048 cover</i>	1
1ah, 5h	504205	Backrest Cover with Two Straps 17W Rigid, 11-14H <i>Contains 003049 cover</i>	1
1ai, 5h	504206	Backrest Cover with Two Straps 17W Rigid, 13-16H <i>Contains 003050 cover</i>	1
1aj, 5h	504207	Backrest Cover with Two Straps 17W Rigid, 15-18H <i>Contains 003051 cover</i>	1
1ak, 5h	504208	Backrest Cover with Two Straps 17W Rigid, 17-20H <i>Contains 003052 cover</i>	1
1al, 5i	504755	Backrest Cover with Two Straps 18W Rigid, 9-12H <i>Contains 003054 cover</i>	1
1am, 5i	504209	Backrest Cover with Two Straps 18W Rigid, 11-14H <i>Contains 003055 cover</i>	1
1an, 5i	504210	Backrest Cover with Two Straps 18W Rigid, 13-16H <i>Contains 003056 cover</i>	1
1ao, 5i	504211	Backrest Cover with Two Straps 18W Rigid, 15-18H <i>Contains 003057 cover</i>	1
1ap, 5i	504212	Backrest Cover with Two Straps 18W Rigid, 17-20H <i>Contains 003058 cover</i>	1

Item Number	Part Number	Description	Quantity
1aq, 5j	504756	Backrest Cover with Two Straps 19W Rigid, 9-12H <i>Contains 003060 cover</i>	1
1ar, 5j	504213	Backrest Cover with Two Straps 19W Rigid, 11-14H <i>Contains 003061 cover</i>	1
1as, 5j	504214	Backrest Cover with Two Straps 19W Rigid, 13-16H <i>Contains 003062 cover</i>	1
1at, 5j	504215	Backrest Cover with Two Straps 19W Rigid, 15-18H <i>Contains 003063 cover</i>	1
1au, 5j	504216	Backrest Cover with Two Straps 19W Rigid, 17-20H <i>Contains 003064 cover</i>	1
1av, 5k	504757	Backrest Cover with Two Straps 20W Rigid, 9-12H <i>Contains 003066 cover</i>	1
1aw, 5k	504217	Backrest Cover with Two Straps 20W Rigid, 11-14H <i>Contains 003067 cover</i>	1
1ax, 5k	504218	Backrest Cover with Two Straps 20W Rigid, 13-16H <i>Contains 003068 cover</i>	1
1ay, 5k	504219	Backrest Cover with Two Straps 20W Rigid, 15-18H <i>Contains 003069 cover</i>	1
1az, 5k	504220	Backrest Cover with Two Straps 20W Rigid, 17-20H <i>Contains 003070 cover</i>	1
1ba, 5l	504221	Backrest Cover with Two Straps 21W Rigid, 13-16H <i>Contains 003072 cover</i>	1
1bb, 5l	504222	Backrest Cover with Two Straps 21W Rigid, 15-18H <i>Contains 003073 cover</i>	1
1bc, 5l	504223	Backrest Cover with Two Straps 21W Rigid, 17-20H <i>Contains 003074 cover</i>	1
1bd, 5m	504224	Backrest Cover with Two Straps 22W Rigid, 13-16H <i>Contains 003076 cover</i>	1
1be, 5m	504225	Backrest Cover with Two Straps 22W Rigid, 15-18H <i>Contains 003077 cover</i>	1
1bf, 5m	504226	Backrest Cover with Two Straps 22W Rigid, 17-20H <i>Contains 003078 cover</i>	1
4	100973	Back Upholstery D-Ring	4
5a	003483	Back Upholstery Strap 7.3" <i>Used on 10" chair widths</i>	2

Item Number	Part Number	Description	Quantity
5b	003484	Back Upholstery Strap 9.3" <i>Used on 11" chair widths</i>	2
5c	003485	Back Upholstery Strap 11.3" <i>Used on 12" chair widths</i>	2
5d	003486	Back Upholstery Strap 13.3" <i>Used on 13" chair widths</i>	2
5e	003487	Back Upholstery Strap 15.3" <i>Used on 14" chair widths</i>	2
5f	002532	Back Upholstery Strap 17.3" <i>Used on 15" chair widths</i>	2
5g	002533	Back Upholstery Strap 19.3" <i>Used on 16" chair widths</i>	2
5h	002534	Back Upholstery Strap 21.3" <i>Used on 17" chair widths</i>	2
5i	002535	Back Upholstery Strap 23.3" <i>Used on 18" chair widths</i>	2
5j	002536	Back Upholstery Strap 25.3" <i>Used on 19" chair widths</i>	2
5k	002537	Back Upholstery Strap 27.3" <i>Used on 20" chair widths</i>	2
5l	002538	Back Upholstery Strap 29.3" <i>Used on 21" chair widths</i>	2
5m	002539	Back Upholstery Strap 31.3" <i>Used on 22" chair widths</i>	2
7a	001750	Tension Strap Folding Handle 14-15 Wide	1
7b	001751	Tension Strap Folding Handle 16-17 Wide	1
7c	001752	Tension Strap Folding Handle 18-19 Wide	1
7d	001753	Tension Strap Folding Handle 20-22 Wide	1
8	100734	M5 x 16 PHMS ZC	2