


Ethos Wheel Lock Clamps



Item Number	Part Number	Description	Quantity	Retail
1 - 4	003973	Wheel Lock Clamp Assembly, Rigid <i>Used with 1.25" or less rear wheel spacing</i>	1	\$12.87
3 - 7	112430	Wheel Lock Clamp Assembly with Set Screws, Extended, Rigid <i>Replaced part number 004162. Used with 1.5" or greater rear wheel spacing or if Height Adjustable T-Arms are present</i>	1	\$13.89
1	003974	Wheel Lock Clamp Threaded	1	\$5.56
2	003975	Wheel Lock Clamp C'Bored	1	\$5.56
3	101701	Washer 3.315/.350 x .425/.525 x .059/.067 18-8 F/W	1	\$2.58
4	103327	M8 x 20 LHCS Blk Zc	1	\$0.74
5	105015	M5 X 6 DIN914 SSS CONE PT PLN	2	\$2.58
7	004163	Wheel Lock Clamp Threaded Long	1	\$5.74