


# Catalyst 5Ti Swing Away Front Frame



Frame assemblies are based on Seat Depth. If switching to a new Seat Depth, new Height Adjustable T-Arms or Height Adjustable Flip Back T-Arms may be needed. <p>

If chair has Flip Back or the discontinued Height Adjustable Flip Back armrests, select a frame assembly with the Flip Back Arm Modification. This modification does not apply to Height Adjustable Flip Back T-Arms.<p>

If retrofitting from a Fixed Front Frame, new caster housings and attaching hardware is required.

Item Number	Part Number	Description	Quantity	Retail
3, 4, 5, 6	501005	Seat Hook and Support Hardware Kit, Standard, Catalyst <i>This kit contains quantity of 1 each</i>	1	\$17.40

Item Number	Part Number	Description	Quantity	Retail
1a, 2	100270	Frame Assembly, Catalyst Ti Front, Left 14	1	\$388.66
1b, 2	100271	Frame Assembly, Catalyst Ti Front, Left 15	1	\$524.44
1c, 2	100272	Frame Assembly, Catalyst Ti Front, Left 16	1	\$524.44
1d, 2	100273	Frame Assembly, Catalyst Ti Front, Left 17	1	\$524.44
1e, 2	501639	Frame Assembly, Catalyst Ti Front, Left 18	1	\$394.68
1f, 2	501640	Frame Assembly, Catalyst Ti Front, Left 19	1	\$431.27
1g, 2	100274	Frame Assembly, Catalyst Ti Front, Left 20	1	\$524.44
1h, 2	100279	Frame Assembly, Catalyst Ti Front, Right 14	1	\$524.44
1i, 2	100280	Frame Assembly, Catalyst Ti Front, Right 15	1	\$524.44
1j, 2	100281	Frame Assembly, Catalyst Ti Front, Right 16	1	\$524.44
1k, 2	100282	Frame Assembly, Catalyst Ti Front, Right 17	1	\$524.44
1l, 2	501641	Frame Assembly, Catalyst Ti Front, Right 18	1	\$394.68
1m, 2	501642	Frame Assembly, Catalyst Ti Front, Right 19	1	\$431.27
1n, 2	100283	Frame Assembly, Catalyst Ti Front, Right 20	1	\$524.44
1c, 2	114361	Frame Assembly, Catalyst Ti Front with Flip Back Arm Modification, Left 16	1	\$524.44
1d, 2	114362	Frame Assembly, Catalyst Ti Front with Flip Back Arm Modification, Left 17	1	\$524.44
1e, 2	114389	Frame Assembly, Catalyst Ti Front with Flip Back Arm Modification, Left 18	1	\$394.68
1f, 2	114390	Frame Assembly, Catalyst Ti Front with Flip Back Arm Modification, Left 19	1	\$431.27
1g, 2	114363	Frame Assembly, Catalyst Ti Front with Flip Back Arm Modification, Left 20	1	\$524.44
1j, 2	114366	Frame Assembly, Catalyst Ti Front with Flip Back Arm Modification, Right 16	1	\$524.44
1k, 2	114367	Frame Assembly, Catalyst Ti Front with Flip Back Arm Modification, Right 17	1	\$524.44
1l, 2	114391	Frame Assembly, Catalyst Ti Front with Flip Back Arm Modification, Right 18	1	\$394.68
1m, 2	114392	Frame Assembly, Catalyst Ti Front with Flip Back Arm Modification, Right 19	1	\$431.27
1n, 2	114368	Frame Assembly, Catalyst Ti Front with Flip Back Arm Modification, Right 20	1	\$524.44
2	100804	Tube End Plug 1"	1	\$1.58
3	100831	1/2" Washer OD .040" Thick	1	\$1.58
4	000003	Seat Hook, Catalyst	1	\$7.70

<b>Item Number</b>	<b>Part Number</b>	<b>Description</b>	<b>Quantity</b>	<b>Retail</b>
5	000004	Seat Hook Support, Catalyst	1	\$6.55
6	100668	M6 x 20 CLS10.9 BTN SHCS ZC w/Patch	1	\$1.58

<b>Code</b>	<b>Color Name</b>
<b>Matte</b>	
AGM	Army Green
EGM	Emerald Green
GGM	Graphite Gray
IBM	Iridescent Burgundy
INM	Iridescent Navy
RRM	Retro Red
RBM	Royal Blue
SBM	Shadow Black
SWM	Smooth White
<b>Gloss</b>	
BLK	Black
SKY	Blue Skies
PNK	Bubblegum Pink
CBL	Candy Blue
PUR	Candy Purple
RED	Candy Red
COP	Copper
FUS	Fuchsia
GSM	Granny Smith
ORG	Orange
FLM	Graphic Image Transfer
<b>Discontinued</b>	
BLU	Blue
BUR	Burgundy
EBL	Electric Blue
GRN	Forest Green
SAP	Sapphire
GRY	Slate Gray