



# Raz-CAT

## Compact Attendant Tilt

Rehab Shower Commode

The **Raz-CAT** (Compact Attendant Tilt) has all of the features of our Raz-AT but in a more compact size. For pediatrics, adolescents or smaller adults, the Raz-CAT provides a fit that is better suited to physically-smaller clients. To accommodate client growth, the Raz-CAT has adjustable seat depths and is available with 3 back-width options.



# FEATURES

Extensive range of options available

Flip-up padded arms

Unique oblong seat aperture accommodates seat depth adjustment

Exclusive to Raz, Ischial Pelvic Align System (IPAS), provides 2" (5 cm) of fore/aft or 1" (2.5 cm) of lateral seat adjustment

Stainless steel frame with 4" (10 cm) of height adjustment

Dual-locking (swivel and roll) casters



Tension-adjustable back upholstery

Available back frames: 14", 16", 18" (36 cm, 41 cm, 46 cm) wide 15", 18" (38 cm, 46 cm) tall

Backrest assembly adjusts fore/aft for seat depth change

Adjustable back angle

Locking gas struts to facilitate tilting

## ALSO AVAILABLE

- Custom seats
- Shown with optional laterals, hip adductors and heel loops

| SPECIFICATIONS                                     | Raz-CAT (Compact Attendant Tilt)   |
|--|--|
| WEIGHT CAPACITY                                    | 250 lbs (115 kg)   |
| FRAME WIDTH (seat and Back frame)                  | 18" (46 cm)  |
| BACK FRAME WIDTHS                                  | 14" (35 cm), 16" (40 cm), or 18" (46 cm)   |
| OVERALL WIDTH (with standard armpad)               | 21.75" (55 cm), 23" (58 cm) with 18" (46 cm) wide Back frame   |
| SEAT WIDTH (commode seat)                          | 18" (46 cm)  |
| BACK HEIGHTS (top of seat to top of upholstery)    | 15" or 18" (38 cm or 46 cm)  |
| RECLINING BACKS (top of seat to top of upholstery) | 15" or 18" (38 cm or 46 cm)  |
| SEAT DEPTH (front of back post to front of seat)   | 11" – 16" (28 cm – 40 cm)  |
| OVERALL LENGTH WHEN UPRIGHT                        | 32" (81 cm)  |
| CLEARANCE HEIGHT UNDER SEAT                        | 16.75", 17.75", 18.75", 19.75", 20.75" (42.5, 45, 47.5, 50, 53 cm) with pan brackets<br>17.5", 18.5", 19.5", 20.5", 21.5" (44.5, 47, 49.5, 52, 54.5 cm) without pan brackets |
| ARMREST HEIGHT ABOVE SEAT (20" Adjusta-Back)       | 15" (38 cm) High Back: 5.5" – 8" in ½" increments (14 – 20 cm)<br>18" (46 cm) High Back: 6.5" – 9" in ½" increments (16.5 – 23 cm)   |
| SEAT-TO-FLOOR HEIGHT                               | 21.5" – 25.5" in 1" increments (55 – 65 cm)  |
| PRODUCT WEIGHT                                     | 55 lbs (25 kg)   |
| FRAME AND HARDWARE MATERIAL                        | Stainless Steel (Lifetime warranty against corrosion)  |

All measurements taken with standard molded seat and standard arm pad.

Two-year limited warranty covers the frame against manufacturing defects. A one-year limited warranty covers all other components, subject to normal wear and tear exclusions.

Specifications, and options are subject to change without prior notice.

## MANUFACTURED BY:

RAZ Design Inc. 135 Railside Road, Toronto, Ontario Canada, M3A 1B7 T: (877) 720-5678 E: info@razdesigninc.com