


22" MAXX Mag Wheel / Tire / Handrim Kits



Not all components of the Wheel, Tire and Handrim kits will be listed. These kits cannot be used if Drum Brake is present. If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles and/or receivers may be needed. For the Catalyst 4, new axle plates may be needed. See chair specific axle page for axle and receiver utilization.

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|---|----------|
| | 118218 | 22" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount | 1 |

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|---|----------|
| | 118219 | 22" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Close Mount | 1 |
| | 118220 | 22" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip | 1 |
| | 118221 | 22" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip | 1 |
| | 118222 | 22" MAXX Mag, Pneumatic Tire (Black), Plastic Coated Handrim Kit, Standard Mount | 1 |
| | 118223 | 22" MAXX Mag, Pneumatic Tire (Black), Plastic Coated Handrim Kit, Close Mount | 1 |
| | 118224 | 22" MAXX Mag, Pneumatic Tire (Black), Projection Handrim Kit, Standard Mount | 1 |
| | 118225 | 22" MAXX Mag, Pneumatic Tire (Black), Projection Handrim Kit, Close Mount | 1 |
| | 119527 | 22" MAXX Mag Pneumatic Tire (Black) Ergonomic Standard Handrim Omit Thumb Grip Kit | 1 |
| | 119528 | 22" MAXX Mag Pneumatic Tire (Black) Ergonomic Standard Handrim Regular Thumb Grip Kit | 1 |
| | 119529 | 22" MAXX Mag Pneumatic Tire (Black) Ergonomic Standard Handrim Extra Thumb Grip Kit | 1 |
| | 119531 | 22" MAXX Mag Pneumatic Tire (Black) Ergonomic Lite Handrim Regular Thumb Grip Kit | 1 |
| | 119532 | 22" MAXX Mag Pneumatic Tire (Black) Ergonomic Lite Handrim Extra Thumb Grip Kit | 1 |
| | 118122 | 22" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount | 1 |
| | 118123 | 22" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Close Mount | 1 |
| | 118124 | 22" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip | 1 |
| | 118125 | 22" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip | 1 |
| | 118126 | 22" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Plastic Coated Handrim Kit, Standard Mount | 1 |
| | 118127 | 22" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Plastic Coated Handrim Kit, Close Mount | 1 |
| | 118128 | 22" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Projection Handrim Kit, Standard Mount | 1 |
| | 118129 | 22" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Projection Handrim Kit, Close Mount | 1 |
| | 119521 | 22" MAXX Mag Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Omit Thumb Grip Kit | 1 |
| | 119522 | 22" MAXX Mag Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Regular Thumb Grip Kit | 1 |

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|---|----------|
| | 119523 | 22" MAXX Mag Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Extra Thumb Grip Kit | 1 |
| | 119525 | 22" MAXX Mag Pneumatic w/ Airless Insert Tire (Black) Ergonomic Lite Handrim Regular Thumb Grip Kit | 1 |
| | 119526 | 22" MAXX Mag Pneumatic w/ Airless Insert Tire (Black) Ergonomic Lite Handrim Extra Thumb Grip Kit | 1 |
| | 117682 | 22" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount | 1 |
| | 117683 | 22" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Close Mount | 1 |
| | 117684 | 22" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip | 1 |
| | 117685 | 22" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip | 1 |
| | 117686 | 22" MAXX Mag, Full Poly Tire, Plastic Coated Handrim Kit, Standard Mount | 1 |
| | 117687 | 22" MAXX Mag, Full Poly Tire, Plastic Coated Handrim Kit, Close Mount | 1 |
| | 117688 | 22" MAXX Mag, Full Poly Tire, Projection Handrim Kit, Standard Mount | 1 |
| | 117689 | 22" MAXX Mag, Full Poly Tire, Projection Handrim Kit, Close Mount | 1 |
| | 119497 | 22" MAXX Mag Full Poly Tire Ergonomic Standard Handrim Omit Thumb Grip Kit | 1 |
| | 119498 | 22" MAXX Mag Full Poly Tire Ergonomic Standard Handrim Regular Thumb Grip Kit | 1 |
| | 119499 | 22" MAXX Mag Full Poly Tire Ergonomic Standard Handrim Extra Thumb Grip Kit | 1 |
| | 119501 | 22" MAXX Mag Full Poly Tire Ergonomic Lite Handrim Regular Thumb Grip Kit | 1 |
| | 119502 | 22" MAXX Mag Full Poly Tire Ergonomic Lite Handrim Extra Thumb Grip Kit | 1 |
| | 117720 | 22" MAXX Mag, Low Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount | 1 |
| | 117721 | 22" MAXX Mag, Low Poly Tire, Aluminum Anodized Handrim Kit, Close Mount | 1 |
| | 117722 | 22" MAXX Mag, Low Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip | 1 |
| | 117723 | 22" MAXX Mag, Low Poly Tire, Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip | 1 |
| | 117724 | 22" MAXX Mag, Low Poly Tire, Plastic Coated Handrim Kit, Standard Mount | 1 |
| | 117725 | 22" MAXX Mag, Low Poly Tire, Plastic Coated Handrim Kit, Close Mount | 1 |
| | 119503 | 22" MAXX Mag Low Poly Tire Ergonomic Standard Handrim Omit Thumb Grip Kit | 1 |
| | 119504 | 22" MAXX Mag Low Poly Tire Ergonomic Standard Handrim Regular Thumb Grip Kit | 1 |

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|--|----------|
| | 119505 | 22" MAXX Mag Low Poly Tire Ergonomic Standard Handrim Extra Thumb Grip Kit | 1 |
| | 119507 | 22" MAXX Mag Low Poly Tire Ergonomic Lite Handrim Regular Thumb Grip Kit | 1 |
| | 119508 | 22" MAXX Mag Low Poly Tire Ergonomic Lite Handrim Extra Thumb Grip Kit | 1 |
| | 117754 | 22" MAXX Mag, Sentinel Tire, Aluminum Anodized Handrim Kit, Standard Mount | 1 |
| | 117755 | 22" MAXX Mag, Sentinel Tire, Aluminum Anodized Handrim Kit, Close Mount | 1 |
| | 117756 | 22" MAXX Mag, Sentinel Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip | 1 |
| | 117757 | 22" MAXX Mag, Sentinel Tire, Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip | 1 |
| | 117758 | 22" MAXX Mag, Sentinel Tire, Plastic Coated Handrim Kit, Standard Mount | 1 |
| | 117759 | 22" MAXX Mag, Sentinel Tire, Plastic Coated Handrim Kit, Close Mount | 1 |
| | 119509 | 22" MAXX Mag Sentinel Tire Ergonomic Standard Handrim Omit Thumb Grip Kit | 1 |
| | 119510 | 22" MAXX Mag Sentinel Tire Ergonomic Standard Handrim Regular Thumb Grip Kit | 1 |
| | 119511 | 22" MAXX Mag Sentinel Tire Ergonomic Standard Handrim Extra Thumb Grip Kit | 1 |
| | 119513 | 22" MAXX Mag Sentinel Tire Ergonomic Lite Handrim Regular Thumb Grip Kit | 1 |
| | 119514 | 22" MAXX Mag Sentinel Tire Ergonomic Lite Handrim Extra Thumb Grip Kit | 1 |
| | 117800 | 22" MAXX Mag, Mountain Tire, Aluminum Anodized Handrim Kit, Standard Mount | 1 |
| | 117801 | 22" MAXX Mag, Mountain Tire, Aluminum Anodized Handrim Kit, Close Mount | 1 |
| | 117802 | 22" MAXX Mag, Mountain Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip | 1 |
| | 117803 | 22" MAXX Mag, Mountain Tire, Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip | 1 |
| | 117804 | 22" MAXX Mag, Mountain Tire, Plastic Coated Handrim Kit, Standard Mount | 1 |
| | 117805 | 22" MAXX Mag, Mountain Tire, Plastic Coated Handrim Kit, Close Mount | 1 |
| | 117806 | 22" MAXX Mag, Mountain Tire, Projection Handrim Kit, Standard Mount | 1 |
| | 117807 | 22" MAXX Mag, Mountain Tire, Projection Handrim Kit, Close Mount | 1 |
| | 119515 | 22" MAXX Mag Mountain Tire Ergonomic Standard Handrim Omit Thumb Grip Kit | 1 |
| | 119516 | 22" MAXX Mag Mountain Tire Ergonomic Standard Handrim Regular Thumb Grip Kit | 1 |
| | 119517 | 22" MAXX Mag Mountain Tire Ergonomic Standard Handrim Extra Thumb Grip Kit | 1 |
| | 119519 | 22" MAXX Mag Mountain Tire Ergonomic Lite Handrim Regular Thumb Grip Kit | 1 |

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|---|----------|
| | 119520 | 22" MAXX Mag Mountain Tire Ergonomic Lite Handrim Extra Thumb Grip Kit | 1 |