



## BATH/TOILET MODULAR SYSTEM- REHAB SHOWER/COMMUNE CHAIR

- Complete Packages for Easier Ordering Available
- Stainless Steel Frame
- Height Adjustable
- Can be used as a Raised Toilet Seat
- Lightweight
- 350 lb. (158kg) Weight Capacity
- Several Seat Options Available
- Seat usable in Four Different Directions

The 101's unique design accommodates perineal access from four directions. The cut-out portion of the seat can be placed over the toilet in the desired direction, and the legs can be adjusted in height while maintaining a stable platform for seating.

With the 12 different seating options and large seating surface the 101 provides excellent support, while the cut-out may offer some additional pressure relief. With these options, the customer can find the size, contour, and shape to fit a wide range of needs. It is lightweight, easy to move out of the way and no tools are required for disassembling the legs for travel or for height adjustment. Easy to use, versatile, sturdy and comfortable makes the 101 an outstanding bathroom aid.

101 OPTIONS

FRAMES



**PRT50000**  
18"Wx18"D  
White Stainless  
Steel Frame



Shown with  
**P61086 Seat**



Shown with  
**P60725 Seat**



Shown with  
**P61063 Seat**

LEGS



**PRT50001**  
Short Legs w/  
Crutch Tips  
(seat height range  
from floor to bottom  
frame-15.5"-18.5")



**PRT50002**  
Short Legs w/  
Suction Tips  
(seat height range  
from floor to bottom  
frame-15.5"-18.5")




**PRT50003**  
Long Legs w/  
Crutch Tips  
(seat height range  
from floor to bottom  
frame-18.5"-21.5")




**PRT50004**  
Long Legs w/  
Suction Tips  
(seat height range  
from floor to bottom  
frame-18.5"-21.5")

Seat Opening  
Measurements



Oval Opening  
Measurements



SEATS

\*Seat works in all 4 positions



**P60928\***  
16"Wx16"D  
Ensolite® 1" foam  
Front/Rear Open Seat  
A-4"xB-6.5" opening



**P60686\***  
16"Wx18"D  
Ensolite® 2" foam  
Front/Rear Open Seat  
A-4.25"xB-7" opening



**P60811\***  
16"Wx15.5"D  
Ensolite® 2" waterfall foam  
Front/Rear Open Seat  
A-3.5"xB-6" opening



**P60722\***  
18"Wx16"D  
Ensolite® 2" foam  
Front/Rear Open Seat  
A-4.5"xB-7.5" opening



**P60911\***  
18"Wx18"D  
Ensolite® 1" foam  
Front/Rear Open Seat  
A-4.5"xB-7.5" opening



**P60913\***  
18"Wx18"D  
Ensolite® 2" waterfall foam  
Front/Rear Open Seat  
A-5.25"xB-8" opening



**P60993\***  
18"Wx18"D  
Comfortuff 1.5" waterfall foam  
Front/Rear Open Seat  
A-5"xB-8" opening




**P60850\***  
18"Wx18"D  
Ensolite® 1" waterfall foam  
Front/Rear Open Seat  
A-5"xB-8" opening



**P61063**  
18"Wx18"D  
Ensolite® 2" foam  
Front Open Contoured Seat  
A-4.25"xB-7.5" opening



**P61058**  
18"Wx18"D  
Ensolite® 2" waterfall foam  
Front Open Contoured Seat  
A-5.5"xB-8.25" opening



**P60725**  
18"Wx16"D  
Ensolite® 2" foam  
Oval Seat  
A-7.75"xB-9" opening



**P61086**  
18"Wx18"D  
Ensolite® 2" foam  
Solid Seat

full option details at  
[activeaid.com](http://activeaid.com)

SPECIFICATIONS

Weight Capacity

350 lbs. (158kg)

Seat Width

16" (40cm)

Seat Depth

16" (40cm)

Seat Height Range From Floor to Top of Seat (transfer height)

18.25"-21.25" (46-53cm) short legs, crutch tip  
21"-24" (53-60cm) long legs, crutch tip

Under Seat Height Clearance

15.5"-18.5" (39-46cm) short legs, crutch tip  
18.5"-21.5" (46-54cm) long legs, crutch tip

Under Seat Width Clearance

16" (40cm)

Overall Width

min.- short leg, crutch tip  
min.- short leg, suction cup  
max.- long leg, crutch tip  
max.- long leg, suction cup

21.75" (55cm)  
23.75" (60cm)  
23" (58cm)  
25" (63cm)

Back of Unit to Start of Commode Opening

6.75"-7" (17cm) short legs, crutch tip  
7"-7.25" (17cm) long legs, crutch tip

Seat Used to Generate Specifications

#60928 16"Wx16"Dx1.5"H