


# (Discontinued) Catalyst Depth Adjustable Backrest Stroller Handle Extension

---



*This backrest style and stroller handle extension discontinued February 8th, 2021 Not compatible with Drum Brake. A backrest retrofit to the Catalyst Stroller Handle Back Post with Removable Stroller Handle Extension (Depth Adjustable) is required if retrofitting Drum Brake.*

<b>Item Number</b>	<b>Part Number</b>	<b>Description</b>	<b>Quantity</b>
5-9	502397	Stroller Handle Depth Adj Back Post Assembly, Short CAT	1
5-9	502398	Stroller Handle Depth Adj Back Post Assembly, Tall CAT	1
4	100533	Detent Button	1
6	000034	Sleeve 3/4" x 1 1/2"	1
7	000065	M6 x 23.68 Threaded Barrel	1
8	100667	M6 x 12 BTN SHCS BLKZC w/Patch	2