





RIDICULOUSLY LIGHT: 32 LBS

One of the world's lightest carbon fiber folding power wheelchairs.

MODERN DESIGN

Sleek, ergonomic design with premium materials – the epitome of comfort and style. Available in 5 accent colors: Blue, Orange, Purple, Red, and Yellow.

GO UP TO 15 MILES

Bring an (optional) second battery along for the ride to double your range to up to 15 miles!

EFFORTLESS FOLDING

Flip the clip & push the backrest down. You don't even need to remove the battery.

COMPACT FOLDED SIZE

Folds to a compact size of 30.2" length x 22.4" width x 14.6" height, to fit small car trunks, buses, and closets.

SMOOTH RIDE

Its carbon fiber frame can naturally flex over uneven surfaces, providing unique shock-absorbing characteristics to deliver one of the smoothest, most efficient rides around.

LIGHT BUT INCREDIBLY STRONG

Can support a user weight of up to 300 lbs.

NARROW 22.4" WIDTH

Its slim width easily fits through narrow doorways and navigates busy rooms.

RELIABLE MOBILITY

250W high-efficiency brushless motors provide the power and reliability you need.





TECHNICAL DATA	
36.5" Long x 22.4" Wide	
17" Wide X 16.5" Deep	
18.3"	
7.5 miles (one battery), 15 miles (with optional additional battery)	
3.7 MPH	
10Ah lithium-ion	
Rear	
8" Drive Wheel; 6.5" Casters	
300 lbs.	
Chair (without batteries): from 32 lbs., Battery (single): 3 lbs.	

Refer to order form as specifications may vary based on configuration.

















