

INVACARE® MATTRESSES



Now you can offer your customers the comfort, protection and durability of superior **Homecare Mattresses** by Invacare, along with **Invacare® Mattress Covers**. Choose from innerspring or foam mattresses. Mattress covers will keep your Invacare® mattress clean and fresh, and protect it from stains, odors, dust and incontinence. Covers are made of non-allergenic waterproof vinyl and are available in zippered and contoured styles.

Invacare beds, mattresses and rails, excluding therapeutic support surfaces and bariatric beds, when used as a system and sold new after August 1, 2007, meet the voluntary FDA Guidelines for bed rail entrapment. Visit www.invacare.com/bedsafety for additional safety and entrapment information.

Invacare® Innerspring Mattress

Model # 5185 , 5185XL

High quality .312 innerspring design

Durable 11-ounce vinyl cover is anti-bacterial, anti-static, acid-resistant and waterproof for easy use and care

5185XL (not shown) is available for use with extender kit

Invacare® Economy Foam Mattress

Model # 5180 , 5184

High quality 1.5 lb. density foam provides comfort and durability in a one-piece construction

Invacare® Bariatric Foam Mattress

Model # BARMATT(39/42/48/EXP)

Constructed of latex-free materials

Durable 11-ounce vinyl cover is anti-bacterial and waterproof for easy use and care

Solid inner core

Mattress is reversible and flippable

Constructed of high resilience foam for extended life and service for the homecare provider

Easy to clean

HCPCS Code: E0271

HCPCS Code: E01272



5185



BARMATT39

- 1. Invacare® MC0195-1 zippered mattress cover comes individually wrapped; 12 per case
- 2. Invacare® MC0065-1 contoured mattress cover comes individually wrapped; 12 per case

Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2012
Invacare Corporation.
All rights reserved.
Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.
Form No. 11-219
111225



| Model # | Dimensions LxWxD | Gauge | Regulation Compliance | Limited Warranty |
|--------------------|------------------|-------|-----------------------|------------------|
| MC0165-1, MC0195-1 | 80"x36"x6" | 4" | N/A | 30 Days |
| 5180 | 80"x36"x5.5" | N/A | 16CFR1632/33 | 1 Year |
| 5184 | 84"x36"x5.5" | N/A | 16CFR1632/33 | 1 Year |
| 5185 | 80"x16"x6" | N/A | 16CFR1632/33 | 1 Year |
| 5185XL | 84"x36"x6" | N/A | 16CFR1632/33 | 1 Year |
| BARMATT39 | 80"x39"x6.5" | N/A | 16CFR1632/33 | 1 Year |
| BARMATT42 | 80"x42"x6.5" | N/A | 16CFR1632/33 | 1 Year |
| BARMATT48 | 80"x48"x6.5" | N/A | 16CFR1632/33 | 1 Year |
| BARMATTEXP | 80"x39"-48"x6.5" | N/A | 16CFR1632/33 | 1 Year |

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

