


Ki Mobility Belt Mounting Clamp



Item Number	Part Number	Description	Quantity	Retail
1-5	503044	Ki Mobility Belt Mounting Clamp Kit, Pair, 1 3/8" Tubing	1	\$35.00
1	003041	Seat Belt Clamp A Rigid	2	\$9.19
2	003042	Seat Belt Clamp B Rigid	2	\$9.19
3	101781	M6 x 50mm Partial Thread BHCS Blk ZC	2	\$1.58
4	100746	M6 Flat Washer Blk Zn	2	\$1.58
5	100658	M6 Nylock Nut Blk Zn	2	\$1.05