



Raz-SP600

RazBari Self Propel

MSRP – USD – USA
Effective Feb. 1, 2021

Date: _____

Dealer / Location: _____ T: _____ F: _____

PO# / Client: _____ Contact: _____ Email: _____

Raz-SP600 Heavy Duty Mobile Shower Commode Chair with stainless steel frame, fore / aft-adjustable 24" rear wheels, flip-up arm supports, heavy duty foot supports, front dual-locking 5" casters and commode pan. **User weight capacity 600 lbs.** User weight over 450 lbs. requires 6" casters (Z66SP). See our [Evaluation and Justification Form](#).

1. Chair Codes by size – Must Choose One

Chair Code	Seat Frame Width	MSRP
Z26022	22"W Frame	\$2,895
Z26024	24"W Frame	\$3,120
Z26026	26"W Frame	\$3,545
Z26028	28"W Frame	\$3,770
Z26030	30"W Frame – Requires custom seat	\$3,995

2. Seats – Must Choose One

Molded and Visco Foam Interface Seats			
Z4M23	Molded 23"W x 18"D – Front Access	Fits 22" / 24"W ¹ frames only	N/C
Z4M26	Molded 26"W x 18"D – Front Access	Fits 26" / 28"W ¹ frames only	N/C
Z4V23	Visco Foam Interface 23"W x 18"D – Front Access	Fits 22" / 24"W ¹ frames only	\$595
Z4V26	Visco Foam Interface 26"W x 18"D – Front Access	Fits 26" / 28"W ¹ frames only	\$695
Z4SSKASHD	Seat Support Kit for Raz-AP600 – for seats that are narrower than seat frames		\$195

¹ Required for seats that are narrower than seat frames

Custom Seats

For custom seats, please specify seat width and depth. Indicate if the access opening is to be in the Front, Left, Right or Rear. Select Bridged if an access opening is not required. See [Custom Seat Assessment Form](#) – Resources tab at razdesigninc.com

Z4CSHD	Heavy Duty Custom Seat (enter specs below) Select if both width and depth are ≤ 28"	\$1,095
Z4CSXHD	Heavy Duty Custom Seat for widths 28" and 30" (enter specs below) Select if width and/or depth is ≥ 28"	\$1,295

Seat Frame Width: ___ 22" / 24" / 26" / 28" / 30"
 Custom Seat Width: ___ 22" / 24" / 26" / 28" / 30" / 32" (cannot be more than 2" wider than frame width)
 Seat Depth: ___ 16" / 17" / 18" / 19" / 20" / 21" / 22" / 23" / 24" [≥ 21" requires custom longer frame]
 Access: ___ FR – Front / LT – Left / RT – Right / RR – Rear / BR – Bridged
 Aperture Width: ___ 6.5" recommended

- Unless otherwise requested, the width of the aperture (inside width dimension of the commode opening) will be 6.5" for custom seats.
- Custom seat widths can be up to 2" wider than the frame width to provide additional thigh support without increasing overall width of the chair.
- For Custom Seats, the Ischial Pelvic Alignment System (IPAS) allows +/- 1" fore / aft adjustment from specified custom seat depth.
- Access for Rear-Open seats is limited due to rear seat support tube.
- Special modifications on Custom Seats include: leg-length discrepancy, soft spots, pressure-relief areas, contouring, pelvic obliquity accommodation, etc.

3. Back Support Options (Symphony Back replaces standard Tension-Adjustable Vinyl Upholstery – recommended for infection-control environments) All Symphony Back panels are 13”H.

Z5S22	Symphony Back – Recessed Solid Back Support – 22”W	\$595
Z5S24	Symphony Back – Recessed Solid Back Support – 24”W	\$695
Z5S26	Symphony Back – Recessed Solid Back Support – 26”W	\$795
Z5S28	Symphony Back – Recessed Solid Back Support – 28”W	\$895
Z5S30	Symphony Back – Recessed Solid Back Support – 30”W	\$995

4. Rear Wheel, Caster and Anti-Tipper Options

Z614RHD	24” HD Wheels w/ Ergonomic Handrims	Standard	Z66SPR	Rear Anti-Tippers for regular frame (pr)	\$110
Z6WLE	Wheel Lock with Extensions (pr) TEMPORARILY OUT OF STOCK	\$95	Z66AS5	Front Anti-Tippers for 5” casters (pr)	\$245
Z624VHR	Vinyl-Coated Handrims – pr	\$265	Z66AS6	Front Anti-Tippers for 6” casters (pr)	\$245
Z624PHR	Vertical Projection Handrims – pr	\$495	Z66SP	6” Caster Upgrade (2x) – required for user weight > 450 lbs	\$175
Z674HD	User-Actuated Front Directional Caster Locks w/ extension handles (pr)	\$345			

5. Arm Support Accessories

Z710HD	Arm Support Lock Kit – HD (pr)	\$95	Z7PAML	Pivoting Arm Mount – Left (for Z7FL / Z7TR)	\$165
Z7FLL	Long Molded Flat Armpad ¹ – Left (14”L x 3.25”W)	\$95	Z7PAMR	Pivoting Arm Mount – Right (for Z7FL / Z7TR)	\$165
Z7FLR	Long Molded Flat Armpad ¹ – Right (14”L x 3.25”W)	\$95			
Z7TRL	Molded Arm Trough ¹ – Left (18.5”L x 4.5”W)	\$110			
Z7TRR	Molded Arm Trough ¹ – Right (18.5”L x 4.5”W)	\$110			
¹ May require Z7PAM for arm support to flip up					

6. Foot and Leg Supports – Standard Foot Supports (Z709N) are included with the chair. The selection of any Foot / Leg Support Option below automatically replaces the relevant Standard Foot Support. e.g. By selecting Elevating Leg Support – Left (Z7ELRL), the chair will be supplied with one Foot Support (Z709N) on the right side and a Z7ELRL on the left side.

Z709N	Adjustable Footrests HD (16”-20.5”) – pr	Std	Z7ELRL	Elevating Legrest ¹ – Left (incl Z7ADM)	\$475
Z7FET2	Footrest Extensions (2” longer) – pr	\$50	Z7ELRR	Elevating Legrest ¹ – Right (incl Z7ADM)	\$475
Z7ADM	Angle/Depth Adj Footplates – pr	\$190	Z7ELROFL	Forward Offset Elevating Legrest ² – Left (incl Z7ADM)	\$525
Z7FET2A	Footrest Extension Tubes for Z7ADM (pr) – 2” longer	\$50	Z7ELROFR	Forward Offset Elevating Legrest ² – Right (incl Z7ADM)	\$525
Z7RLSL	Universal Residual Limb Support – Left	\$395			
Z7RLSR	Universal Residual Limb Support – Right	\$395			
¹ Z7ELRL/R compatible with adjusted seat depths ≤ 20”D					
² Z7ELROFL/R compatible with adjusted seat depths 20”-22”D					

7. Foot and Leg Support Accessories

Z7FSLK	Foot Support Lock (pr)	\$95	Z708ADJ	Adjustable Heel Loops (pr)	\$60
Z7CS__	Calf Strap, Snap-On – Bodypoint Aeromesh * Enter seat frame width to complete code (22"/24"W)	\$95	Z7NCS	Neoprene Covers for Standard Footplates (pr)	\$90
Z7CS__*	Calf Strap, Snap-On – Bodypoint Aeromesh * Enter seat frame width to complete code (26"/28"/ 30"W)	\$195	Z7NCA	Neoprene Covers for Adj Angle/Depth Molded Footplates (pr)	\$90
Z7CSIC__*	Calf Strap, Snap-On – Infection Control *Enter seat frame width to complete code (22" / 24" / 26" / 28" / 30")	\$140	Z7OFLM	Lateral Offset Foot Support Receiver – moves foot support 2" laterally or medially (ea) Quantity (1 or 2) __	\$160
Z712	H-Strap – for 22"/24"W frames – hook and loop attachment	\$110			

8. Positioning Pads

Z7HSASS	Headrest Pad – Small (10" x 4")	\$495	Z7LGAS2L	Leg Adductor – w/ Large Pad ¹ – Left	\$320
Z7HSASL	Headrest Pad – Large (12" x 6")	\$545	Z7LGAS2R	Leg Adductor – w/ Large Pad ¹ – Right	\$320
Z7LTS2HDL	Lateral Thoracic Support w/ Large Pad ^{1,2} – Left	\$320	Z742L	Lateral/Adductor Extension Plate (5.5" x 3.5") – Left	\$25
Z7LTS2HDR	Lateral Thoracic Support w/ Large Pad ^{1,2} – Right	\$320	Z742R	Lateral/Adductor Extension Plate (5.5" x 3.5") – Right	\$25

9. Pelvic and Chest Belts – A CHEST BELT MUST NEVER BE USED WITHOUT A PELVIC BELT!

Pelvic Belts				Chest Belts			
ZPBLHDT	Pelvic Belt – Large (frame widths 22", 2-piece w/ side-release buckle)	\$65		ZCBL	Chest Belt ² – Large (back widths 22", 1-piece w/ hook and loop closure)	\$45	
ZPBXLHDT	Pelvic Belt – XLarge (seat frame widths 22"-26", 2-piece w/ side-release buckle)	\$75		ZCBXL	Chest Belt – XLarge ² (back frame widths 22"-26", 1-piece w/ hook and loop closure)	\$55	
ZPBXXLHDT	Pelvic Belt – XXL Large (seat frame widths 26"-30", 2-piece w/ side-release buckle)	\$85		ZCBXXL	Chest Belt – XXL Large ² (back frame widths 26"-30", 1-piece w/ hook and loop closure)	\$65	
ZPBICLHDT	Pelvic Belt ¹ – Infection Control, Large 2-piece w/ frame clamps, mounting plates, side-release buckle (frame widths 22")	\$115		ZCBICL	Chest Belt ^{1,2} – Infection Control, Large 1-piece w/ side-release buckle (back widths 22")	\$85	
ZPBICXLHDT	Pelvic Belt ¹ – Infection Control, XLarge 2-piece w/ frame clamps, mounting plates, side-release buckle (frame widths 22"-26")	\$125		ZCBICXL	Chest Belt ^{1,2} – Infection Control, XLarge 1-piece w/ side-release buckle (back frame widths 22"-26")	\$95	
ZPBICXXLHDT	Pelvic Belt ¹ – Infection Control, XXL Large 2-piece w/ frame clamps, mounting plates, side-release buckle (frame widths 26"-30")	\$135		ZCBICXXL	Chest Belt ^{1,2} – Infection Control, XXL Large 1-piece w/ side-release buckle (back frame widths 26"-30")	\$105	
ZPBBPLHDT	Pelvic Belt – Bodypoint Aeromesh, Large	\$195		ZCBBP_*	* Chest Belt – Bodypoint Aeromesh * enter size to complete code (L / XL)	\$205	
ZPBBPXLHDT	Pelvic Belt – Bodypoint Aeromesh, XLarge	\$165					
Recommended sizes for Bodypoint belts	Chest belt	ZCBBPL	22"W back frame	ZCBBPXL	22"-26"W back frame	ZCBBPXXL	26"-30"W back frame
	Pelvic belt	ZPBBPLHDT	22"W base frame	ZPBBPXLHDT	22"-26"W base frame	ZPBBPXXLHDT	26"-30"W base frame
¹ Non-absorbent belts – recommended for infection-control environments.							
² Tension-Adjustable Upholstery includes loops to retain a chest belt (ZCB). When the chair is supplied with a Symphony Back, Chest Belt loops are provided.							

10. Other

ZPSK1	Commode Pan Spacer Kit ¹ – 1" lower	\$200	ZIHW	Whizard – winged urine deflector ²	\$67
ZPSK2	Commode Pan Spacer Kit ¹ – 2" lower	\$215	ZBUDBH_*	Bard Urine Drainage Bag Holder * Enter size to complete code (22"/24"/26"/28"/30")	\$155
ZPSK3	Commode Pan Spacer Kit ¹ – 3" lower	\$230			
ZPSK4	Commode Pan Spacer Kit ¹ – 4" lower	\$245			
¹ Commode Pan Spacer Kits provide extra clearance for perianal care with pan in place. Reduces toilet bowl clearance height.			² Has patent-pending lateral wings that are positioned under the medial aspect of the thighs to keep the Whizard in place.		

11. Credits and Omissions

Z905	Omit Commode Pan and Mounting Brackets	-\$20	Z907HD	Omit HD Foot Supports – if no other Foot / Leg Support Option is selected	-\$45
------	--	-------	--------	---	-------

Some assembly of Raz products is required.

12. Notes:

Specifications

Raz-SP600 (Self Propel)

Weight Capacity	600 lbs	
Overall Width	Frame width plus 7.5"	
Overall Length (with foot supports in middle position)	44" (with rear wheel in middle position)	
Seat Depth	18"-20"	
Clearance Height (under seat)	with 5" casters	with 6" casters
- with pan brackets	14.35"-18.35"	15.35"-19.35"
- without pan brackets	15.12"-19.12"	16.12"-20.12"
Arm Support Height - above molded seat	10"	
Seat-to-Floor Height (with molded seat)	with 5" casters	with 6" casters
	20" and 21"	20", 21", 24"
Back Height (Molded Seat to top of upholstery)	20"	
Product Weight (with 22"W seat frame)	68 lbs	
Frame and Hardware Material	Stainless Steel (lifetime warranty against corrosion)	