

X'Cape Growth Seat Sling Order Matrix													
Frame Depth	Available Seat Sling Depths												
	8	9	10	11	12	13	14	15	16	17	18	19	20
8	8	9	10	11									
9	8	9	10	11	12								
10	8	9	10	11	12	13							
11	8	9	10	11	12	13	14						
12		9	10	11	12	13	14	15					
13			10	11	12	13	14	15	16				
14				11	12	13	14	15	16	17			
15					12	13	14	15	16	17	18		
16						13	14	15	16	17	18	19	
17							14	15	16	17	18	19	20
18								15	16	17	18	19	20

Green seat sling depths require a growth seat sling.

To achieve this, have the chairs frame depth match the seat sling depth (may need new seat sling) and order a growth seat sling. This will allow enough growth to accommodate for the new seat sling depth.

White seat sling depths can ether order a new seat sling to match new depth or keep current seat sling and order a growth seat sling.

To change width, a new seat sling will need to be ordered.