


Discontinued 2 Catalyst Tension Adjustable Back Upholstery



Tension Straps (7a-7d) are used with the combination of Fold Down Push Handle and Tension Adjustable Back Upholstery. Tension straps (5a-5d and 7a-7d) are based on chair width. Backrest cover kits with two straps and Tension Adjustable assemblies include Installation Instructions.

Item Number	Part Number	Description	Quantity
1al, 3a, 4a, 5d, 5m	504159	Back Upholstery Assembly Tension Adjustable 21W 14.5-15.5H, Standard	1
1al, 3c, 4c, 5d, 5m	504597	Back Upholstery Assembly Tension Adjustable 21W 14.5-15.5H, Standard, Heavy Duty	1
1ao, 3a, 4a, 5d, 5n	504148	Back Upholstery Assembly Tension Adjustable 22W 14.5-15.5H, Standard	1
1ao, 3c, 4c, 5d, 5n	504632	Back Upholstery Assembly Tension Adjustable 22W 14.5-15.5H, Standard, Heavy Duty	1
1a, 5f	504254	Backrest Cover with Two Straps 14W Cat Std, 9-12H <i>Contains 002497 cover and two 003485 straps</i>	1
1b, 5f	504255	Backrest Cover with Two Straps 14W Cat Std, 11-14H <i>Contains 002498 cover and two 003485 straps</i>	1
1c, 5f	504256	Backrest Cover with Two Straps 14W Cat Std, 13-16H <i>Contains 002499 cover and two 003485 straps</i>	1
1d, 5f	504257	Backrest Cover with Two Straps 14W Cat Std, 15-18H <i>Contains 002500 cover and two 003485 straps</i>	1
1e, 5f	504258	Backrest Cover with Two Straps 14W Cat Std, 17-20H <i>Contains 003043 cover and two 003485 straps</i>	1
1c, 5e	504227	Backrest Cover with Two Straps 14W Cat Angle Adjustable, 13-16H <i>Contains 002499 cover and two 003484 straps</i>	1
1d, 5e	504228	Backrest Cover with Two Straps 14W Cat Angle Adjustable, 15-18H <i>Contains 002500 cover and two 003484 straps</i>	1
1e, 5e	504229	Backrest Cover with Two Straps 14W Cat Angle Adjustable, 17-20H <i>Contains 003045 cover and two 003484 straps</i>	1
1f, 5g	504259	Backrest Cover with Two Straps 15W Cat Std, 9-12H <i>Contains 002501 cover and two 003486 straps</i>	1
1g, 5g	504260	Backrest Cover with Two Straps 15W Cat Std, 11-14H <i>Contains 002502 cover and two 003486 straps</i>	1
1h, 5g	504261	Backrest Cover with Two Straps 15W Cat Std, 13-16H <i>Contains 002503 cover and two 003486 straps</i>	1
1i, 5g	504262	Backrest Cover with Two Straps 15W Cat Std, 15-18H <i>Contains 002504 cover and two 003486 straps</i>	1
1j, 5g	504263	Backrest Cover with Two Straps 15W Cat Std, 17-20H <i>Contains 003044 cover and two 003486 straps</i>	1
1h, 5f	504230	Backrest Cover with Two Straps 15W Cat Angle Adjustable, 13-16H <i>Contains 003485 cover and two 002503 straps</i>	1

Item Number	Part Number	Description	Quantity
1i, 5f	504231	Backrest Cover with Two Straps 15W Cat Angle Adjustable, 15-18H <i>Contains 003485 cover and two 002504 straps</i>	1
1j, 5f	504232	Backrest Cover with Two Straps 15W Cat Angle Adjustable, 17-20H <i>Contains 003044 cover and two 003485 straps</i>	1
1k, 5h	504264	Backrest Cover with Two Straps 16W Cat Std, 9-12H <i>Contains 002505 cover and two 003487 straps</i>	1
1l, 5h	504265	Backrest Cover with Two Straps 16W Cat Std, 11-14H <i>Contains 002506 cover and two 003487 straps</i>	1
1m, 5h	504266	Backrest Cover with Two Straps 16W Cat Std, 13-16H <i>Contains 002507 cover and two 003487 straps</i>	1
1n, 5h	504267	Backrest Cover with Two Straps 16W Cat Std, 15-18H <i>Contains 002508 cover and two 003487 straps</i>	1
1o, 5h	504268	Backrest Cover with Two Straps 16W Cat Std, 17-20H <i>Contains 003045 cover and two 003487 straps</i>	1
1m, 5g	504233	Backrest Cover with Two Straps 16W Cat Angle Adjustable, 13-16H <i>Contains 002507 cover and two 003486 straps</i>	1
1n, 5g	504234	Backrest Cover with Two Straps 16W Cat Angle Adjustable, 15-18H <i>Contains 002508 cover and two 003486 straps</i>	1
1o, 5g	504235	Backrest Cover with Two Straps 16W Cat Angle Adjustable, 17-20H <i>Contains 003045 cover and two 003486 straps</i>	1
1p, 5i	504270	Backrest Cover with Two Straps 17W Cat Std, 9-12H <i>Contains 002509 cover and two 002532 straps</i>	1
1q, 5i	504271	Backrest Cover with Two Straps 17W Cat Std, 11-14H <i>Contains 002510 cover and two 002532 straps</i>	1
1r, 5i	504272	Backrest Cover with Two Straps 17W Cat Std, 13-16H <i>Contains 002511 cover and two 002532 straps</i>	1
1s, 5i	504273	Backrest Cover with Two Straps 17W Cat Std, 15-18H <i>Contains 002512 cover and two 002532 straps</i>	1
1t, 5i	504274	Backrest Cover with Two Straps 17W Cat Std, 17-20H <i>Contains 003046 cover and two 002532 straps</i>	1
1r, 5h	504236	Backrest Cover with Two Straps 17W Cat Angle Adjustable, 13-16H <i>Contains 002511 cover and two 003487 straps</i>	1
1s, 5h	504237	Backrest Cover with Two Straps 17W Cat Angle Adjustable, 15-18H <i>Contains 003487 cover and two 002512 straps</i>	1
1t, 5h	504238	Backrest Cover with Two Straps 17W Cat Angle Adjustable, 17-20H <i>Contains 003046 cover and two 003487 straps</i>	1

Item Number	Part Number	Description	Quantity
1u, 5j	504276	Backrest Cover with Two Straps 18W Cat Std, 9-12H <i>Contains 003048 cover and two 002533 straps</i>	1
1v, 5j	504277	Backrest Cover with Two Straps 18W Cat Std, 11-14H <i>Contains 003049 cover and two 002533 straps</i>	1
1w, 5j	504278	Backrest Cover with Two Straps 18W Cat Std, 13-16H <i>Contains 003050 cover and two 002533 straps</i>	1
1x, 5j	504279	Backrest Cover with Two Straps 18W Cat Std, 15-18H <i>Contains 003051 cover and two 002533 straps</i>	1
1y, 5j	504280	Backrest Cover with Two Straps 18W Cat Std, 17-20H <i>Contains 003052 cover and two 002533 straps</i>	1
1w, 5i	504239	Backrest Cover with Two Straps 18W Cat Angle Adjustable, 13-16H <i>Contains 003050 cover and two 002532 straps</i>	1
1x, 5i	504240	Backrest Cover with Two Straps 18W Cat Angle Adjustable, 15-18H <i>Contains 003051 cover and two 002532 straps</i>	1
1y, 5i	504241	Backrest Cover with Two Straps 18W Cat Angle Adjustable, 17-20H <i>Contains 003052 cover and two 002532 straps</i>	1
1z, 5k	504282	Backrest Cover with Two Straps 19W Cat Std, 9-12H <i>Contains 003054 cover and two 002534 straps</i>	1
1aa, 5k	504283	Backrest Cover with Two Straps 19W Cat Std, 11-14H <i>Contains 003055 cover and two 002534 straps</i>	1
1ab, 5k	504284	Backrest Cover with Two Straps 19W Cat Std, 13-16H <i>Contains 003056 cover and two 002534 straps</i>	1
1ac, 5k	504285	Backrest Cover with Two Straps 19W Cat Std, 15-18H <i>Contains 003057 cover and two 002534 straps</i>	1
1ad, 5k	504286	Backrest Cover with Two Straps 19W Cat Std, 17-20H <i>Contains 003058 cover and two 002534 straps</i>	1
1ab, 5j	504242	Backrest Cover with Two Straps 19W Cat Angle Adjustable, 13-16H <i>Contains 003056 cover and two 002533 straps</i>	1
1ac, 5j	504243	Backrest Cover with Two Straps 19W Cat Angle Adjustable, 15-18H <i>Contains 003057 cover and two 002533 straps</i>	1
1ad, 5j	504244	Backrest Cover with Two Straps 19W Cat Angle Adjustable, 17-20H <i>Contains 003058 cover and two 002533 straps</i>	1
1ae, 5l	504288	Backrest Cover with Two Straps 20W Cat Std, 9-12H <i>Contains 003060 cover and two 002535 straps</i>	1
1af, 5l	504289	Backrest Cover with Two Straps 20W Cat Std, 11-14H <i>Contains 003061 cover and two 002535 straps</i>	1

Item Number	Part Number	Description	Quantity
1ag, 5l	504290	Backrest Cover with Two Straps 20W Cat Std, 13-16H <i>Contains 003062 cover and two 002535 straps</i>	1
1ah, 5l	504291	Backrest Cover with Two Straps 20W Cat Std, 15-18H <i>Contains 003063 cover and two 002535 straps</i>	1
1ai, 5l	504292	Backrest Cover with Two Straps 20W Cat Std, 17-20H <i>Contains 003064 cover and two 002535 straps</i>	1
1ag, 5k	504245	Backrest Cover with Two Straps 20W Cat Angle Adjustable, 13-16H <i>Contains 003062 cover and two 002534 straps</i>	1
1ah, 5k	504246	Backrest Cover with Two Straps 20W Cat Angle Adjustable, 15-18H <i>Contains 003063 cover and two 002534 straps</i>	1
1ai, 5k	504247	Backrest Cover with Two Straps 20W Cat Angle Adjustable, 17-20H <i>Contains 003064 cover and two 002534 straps</i>	1
1aj, 5m	504294	Backrest Cover with Two Straps 21W Cat Std, 9-12H <i>Contains 003066 cover and two 002536 straps</i>	1
1ak , 5m	504295	Backrest Cover with Two Straps 21W Cat Std, 11-14H <i>Contains 003067 cover and two 002536 straps</i>	1
1al, 5m	504296	Backrest Cover with Two Straps 21W Cat Std, 13-16H <i>Contains 003068 cover and two 002536 straps</i>	1
1am, 5m	504297	Backrest Cover with Two Straps 21W Cat Std, 15-18H <i>Contains 003069 cover and two 002536 straps</i>	1
1an, 5m	504298	Backrest Cover with Two Straps 21W Cat Std, 17-20H <i>Contains 003070 cover and two 002536 straps</i>	1
1al, 5l	504248	Backrest Cover with Two Straps 21W Cat Angle Adjustable, 13-16H <i>Contains 003068 cover and two 002535 straps</i>	1
1am, 5l	504249	Backrest Cover with Two Straps 21W Cat Angle Adjustable, 15-18H <i>Contains 003069 cover and two 002535 straps</i>	1
1an, 5l	504250	Backrest Cover with Two Straps 21W Cat Angle Adjustable, 17-20H <i>Contains 003070 cover and two 002535 straps</i>	1
1ao, 5n	504300	Backrest Cover with Two Straps 22W Cat Std, 13-16H <i>Contains 003072 cover and two 002537 straps</i>	1
1ap, 5n	504301	Backrest Cover with Two Straps 22W Cat Std, 15-18H <i>Contains 003073 cover and two 002537 straps</i>	1
1aq, 5n	504302	Backrest Cover with Two Straps 22W Cat Std, 17-20H <i>Contains 003074 cover and two 002537 straps</i>	1
1ao, 5m	504251	Backrest Cover with Two Straps 22W Cat Angle Adjustable, 13-16H <i>Contains 003072 cover and two 002536 straps</i>	1

Item Number	Part Number	Description	Quantity
1ap, 5m	504252	Backrest Cover with Two Straps 22W Cat Angle Adjustable, 15-18H <i>Contains 003073 cover and two 002536 straps</i>	1
1q, 5m	504253	Backrest Cover with Two Straps 22W Cat Angle Adjustable, 17-20H <i>Contains 003074 cover and two 002536 straps</i>	1
5e	003484	Back Upholstery Strap 9.3" <i>Used on 14" chair width with Angle Adjustable backposts</i>	2
5f	003485	Back Upholstery Strap 11.3" <i>Used on 14" chair width with 8° Bend, Straight, or Straight with Push Handle backposts and 15" chair width with Angle Adjustable backposts</i>	2
5g	003486	Back Upholstery Strap 13.3" <i>Used on 15" chair width with 8° Bend, Straight, or Straight with Push Handle backposts and 16" chair width with Angle Adjustable backposts</i>	2
5h	003487	Back Upholstery Strap 15.3" <i>Used on 16" chair width with 8° Bend, Straight, or Straight with Push Handle backposts and 17" chair width with Angle Adjustable backposts</i>	2
5i	002532	Back Upholstery Strap 17.3" <i>Used on 17" chair width with 8° Bend, Straight, or Straight with Push Handle backposts and 18" chair width with Angle Adjustable backposts</i>	2
5j	002533	Back Upholstery Strap 19.3" <i>Used on 18" chair width with 8° Bend, Straight, or Straight with Push Handle backposts and 19" chair width with Angle Adjustable backposts</i>	2
5k	002534	Back Upholstery Strap 21.3" <i>Used on 19" chair width with 8° Bend, Straight, or Straight with Push Handle backposts and 20" chair width with Angle Adjustable backposts</i>	2
5l	002535	Back Upholstery Strap 23.3" <i>Used on 20" chair width with 8° Bend, Straight, or Straight with Push Handle backposts and 21" chair width with Angle Adjustable backposts</i>	2
5m	002536	Back Upholstery Strap 25.3" <i>Used on 21" chair width with 8° Bend, Straight, or Straight with Push Handle backposts and 22" chair width with Angle Adjustable backposts</i>	2
5n	002537	Back Upholstery Strap 27.3" <i>Used on 22" chair width with 8° Bend, Straight, or Straight with Push Handle backposts</i>	2
6	100734	M5 x 16 PHMS ZC	2
7a	001750	Tension Strap Folding Handle 14-15 Wide	1
7b	001751	Tension Strap Folding Handle 16-17 Wide	1

Item Number	Part Number	Description	Quantity
7c	001752	Tension Strap Folding Handle 18-19 Wide	1
7d	001753	Tension Strap Folding Handle 20-22 Wide	1