



# Focus / Flip / Arc Heel Loops



Heel Loop							
Part No.	Used on Chair Widths						Strap Length
	Standard Hanger			Offset Hanger			
	Angle Adjustable / Multi-Angle	Composite	Composite Angle Adjustable	Angle Adjustable / Multi-Angle	Composite	Composite Angle Adjustable	
100591	N/A	14" - 15"	N/A	12" - 13"	12" - 13"	N/A	9.4'
100592	14" - 15"	16" - 17"	15" - 16"	14" - 15"	14" - 15"	13" - 14"	10.4'
100593	16" - 24"	18" - 24"	17" - 24"	16" - 24"	16" - 24"	15" - 24"	11.4'

Adjustable Heel Loop							
Part No.	Used on Chair Widths						Strap Length
	Standard Hanger			Offset Hanger			
	Angle Adjustable / Multi-Angle	Composite	Composite Angle Adjustable	Angle Adjustable / Multi-Angle	Composite	Composite Angle Adjustable	
112009	10" - 13"	N/A	N/A	10" - 11"	N/A	N/A	17"
122956	10" - 13"	14" - 15"	N/A	10" - 11"	12" - 13"	N/A	15.5"
122957	14" - 15"	16" - 17"	15" - 16"	12" - 13"	14" - 15"	13" - 14"	16.4"
103364	14" - 17"	14" - 18"	15" - 17"	12" - 15"	12" - 16"	13" - 15"	19"
122958	16" - 24"	18" - 24"	17" - 24"	14" - 24"	16" - 24"	15" - 24"	17.6"
103365	18" - 20"	19" - 24"	18" - 20"	16" - 18"	17" - 24"	16" - 18"	21"
103366	21" - 24"	N/A	21" - 24"	19" - 24"	N/A	19" - 24"	25"

Item Number	Part Number	Description	Quantity
1a	100591	Heel Loop Strap Short	1
1b	100592	Heel Loop Strap Medium	1
1c	100593	Heel Loop Strap Long	1
1d	122956	Adjustable Heel Loop Small	1
1e	122957	Adjustable Heel Loop Medium	1
1f	122958	Adjustable Heel Loop Large	1
2	100551	Heel Loop Peg	1
3	101471	10-32 x 1-1/4 PH Pan HD T-B T/S ZC	1