


Rogue Transit



This Rogue Transit is effective starting with serial number RG0016612. <p>

For all mounting hardware, please see Rogue Backrest Mount and Hardware page. <p>

Rogue Transit cannot be retrofitted. If needing to add transit, please see the Rogue Rigid Frame page or call Customer Service.

Item Number	Part Number	Description	Quantity	Retail
2-6	112551	Kit Rogue Transit	1	\$102.45
2	003232	M6 x 55 BHCS Blk Zc	4	\$8.94
3	003715	Back Saddle 1 3/8" Tube	2	\$2.43
4	003716	Saddle 1 3/8" Tube	4	\$1.05
5	004434	Back Bracket with Transit, Black, Rigid	2	\$26.72
6	100658	M6 Nylock Nut Blk Zn	4	\$1.05