



## Intelect Legend XT System

- 2763** – 2-Channel Electrotherapy
- 2786K** – 4-Channel Electrotherapy
- 2760K** – 2-Channel Combination
- 2788K** – 4-Channel Combination

## Standard Accessories:

- 27465** – Patient Data Cards (5x)
- 27469** – Patient Interrupt Switch
- 27383** – 5 cm<sup>2</sup> Sound Head Applicator
- 27312** – Electrotherapy Leadwire Channel 1
- 27313** – Electrotherapy Leadwire Channel 2
- 42044** – DURA-STICK® II 2.75" (7 cm) Electrodes (4x)
- 27852** – Black Rubber Carbon 3" (8 cm) Electrodes (2x)
- 27853** – Red Rubber Carbon 3" (8 cm) Electrodes (2x)
- 10648** – Nylatex® Wrap (2x)
- 4248** – Conductor™ Ultrasound Gel - 8.5 oz (250 ml) Bottle (1x) (Combination Systems only)

## Intelect Legend XT System with Cart

- 2793K** – 2-Channel Electrotherapy
- 2797K** – 4-Channel Electrotherapy
- 2791K** – 2-Channel Combination
- 2795K** – 4-Channel Combination

## Optional Accessories

- 2758** – Light Therapy Module
- 2781** – Channel 3 and 4 Electrotherapy Module
- 2780ASY** – Intelect Legend XT Therapy Cart
- 2767** – Battery Module (Self-contained Nickel Metal Hydride battery automatically and continually charges when unit is plugged into mains power supply.)
- 27465** – Patient Data Cards (25 /pack)
- 2790** – Intelect Legend XT Patient Data Management System
- 27508** – Operator Remote Channels 1 and 2
- 27079** – Operator Remote Channels 3 and 4
- 27381** – 1 cm<sup>2</sup> Sound Head Applicator
- 27382** – 2 cm<sup>2</sup> Sound Head Applicator
- 27384** – 10 cm<sup>2</sup> Sound Head Applicator
- 27133** – Carry Bag
- 28205** – 19 Diode Cluster Applicator, 325mW Total: no Lasers, 6 x 670nm 10mW LED, 7 x 880nm 25mW LED, 6 x 950nm 15mW LED 33 Diode Cluster
- 28206** – 33 Diode Cluster Applicator, 565mW Total: no Lasers, 12 x 670nm 10mW LED, 13 x 880nm 25mW LED, 8 x 950nm 15mW LED

## Intelect® Legend XT

The Intelect® Legend XT is one of the most advanced therapy systems available. Beneath its award-winning, ergonomic design is a sophisticated, versatile, user-friendly system. Engineered to meet the electrotherapy and ultrasound clinician's needs, the Intelect Legend XT delivers a comprehensive range of therapies in a single, state-of-the-art unit.

# Intelect® LEGENDXT

- Patented award-winning, ergonomic modular design
- Easy-to-add modularity allows the addition of two additional Channels of Electrotherapy, Light Therapy or a Battery Module
- High contrast 5" (13 cm) FSTN LCD monochromatic user interface
- Two independent Electrotherapy channels, expandable to four independent channels
- Multiple waveform electrotherapy- six clinical waveforms: Interferential, Premodulated, High Volt, Russian, Symmetrical Biphasic and Microcurrent
- Documentation of treatment data with Patient Data Cards
- User defined protocols for your specific needs
- Dual frequency Ultrasound at 1 or 3.3 MHz
  - » Pulsed or Continuous Duty Cycles (10%, 20%, 50% & 100%)
  - » Four ergonomically designed sound head applicator sizes: 1 cm<sup>2</sup>, 2 cm<sup>2</sup>, 5 cm<sup>2</sup> & 10 cm<sup>2</sup>
  - » Soundheads are watertight for use in underwater therapy
- Therapy System Cart maximizes storage of clinical supplies and organizes cables
- Optional Operator Remotes allow treatment stop, start and pause, as well as change in amplitude, at any time during a treatment session. Each remote is color-coded for designation between channels



## PATENTED TECHNOLOGY

The Intelect Legend XT Therapy System is protected by one or more of the following patents:

D503,807	D511,008
D497,462	D513,615
D506,260	D504,593
D505,726	D508,998

Additional Patents Pending