


Catalyst E 22" Spoke Wheel / Tire / Handrim Kits



Not all components of the Wheel, Tire and Handrim kits will be listed. If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles may be needed.

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|---|----------|
| | 118602 | 22" Spoke, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, 5 Rivet Nut, Standard Mount | 1 |
| | 118603 | 22" Spoke, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, 5 Rivet Nut, Close Mount | 1 |
| | 118606 | 22" Spoke, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit w/ Micro Grip, 5 Rivet Nut, Standard Mount | 1 |
| | 118607 | 22" Spoke, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit w/ Micro Grip, 5 Rivet Nut, Close Mount | 1 |
| | 118612 | 22" Spoke, Pneumatic Tire (Black), Plastic Coated Handrim Kit, 6 Rivet Nut, Standard Mount | 1 |
| | 118613 | 22" Spoke, Pneumatic Tire (Black), Plastic Coated Handrim Kit, 6 Rivet Nut, Close Mount | 1 |
| | 118614 | 22" Spoke, Pneumatic Tire (Black), Projection Handrim Kit, Standard Mount | 1 |
| | 118615 | 22" Spoke, Pneumatic Tire (Black), Projection Handrim Kit, Close Mount | 1 |

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|--|----------|
| | 119766 | 22" Spoke, Pneumatic Tire (Black) Ergonomic Standard Handrim Regular Thumb Grip Kit | 1 |
| | 119769 | 22" Spoke, Pneumatic Tire (Black) Ergonomic Lite Handrim Regular Thumb Grip Kit | 1 |
| | 118622 | 22" Spoke, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount | 1 |
| | 118623 | 22" Spoke, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Close Mount | 1 |
| | 118624 | 22" Spoke, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit w/ Micro Grip, Standard Mount | 1 |
| | 118625 | 22" Spoke, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit w/ Micro Grip, Close Mount | 1 |
| | 118626 | 22" Spoke, Pneumatic w/ Airless Insert Tire (Black), Plastic Coated Handrim Kit, Standard Mount | 1 |
| | 118627 | 22" Spoke, Pneumatic w/ Airless Insert Tire (Black), Plastic Coated Handrim Kit, Close Mount | 1 |
| | 118628 | 22" Spoke, Pneumatic w/ Airless Insert Tire (Black), Projection Handrim Kit, Standard Mount | 1 |
| | 118629 | 22" Spoke, Pneumatic w/ Airless Insert Tire (Black), Projection Handrim Kit, Close Mount | 1 |
| | 119772 | 22" Spoke, Pneumatic w/ Airless Insert Tire (Black) Ergonomic Lite Handrim Regular Thumb Grip Kit | 1 |
| | 119775 | 22" Spoke, Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Regular Thumb Grip Kit | 1 |
| | 118636 | 22" Spoke, Low Poly Tire, Aluminum Anodized Handrim Kit, 5 Rivet Nut, Standard Mount | 1 |
| | 118637 | 22" Spoke, Low Poly Tire, Aluminum Anodized Handrim Kit, 5 Rivet Nut, Close Mount | 1 |
| | 118640 | 22" Spoke, Low Poly Tire, Aluminum Anodized Handrim Kit w/ Micro Grip, 5 Rivet Nut, Standard Mount | 1 |
| | 118641 | 22" Spoke, Low Poly Tire, Aluminum Anodized Handrim Kit w/ Micro Grip, 5 Rivet Nut, Close Mount | 1 |
| | 118646 | 22" Spoke, Low Poly Tire, Plastic Coated Handrim Kit, 6 Rivet Nut, Standard Mount | 1 |
| | 118647 | 22" Spoke, Low Poly Tire, Plastic Coated Handrim Kit, 6 Rivet Nut, Close Mount | 1 |
| | 119778 | 22" Spoke, Low Poly Tire Ergonomic Standard Handrim Regular Thumb Grip Kit | 1 |
| | 119781 | 22" Spoke, Low Poly Tire Ergonomic Lite Handrim Regular Thumb Grip Kit | 1 |
| | 118654 | 22" Spoke, Full Poly Tire, Aluminum Anodized Handrim Kit, 5 Rivet Nut, Standard Mount | 1 |
| | 118655 | 22" Spoke, Full Poly Tire, Aluminum Anodized Handrim Kit, 5 Rivet Nut, Close Mount | 1 |

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|---|----------|
| | 118658 | 22" Spoke, Full Poly Tire, Aluminum Anodized Handrim Kit w/ Micro Grip, 5 Rivet Nut, Standard Mount | 1 |
| | 118659 | 22" Spoke, Full Poly Tire, Aluminum Anodized Handrim Kit w/ Micro Grip, 5 Rivet Nut, Close Mount | 1 |
| | 118664 | 22" Spoke, Full Poly Tire, Plastic Coated Handrim Kit, 6 Rivet Nut, Standard Mount | 1 |
| | 118665 | 22" Spoke, Full Poly Tire, Plastic Coated Handrim Kit, 6 Rivet Nut, Close Mount | 1 |
| | 118666 | 22" Spoke, Full Poly Tire, Projection Handrim Kit, Standard Mount | 1 |
| | 118667 | 22" Spoke, Full Poly Tire, Projection Handrim Kit, Close Mount | 1 |
| | 119784 | 22" Spoke Full Poly Tire Ergonomic Standard Handrim Regular Thumb Grip Kit | 1 |
| | 119787 | 22" Spoke, Full Poly Tire Ergonomic Lite Handrim Regular Thumb Grip Kit | 1 |
| | 118692 | 22" Spoke, Sentinel Tire, Aluminum Anodized Handrim Kit, 5 Rivet Nut, Standard Mount | 1 |
| | 118693 | 22" Spoke, Sentinel Tire, Aluminum Anodized Handrim Kit, 5 Rivet Nut, Close Mount | 1 |
| | 118696 | 22" Spoke, Sentinel Tire, Aluminum Anodized Handrim Kit w/ Micro Grip, 5 Rivet Nut, Standard Mount | 1 |
| | 118697 | 22" Spoke, Sentinel Tire, Aluminum Anodized Handrim Kit w/ Micro Grip, 5 Rivet Nut, Close Mount | 1 |
| | 118702 | 22" Spoke, Sentinel Tire, Plastic Coated Handrim Kit, 6 Rivet Nut, Standard Mount | 1 |
| | 118703 | 22" Spoke, Sentinel Tire, Plastic Coated Handrim Kit, 6 Rivet Nut, Close Mount | 1 |
| | 119796 | 22" Spoke Sentinel Tire Ergonomic Standard Handrim Regular Thumb Grip Kit | 1 |
| | 119799 | 22" Spoke Sentinel Tire Ergonomic Lite Handrim Regular Thumb Grip Kit | 1 |
| | 118710 | 22" Spoke, High Pressure Tire, Aluminum Anodized Handrim Kit, 5 Rivet Nut, Standard Mount | 1 |
| | 118711 | 22" Spoke, High Pressure Tire, Aluminum Anodized Handrim Kit, 5 Rivet Nut, Close Mount | 1 |
| | 118714 | 22" Spoke, High Pressure Tire, Aluminum Anodized Handrim Kit w/ Micro Grip, 5 Rivet Nut, Standard Mount | 1 |
| | 118715 | 22" Spoke, High Pressure Tire, Aluminum Anodized Handrim Kit w/ Micro Grip, 5 Rivet Nut, Close Mount | 1 |
| | 118720 | 22" Spoke, High Pressure Tire, Plastic Coated Handrim Kit, 6 Rivet Nut, Standard Mount | 1 |
| | 118721 | 22" Spoke, High Pressure Tire, Plastic Coated Handrim Kit, 6 Rivet Nut, Close Mount | 1 |
| | 119802 | 22" Spoke High Pressure Tire Ergonomic Standard Handrim Regular Thumb Grip Kit | 1 |
| | 119805 | 22" Spoke High Pressure Tire Ergonomic Lite Handrim Regular Thumb Grip Kit | 1 |