


Catalyst Upper Pivot Assembly - Growth



Cross Brace Overall Width	
Seat Width	Cross Brace Measurement
14"-16"	Approx. 1" less than seat width
17"	Approx. 1.25" less than seat width
18"-21"	Approx. 1.5" less than seat width
22"	Approx. 1.75" less than seat width

Cross Brace Overall Depth	
Seat Depth	Cross Brace Measurement
14"-20"	Approx. 0.375" less than seat depth

Item Number	Part Number	Description	Quantity
1a, 2a, 3-5	500024	Upper Pivot Assembly 14" Width <i>Measures Approximately 4.5"</i>	2
1b, 2b, 3-5	500025	Upper Pivot Assembly 15" Width <i>Measures Approximately 5"</i>	2
1c, 2b, 3-5	500026	Upper Pivot Assembly 16" Width <i>Measures Approximately 5.5"</i>	2
1d, 2b, 3-5	500027	Upper Pivot Assembly 17" Width <i>Measures Approximately 6"</i>	2
1e, 2b, 3-5	500028	Upper Pivot Assembly 18" Width <i>Measures Approximately 6.5"</i>	2

Item Number	Part Number	Description	Quantity
1f, 2b, 3-5	500029	Upper Pivot Assembly 19" Width <i>Measures Approximately 7"</i>	2
1g, 2b, 3-5	500030	Upper Pivot Assembly 20" Width <i>Measures Approximately 7.25"</i>	2
1h, 2b, 3-5	501479	Upper Pivot Assembly 21/22" Width <i>Measures Approximately 8.25"</i>	2