


# Discontinued 1 Rogue XP Tension Adjustable Back Upholstery 17" - 18" Wide



Back Height	Strap Qty (Item 6)	D-Ring Qty (Item 7)
9"-12"	3	6
11"-14"	3	6
13"-16"	4	8
15"-18"	5	10

*This Tension Adjustable Back Upholstery was effective prior to November 2016. Please see Rogue XP Tension Adjustable Back Upholstery 17" - 18" Wide page for current assemblies.*

<b>Item Number</b>	<b>Part Number</b>	<b>Description</b>	<b>Quantity</b>
7	100973	Back Upholstery D-Ring	6
8	100734	M5 x 16 PHMS ZC	2
9a	001751	Tension Strap Folding Handle 16-17 Wide	1
9b	001752	Tension Strap Folding Handle 18-19 Wide	1