


LECKEY[®]

Freestander



Our child centred approach is not just about sizing equipment to fit the child, it's about understanding and designing for their needs at each level: infants, kids and teens.



Designed to promote weight bearing and mid line symmetry, the Leckey Freestander has been invaluable to carers and therapists working with children and adults for many years.





The Totstander is our pre-school version which has been developed to suit the needs of children at the early stages of standing. (Only available in Meadow Fabric)

Fabric Options

The Freestander comes in 2 fabrics. Size 1 comes as standard in Meadow fabric to suit the nursery environment. Size 2 and 3 in Mondeo blue fabric.



Standard Features



Height and width adjustable flexible lateral supports.



Padded Chest plate adjusts in height and angle for multi positioning.



Multi adjustable pelvic band.



Height and angle adjustable kneeplate with laterally adjustable padded kneecups.



Low base board for easy "step up" during assisted transfer.



Multi adjustable sandals.



Height and angle adjustable activity tray, with bowl and cover.



Four lockable swivel castors.

Accessories



LFS/#/C
Suspension Rail/
Grab Rail



LFS/#/L
Solid Pelvic Plate



LFS/#/N
Grab Posts



Interface frame - LIF/1
Contoured head support - LIF/2/A, B, C
Hyper extension pads - LIF/3/A, B



Whitmyer Pro Series
Head Support System



Whitmyer Pro Series
Plush Headrest System



Whitmyer Contoured cradle

Contact the Leckey Customer Service Centre for further details on the Whitmyer range.

Leckey
Kilwee Business Park
Dunmurry
BT17 OHD
Northern Ireland
United Kingdom

T: 028 9060 2277
F: 028 9060 0795
E: info@leckey.com
W: www.leckey.com

Freephone:
United Kingdom 0800 318265
Ireland 1800 626020

Product sizing						
Size	1		2		3	
Age	1 - 5		4 - 9 1/2		9 - 17	
Codes	LFS/1		LFS/2		LFS/3	
	mm	inches	mm	inches	mm	inches
Freestander						
Chest Plate Height						
min	470	18.5	680	27	975	38
max	680	27	975	38	1350	53
Base						
length	650	26	750	30	750	30
width	650	26	750	30	750	30
Activity Tray						
depth	400	16	450	18	500	20
width	525	21	575	23	625	25
User Weight						
max	30kg	66lbs	60kg	132lbs	100kg	220lbs
Totstander						
Chest Band Height						
min	470	18.5				
max	640	26				
Base						
depth	450	18				
width	450	18				