


# 1" Auto Buckle Pelvic Positioning Belt



| Item Number | Part Number | Description  | Quantity |
|-------------|-------------|--|----------|
| 1, 2a, 3, 4 | 111381      | Belt Pelvic Positioning 1" Auto Buckle 51", Catalyst/Spark                         | 1        |
| 1, 2b, 3, 4 | 111382      | Belt Pelvic Positioning 1" Auto Buckle 51", Focus                                  | 1        |
| 1, 2c, 3, 4 | 111383      | Belt Pelvic Positioning 1" Auto Buckle 51", RG ALX/Rogue/LW XP                     | 1        |
| 1, 2d, 3, 4 | 111384      | Belt Pelvic Positioning 1" Auto Buckle 51", Clik/Ethos                             | 1        |
| 1, 2e, 3, 4 | 111385      | Belt Pelvic Positioning 1" Auto Buckle 51", Flip/Liberty                           | 1        |
| 1           | 100614      | Belt Pelvic Positioning 1" Auto Buckle 51"   | 1        |
| 2a          | 100681      | M6 x 40 CLS8.8 DIN931 HHCS ZC<br><i>Used on Catalyst and Spark</i>                 | 2        |
| 2b          | 101125      | M6 x 55 BTN SHCS BLKZC<br><i>Used on Focus</i>                                     | 2        |
| 2c          | 101781      | M6 x 50mm Partial Thread BHCS Blk ZC<br><i>Used on Tsunami, Rogue and Rogue XP</i> | 2        |

| Item Number | Part Number | Description  | Quantity |
|-------------|-------------|--|----------|
| 2d          | 100670      | M6 x 45 BTN SHCS BLK ZC<br><i>Used on Clik and Ethos</i> | 2        |
| 2e          | 101811      | M6 x 35 BHCS Blk Zn<br><i>Used on Flip and Liberty</i>   | 2        |
| 3           | 100746      | M6 Flat Washer Blk Zn                                    | 2        |
| 4           | 100658      | M6 Nylock Nut Blk Zn                                     | 2        |