


# Technical Data

## FITTING SHOULDER HARNESSSES AND ANTERIOR TRUNK SUPPORTS




<b>XS</b>	9-1/2"-11" (24cm-28cm)
<b>S</b>	11"-13" (28cm-33cm)
<b>M</b>	13"-16" (33cm-41cm)
<b>L</b>	16"-19" (41cm-48cm)
<b>XL</b>	19"-22" (48cm-56cm)

Measure shoulder width. If the shoulder width falls between two sizes, consider other factors such as growth, weight changes, and clothing.


## SPECIFICATIONS: Shoulder Harness pad length

	Standard	Trimline
<b>S</b>	12" (31cm)	9-1/4" (24cm)
<b>M</b>	13-1/2" (34cm)	11-1/2" (29cm)
<b>L</b>	16" (41cm)	13-1/2" (34cm)
<b>XL</b>	18-1/2" (47cm)	

## SPECIFICATIONS: Stayflex™



	Standard	Size	A	B	C	D	E
<b>Standard Stayflex (w/zipper)</b>	SH310S	S	9 (229mm)	5 1/2 (140mm)	6 1/2 (165mm)	1 3/4 (44mm)	4 3/4 (121mm)
	SH310M	M	12 (305mm)	7 (187mm)	8 (203mm)	1 3/4 (44mm)	5 (127mm)
	SH310L	L	15 (381mm)	8 1/2 (216mm)	9 (229mm)	2 1/4 (57mm)	5 1/4 (133mm)
	SH310XL	XL	18 (457mm)	10 (254mm)	10 1/2 (267mm)	2 1/4 (57mm)	6 1/2 (165mm)
<b>Standard One Piece (w/o zipper)</b>	SH311XS	XS	8 1/4 (210mm)	5 (127mm)	5 1/4 (133mm)	1 1/2 (38mm)	3 3/4 (95mm)
	SH311S	S	8 1/4 (210mm)	5 3/4 (146mm)	6 1/2 (165mm)	1 3/4 (44mm)	4 (102mm)
	SH311M	M	11 1/4 (286mm)	7 1/4 (184mm)	8 (203mm)	1 3/4 (44mm)	4 1/4 (108mm)
	SH311L	L	14 1/4 (362mm)	8 3/4 (222mm)	9 (229mm)	2 1/4 (57mm)	4 1/2 (114mm)
	SH311XL	XL	17 (432mm)	10 1/4 (260mm)	10 1/2 (267mm)	2 1/4 (57mm)	5 1/2 (140mm)



	Narrow	Size	A	B	C	D	E
<b>Narrow Stayflex (w/zipper)</b>	SH330M	M	12 (305mm)	7 (178mm)	8 (203mm)	1 3/4 (44mm)	3 1/2 (89mm)
	SH330L	L	15 (381mm)	8 1/2 (216mm)	9 (229mm)	2 1/4 (57mm)	3 1/2 (89mm)
	SH330XL	XL	18 (457mm)	10 (254mm)	10 1/2 (267mm)	2 1/4 (57mm)	3 3/4 (95mm)
<b>Narrow One Piece (w/o zipper)</b>	SH331M	M	11 1/4 (286mm)	7 1/4 (184mm)	8 (203mm)	1 3/4 (44mm)	2 3/4 (70mm)
	SH331L	L	14 1/4 (362mm)	8 3/4 (222mm)	9 (229mm)	2 1/4 (57mm)	2 3/4 (70mm)
	SH331XL	XL	17 (432mm)	10 1/4 (260mm)	10 1/2 (267mm)	2 1/4 (57mm)	3 (76mm)