


(Discontinued 1) Rogue Backrest Mount and Hardware



Adjustment Tool required when ordering Spinergy Back Release Cable

This Rogue Backrest Mount is discontinued and was used on serial numbers prior to RG0016612 This backrest mount and hardware is used on the Discontinued Rogue with XP option frame. If you have a Rogue XP, use the Rogue XP Backrest Mount pages for replacement parts.

Item Number	Part Number	Description	Quantity
1b, 2b, 6-8, 10, 11a, 12-19, 21	500759_ _ _	Color Anodized Backrest Plate Kit	1
1b	102600_ _ _	Color Anodized Backrest Mounting Plate Left RG	1
2a, 6	000958R	Backrest Mounting Plate w/ Pin Right Black Anodized RG	1

Item Number	Part Number	Description	Quantity
2b	102601_ _ _ _	Color Anodized Backrest Mounting Plate Right RG	1
3a	000910	Backrest Pin Lock Housing Black Anodized RG	2
4	000909	Backrest Lockout Pin RG	2
5	000965	Backrest Cable Release Spring RG	2
6	000922	Backrest Stop RG	2
7	000248	Saddle 1.375"	2
8	000055	Saddle 1"	4
9a	200260	Spinergy Back Release Cable 12" Wide Black	1
9b	200261	Spinergy Back Release Cable 12" Wide Blue	1
9c	200262	Spinergy Back Release Cable 12" Wide Red	1
9d	200263	Spinergy Back Release Cable 12" Wide Yellow	1
9e	200264	Spinergy Back Release Cable 12" Wide Green	1
9f	200265	Spinergy Back Release Cable 12" Wide Pink	1
9g	200266	Spinergy Back Release Cable 12" Wide Orange	1
9h	200267	Spinergy Back Release Cable 12" Wide White	1
9q	200268	Spinergy Back Release Cable 13" Wide Black	1
9r	200269	Spinergy Back Release Cable 13" Wide Blue	1
9s	200270	Spinergy Back Release Cable 13" Wide Red	1
9t	200271	Spinergy Back Release Cable 13" Wide Yellow	1
9u	200272	Spinergy Back Release Cable 13" Wide Green	1
9v	200273	Spinergy Back Release Cable 13" Wide Pink	1
9w	200274	Spinergy Back Release Cable 13" Wide Orange	1
9x	200275	Spinergy Back Release Cable 13" Wide White	1
9ag	200276	Spinergy Back Release Cable 14" Wide Black	1
9ah	200277	Spinergy Back Release Cable 14" Wide Blue	1
9ai	200278	Spinergy Back Release Cable 14" Wide Red	1
9aj	200279	Spinergy Back Release Cable 14" Wide Yellow	1
9ak	200280	Spinergy Back Release Cable 14" Wide Green	1

Item Number	Part Number	Description	Quantity
9al	200281	Spinergy Back Release Cable 14" Wide Pink	1
9am	200282	Spinergy Back Release Cable 14" Wide Orange	1
9an	200283	Spinergy Back Release Cable 14" Wide White	1
9aw	200284	Spinergy Back Release Cable 15" Wide Black	1
9ax	200285	Spinergy Back Release Cable 15" Wide Blue	1
9ay	200286	Spinergy Back Release Cable 15" Wide Red	1
9az	200287	Spinergy Back Release Cable 15" Wide Yellow	1
9ba	200288	Spinergy Back Release Cable 15" Wide Green	1
9bb	200289	Spinergy Back Release Cable 15" Wide Pink	1
9bc	200290	Spinergy Back Release Cable 15" Wide Orange	1
9bd	200291	Spinergy Back Release Cable 15" Wide White	1
9bm	200292	Spinergy Back Release Cable 16" Wide Black	1
9bn	200293	Spinergy Back Release Cable 16" Wide Blue	1
9bo	200294	Spinergy Back Release Cable 16" Wide Red	1
9bp	200295	Spinergy Back Release Cable 16" Wide Yellow	1
9bq	200296	Spinergy Back Release Cable 16" Wide Green	1
9br	200297	Spinergy Back Release Cable 16" Wide Pink	1
9bs	200298	Spinergy Back Release Cable 16" Wide Orange	1
9bt	200299	Spinergy Back Release Cable 16" Wide White	1
9cc	200300	Spinergy Back Release Cable 17" Wide Black	1
9cd	200301	Spinergy Back Release Cable 17" Wide Blue	1
9ce	200302	Spinergy Back Release Cable 17" Wide Red	1
9cf	200303	Spinergy Back Release Cable 17" Wide Yellow	1
9cg	200304	Spinergy Back Release Cable 17" Wide Green	1
9ch	200305	Spinergy Back Release Cable 17" Wide Pink	1
9ci	200306	Spinergy Back Release Cable 17" Wide Orange	1
9cj	200307	Spinergy Back Release Cable 17" Wide White	1
9cs	200308	Spinergy Back Release Cable 18" Wide Black	1

Item Number	Part Number	Description	Quantity
9ct	200309	Spinergy Back Release Cable 18" Wide Blue	1
9cu	200310	Spinergy Back Release Cable 18" Wide Red	1
9cv	200311	Spinergy Back Release Cable 18" Wide Yellow	1
9cw	200312	Spinergy Back Release Cable 18" Wide Green	1
9cx	200313	Spinergy Back Release Cable 18" Wide Pink	1
9cy	200314	Spinergy Back Release Cable 18" Wide Orange	1
9cz	200315	Spinergy Back Release Cable 18" Wide White	1
9di	200316	Spinergy Back Release Cable 19" Wide Black	1
9dy	200324	Spinergy Back Release Cable 20" Wide Black	1
9eo	200332	Spinergy Back Release Cable 21" Wide Black	1
9fe	200340	Spinergy Back Release Cable 22" Wide Black	1
10	100670	M6 x 45 BTN SHCS BLK ZC	2
11a	003232	M6 x 55 BHCS Blk Zc <i>Used on chairs with 0" or positive center of gravity</i>	2
11b	003241	M6 x 60 BHCS P/THD Custom Blk Zc <i>Used on chairs with negative center of gravity</i>	2
12	101716	M5 x 16 SHCS Blk Zn	2
13	101121	M5 x 40 CLS8.8 DIN7984 Low SHCS BLKZC	2
14	100797	M6 x 40 BHCS BLK ZC	2
15	100658	M6 Nylock Nut Blk Zn	4
16	101856	Flat Washer, M5, Blk Zn	2
17	100657	M5 Nylock Nut Blk Zn	2
18	101456	M6 Ovl Top L/N ZC	2
19	102072	Washer Nylon .234 ID x .406 OD x .062 Thick	2
20	101991	Spinergy Release Cable/SPOX Adjustment Tool <i>No charge when ordered with Spinergy Release Cable</i>	1
21	000369	Flat Washer .5 x .257 x .0032 Plastic	2

Color Anodizing Samples

