

Rogue ALX Wheels - SpinerGY Spox (Formerly Tsunami)



Recommended to order the SpinerGY Spox Replacement Spoke Tool Kit (110118) when replacing spokes. Replacement SpinerGY Spokes may have an extended lead time. If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel style, new axles may be needed. Please note tab mounted handrims must be installed prior to the tire being installed.

Item Number	Part Number	Description	Quantity
1a	100922	22" x 1/2" Spinergy Spox Black Spoke	1
1b	100924	22" x 1/2" Spinergy Spox Red Spoke	1
1c	100926	22" x 1/2" Spinergy Spox White Spoke	1
1d	100923	22" x 1/2" Spinergy Spox Blue Spoke	1
1f	100862	24" x 1/2" Spinergy Spox Black Spoke	1
1f	118343	24" x 1/2" Spinergy Spox Black Spoke w/ Pneumatic w/ Airless Insert Tire (Black)	1
1g	100928	24" x 1/2" Spinergy Spox Red Spoke	1
1g	118344	24" x 1/2" Spinergy Spox Red Spoke w/ Pneumatic w/ Airless Insert Tire (Black)	1
1h	100929	24" x 1/2" Spinergy Spox White Spoke	1
1h	118345	24" x 1/2" Spinergy Spox White Spoke w/ Pneumatic w/ Airless Insert Tire (Black)	1
1i	100927	24" x 1/2" Spinergy Spox Blue Spoke	1
1i	118347	24" x 1/2" Spinergy Spox Blue Spoke w/ Pneumatic w/ Airless Insert Tire (Black)	1
1k	100989	25" x 1/2" Spinergy Spox Black Spoke	1
1l	100991	25" x 1/2" Spinergy Spox Red Spoke	1
1m	100993	25" x 1/2" Spinergy Spox White Spoke	1
1n	100990	25" x 1/2" Spinergy Spox Blue Spoke	1
1p	100930	26" x 1/2" Spinergy Spox Black Spoke	1
1q	100932	26" x 1/2" Spinergy Spox Red Spoke	1
1r	100934	26" x 1/2" Spinergy Spox White Spoke	1
1s	100931	26" x 1/2" Spinergy Spox Blue Spoke	1
2	101991	Spinergy Release Cable/SPOX Adjustment Tool <i>No charge when ordered with Spinergy wheel</i>	1
	110118	Spinergy Spox Replacement Spoke Tool Kit <i>This tool kit contains: Spinergy Release Cable/SPOX Adjustment Tool, Spinergy Nipple Hex T-Wrench and Spinergy Spoke T-Wrench</i>	1
3a	101909	Replacement Spoke Spinergy SPOX 22", Black (511004-8.75BLA)	1
3b	101910	Replacement Spoke Spinergy SPOX 22", Red (511004-8.75RED)	1
3c	101911	Replacement Spoke Spinergy SPOX 22", White (511004-8.75WHT)	1

Item Number	Part Number	Description	Quantity
3d	101912	Replacement Spoke Spinergy SPOX 22", Blue (511004-8.75BLU)	1
3f	101914	Replacement Spoke Spinergy SPOX 24", Black (511004-9.5BLA)	1
3g	101915	Replacement Spoke Spinergy SPOX 24", Red (511004-9.5RED)	1
3h	101916	Replacement Spoke Spinergy SPOX 24", White (511004-9.5WHT)	1
3i	101917	Replacement Spoke Spinergy SPOX 24", Blue (511004-9.5BLU)	1
3k	101919	Replacement Spoke Spinergy SPOX 25", Black (511004-9.84BLA)	1
3l	101920	Replacement Spoke Spinergy SPOX 25", Red (511004-9.84RED)	1
3m	101921	Replacement Spoke Spinergy SPOX 25", White (511004-9.84WHT)	1
3n	101922	Replacement Spoke Spinergy SPOX 25", Blue (511004-9.84BLU)	1
3p	101924	Replacement Spoke Spinergy SPOX 26", Black (511004-10.47BLA)	1
3q	101925	Replacement Spoke Spinergy SPOX 26", Red (511004-10.47RED)	1
3r	101926	Replacement Spoke Spinergy SPOX 26", White (511004-10.47WHT)	1
3s	101927	Replacement Spoke Spinergy SPOX 26", Blue (511004-10.47BLU)	1

Handrim Hardware Chart													
Wheel		Handrim Connection				Handrim							
Wheel	Wheel Part No.	Connection Points	Tab	Spacer*	Screw*	Aluminum Anodized	Superlight	Plastic Coated	Projection	Natural Fit Standard	Natural Fit LT	Flex Rim	
18" MAXX Spoke	109285	3	100698 (Not used w/ Superlight Handrim)	100653 (Not used w/ Superlight Handrim)	Aluminum, Plastic Ctd, Projection: 100654	100206		101106					
20" MAXX Spoke	108244	6			Natural Fit: 100835	200536		200542					
22" MAXX Spoke	105135				Superlight: Screw: 100669 Nut: 100657	100560		100576	100569	200538	200201		
24" MAXX Spoke	105136				200349	101161	100577	200547	100793	200202			
25" MAXX Spoke	107436				200350		101870	200548	200539	200540			
26" MAXX Spoke	107437				200351		100578	200549	100907	101454			
20" MAXX Mag	111853		6		100629	Aluminum, Plastic Ctd, Projection: 103545	200536		112819				
22" MAXX Mag	111854	100560					112820	112824	200538	200201			
24" MAXX Mag	111855	Natural Fit: 101893				105137		112821	112825	100793	200202		
24" Superlight	101159	6			100536	100754	101161	100836		100830	100828		
25" Superlight	101160						101160	101091		101464	101460		
20" Spinerger Spox	See Spinerger Spox Page	6			Screw: 100669 Nut: 100657	103125		103179					
22" Spinerger Spox						100827		100808		100889	100888		
24" Spinerger Spox						100766	101161	100615		100830	100828		
25" Spinerger Spox						100767	101160	100765		101464	101460		
26" Spinerger Spox						101477		101148		200200	100950		
22" Spinerger LX	See Spinerger LX Page	6			Screw: 100669 Nut: 100657	100827		100808		100889	100888		
24" Spinerger LX						100766	101161	100615		100830	100828	See Spinerger LX Page	
25" Spinerger LX						100767	101160	100765		101464	101460		
26" Spinerger LX						101477		101148		200200	100950		

* Spacer and screw part numbers listed in the chart are for standard handrim mounting. If using close mount handrim mounting on aluminum anodized, plastic coated or projection handrims, use part number 100792 for the spacer, part number 102616 for the screw on Ki Spoke and MAXX Spoke wheels, part number 100666 for the screw on 5-Spoke X Core wheels, or part number 100654 for the screw on MAXX Mag wheels.

See Page 2 for Discontinued Wheels.

Effective 9/6/2022

Handrim Hardware Chart (Discontinued Wheels)

Wheel		Handrim Connection				Handrim							
Wheel	Wheel Part No.	Connection Points	Tab	Spacer*	Screw*	Aluminum Anodized	Superlight	Plastic Coated	Projection	Natural Fit Standard	Natural Fit LT	Flex Rim	
18" Ki Spoke	200529	3	100698	100653	Aluminum, Plastic Ctd, Projection: 100654	100206		101106					
20" Ki Spoke	200530	6				200536		200542					
22" Ki Spoke	200531					100560		100576	100569	200538	200201		
24" Ki Spoke	200532					200349		100577	200547	100793	200202		
25" Ki Spoke	200533				Natural Fit: 100835	200350		101870	200548	200539	200540		
26" Ki Spoke	200534					200351		100578	200549	100907	101454		
20" 5-Spoke X Core	101961		5		Aluminum, Projection, Natural Fit: 100629	100724	101897		101963				
22" 5-Spoke X Core	101962			101898				101964					
24" 5-Spoke X Core	100960			100975				100976	200546				
24" 5-Spoke X Core	100960			Plastic Ctd: 101756	101893				100768	100769			

* Spacer and screw part numbers listed in the chart are for standard handrim mounting. If using close mount handrim mounting on aluminum anodized, plastic coated or projection handrims, use part number 100792 for the spacer, part number 102616 for the screw on Ki Spoke and MAXX Spoke wheels, part number 100666 for the screw on 5-Spoke X Core wheels, or part number 100654 for the screw on MAXX Mag wheels.

Tire Component Chart

Wheel		1 3/8" Street Pneumatic			1" High Pressure			1 3/8" Euro Trek Knobby			1" Schwalbe				1 3/8" Schwalbe HS-127			2" Mountain			
Wheel	Wheel Part #	Tire	Inner Tube	Rim Strip	Tire	Inner Tube	Rim Strip	Tire	Inner Tube	Rim Strip	Mara-thon Tire	One Tire	Inner Tube	Rim Strip	Tire	Inner Tube	Rim Strip	Tire	Inner Tube	Rim Strip	
18" MAXX Spoke	109285	101587	101588	101137																	
20" MAXX Spoke	108244	100911	100912	101137	102957	102844	101137														
22" MAXX Spoke	105135	100634	100588	101908	102961	100609	101908											101475	100588	101908	
24" MAXX Spoke	105136	111818	100589	100596	102958	100580	100596	102982	100589	100596	100864	103176	100580	100596	102920	100589	100596	101476	101188	100596	
25" MAXX Spoke	107436				102959	101457	101452				101357	103177	101457	101452				100166	100208	101452	
26" MAXX Spoke	107437	100636	100590	100608	102960	100582	100597				100865		100582	100597							
20" MAXX Mag	111853	100911	100912	101137																	
22" MAXX Mag	111854	100634	100588	101908														101475	100588	101908	
24" MAXX Mag	111855	111818	100589	100596				102982	100589	100596					102920	100589	100596	101476	101188	100596	
24" Superlight	101159	111818	100589	100596	102958	100580	100596	102982	100589	100596	100864	103176	100580	100596	102920	100589	100596				
25" Superlight	101158				102959	101457	101452				101357	103177	101457	101452							
20" Spinery Spox	See Spinery Spox Page	100911	100912	101137	102957	102844	101137														
22" Spinery Spox		100634	100588	101908	102961	100609	101908											101475	100588	101908	
24" Spinery Spox		111818	100589	100596	102958	100580	100596	102982	100589	100596	100864	103176	100580	100596	102920	100589	100596	101476	101188	100596	
25" Spinery Spox					102959	101457	101452					101357	103177	101457	101452				100166	100208	10145
26" Spinery Spox		100636	100590	100608	102960	100582	100597					100865		100582	100597						
22" Spingery LX	See Spinery LX Page	100634	100588	101908	102961	100609	101908											101475	100588	101908	
24" Spingery LX		111818	100589	100596	102958	100580	100596	102982	100589	100596	100864	103176	100580	100596	102920	100589	100596	101476	101188	100596	
25" Spingery LX					102959	101457	101452					101357	103177	101457	101452				100166	100208	101452
26" Spingery LX		100636	100590	100608	102960	100582	100597					100865		100582	100597						

See Page 2 for Discontinued Wheels.

Tire Component Chart (Discontinued Wheels)

Wheel		1 3/8" Street Pneumatic			1" High Pressure			1 3/8" Euro Trek Knobby			1" Schwalbe				1 3/8" Schwalbe HS-127			2" Mountain			
Wheel	Wheel Part #	Tire	Inner Tube	Rim Strip	Tire	Inner Tube	Rim Strip	Tire	Inner Tube	Rim Strip	Mara-thon Tire	One Tire	Inner Tube	Rim Strip	Tire	Inner Tube	Rim Strip	Tire	Inner Tube	Rim Strip	
18" Ki Spoke	200529	101587	101588	101137																	
20" Ki Spoke	200530	100911	100912	101137	102957	102844	101137														
22" Ki Spoke	200531	100634	100588	101908	102961	100609	101908											101475	100588	101908	
24" Ki Spoke	200532	111818	100589	100596	102958	100580	100596	102982	100589	100596	100864	103176	100580	100596	102920	100589	100596	101476	101188	100596	
25" Ki Spoke	200533				102959	101457	101452				101357	103177	101457	101452				100166	100208	101452	
26" Ki Spoke	200534	100636	100590	100608	102960	100582	100597				100865		100582	100597							
20" 5-Spoke X Core	101961	100911	100912	101137																	
22" 5-Spoke X Core	101962	100634	100588	101908																	
24" 5-Spoke X Core	100960	111818	100589	100596				102982	100589	100596					102920	100589	100596				

Tire Compatibility Chart

Wheel		Tire											
Wheel	Wheel Part #	1 3/8" Street Pneumatic	1 3/8" Full Profile Poly	1" Low Profile Poly	1 3/8" Pneumatic w/ Airless Insert (Tire / Insert)	1" Solid Sentinel	1" Solid Shox	1" High Pressure	1 3/8" Euro Trek Knobby	1" Schwalbe Marathon	1" Schwalbe One	1 3/8" Schwalbe HS-127	2" Mountain
18" MAXX Spoke	109285	101587	101992	100918	101587 / 101092								
20" MAXX Spoke	108244	100911	100917	102085	100911 / 100913	102974	102717	102957					
22" MAXX Spoke	105135	100634	100621	100632	100634 / 100503	103165	102718	102961					101475
24" MAXX Spoke	105136	111818	100622	100633	111818 / 100504	102976	100640	102958	102982	100864	103176	102920	101476
25" MAXX Spoke	107436					102977	101453	102959		101357	103177		100166
26" MAXX Spoke	107437	100636			100636 / 100505	102978	100884	102960		100865			
12" Mag ¹	N/A	001424	001426		001425								
16" Mag ¹	N/A	001159	001161	003453	001160								
20" MAXX Mag	111853	100911	100917	102085	100911 / 100913	102974							
22" MAXX Mag	111584	100634	100621	100632	100634 / 100503	103165							101475
24" MAXX Mag	111855	111818	100622	100633	111818 / 100504	102976			102982			102920	101476
24" Superlight	101159	100635					100894	102958	102982	100864	103176	102920	
25" Superlight	101158						101146	102959		101357	103177		
20" Spinerger Spox	See Spinerger Spox Page	100911	100917	102085	100911 / 100913	102974		102957					
22" Spinerger Spox		100634	100621	100632	100634 / 100503	103165		102961					101475
24" Spinerger Spox		111818	100622	100633	100818 / 100504	102976	100894	102958	102982	100864	103176	102920	101476
25" Spinerger Spox						102977	101146	102959		101357	103177		100166
26" Spinerger Spox		100636				100636 / 100505	102978	100895	102960		100865		
22" Spinerger LX	See Spinerger LX Page	100634	100621	100632	100634 / 100503			102961					101475
24" Spinerger LX		111818	100622	100633	111818 / 100504	102976	100894	102958	102982	100864	103176	102920	101476
25" Spinerger LX						102977	101146	102959		101357	103177		100166
26" Spinerger LX		100636				100636 / 100505	102978	100895	102960		100865		

¹ 12" and 16" Mag Wheels are sold as a combination Wheel and Tire Assembly

See Page 2 for Discontinued Wheels

Tire Compatibility Chart (Discontinued Wheels)

Wheel		Tire											
Wheel	Wheel Part #	1 3/8" Street Pneumatic	1 3/8" Full Profile Poly	1" Low Profile Poly	1 3/8" Pneumatic w/ Airless Insert (Tire / Insert)	1" Solid Sentinel	1" Solid Shox	1" High Pressure	1 3/8" Euro Trek Knobby	1" Schwalbe Marathon	1" Schwalbe One	1 3/8" Schwalbe HS-127	2" Mountain
18" Ki Spoke	200529	101587	101992	100918	101587 / 101092								
20" Ki Spoke	200530	100911	100917	102085	100911 / 100913	102974	102717	102957					
22" Ki Spoke	200531	100634	100621	100632	100634 / 100503	103165	102718	102961					101475
24" Ki Spoke	200532	111818	100622	100633	111818 / 100504	102976	100640	102958	102982	100864	103176	102920	101476
25" Ki Spoke	200533					102977	101453	102959		101357	103177		100166
26" Ki Spoke	200534	100636			100636 / 100505	102978	100884	102960		100865			
20" 5-Spoke X Core	101961	100911	100917	102085	100911 / 100913		102717						
22" 5-Spoke X Core	101962	100634	100621	100632	100634 / 100503		102718						
24" 5-Spoke X Core	100960	111818	100622	100633	111818 / 100504		100640		102982			102920	

Effective 11/2/2022