


Axiom Fixed Laterals



Assembled View

Item Number	Part Number	Description	Quantity	Retail
1-2, 3c, 3d, 4a-4d, 5-8, 9a	LF3X4L	3 x 4 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Lateral and Axiom Lateral Backs (Pair)	1	\$313.00
1-2, 3a, 3b, 4a-4d, 5-8, 9b	LF3X4P	3 x 4 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Back (Pair)	1	\$313.00
1-2, 3c, 3d, 4a-4d, 5-8, 9b	LF4X5L	4 x 5 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Lateral and Axiom Lateral Backs (Pair)	1	\$313.00
1-2, 3a, 3b, 4a-4d, 5-8, 9b	LF4X5P	4 x 5 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Back (Pair)	1	\$313.00
1	100558	M6 NYL INS Domed Cap L/N ZC	3	\$1.10
2	100746	M6 Flat Washer Blk Zn	5	\$1.66
3a	004880	Axiom PB Back Lateral Mount, Fixed, Right	1	\$26.78
3b	004881	Axiom PB Back Lateral Mount, Fixed, Left	1	\$26.78
3c	004882	Axiom PL Back Lateral Mount, Fixed, Right	1	\$37.49

Item Number	Part Number	Description	Quantity	Retail
3d	004883	Axiom PL Back Lateral Mount, Fixed, Left	1	\$37.49
4a	104017	1/4-20 X 1/2 Stainless Button Head Cap Screw <i>Used when no spacers are used.</i>	2	\$3.21
4b	104028	1/4-20 X 3/4 Stainless Button Head Cap Screw <i>Used when the 1/4" spacer is used alone.</i>	2	\$5.36
4c	104014	BHCS, 1/4"-20 X 1.00, 18-8, BLK OX, ASME B18.3 <i>Used when the 1/2" spacer is used alone.</i>	2	\$5.36
4d	104024	1/4-20 X 1-1/4 Stainless Button Head Cap Screw <i>Used when the 1/4" and 1/2" spacers are used together.</i>	2	\$7.50
5	102803	Washer Flat .250 x .500 x .027 HX BLZ	5	\$1.66
6	004889	Axiom Back Lateral Stud Plate	1	\$16.07
7	107555	Axiom Back Lateral Spacer 1/4"	1	\$32.13
8	107556	Axiom Back Lateral Spacer 1/2"	1	\$32.13
9a	108261	3x4 Lateral Pads (Pair)	1	\$192.78
9a	108262	3x4 Lateral Pad (Each)	1	\$96.39
9b	108259	4x5 Lateral Pads (Pair)	1	\$192.78
9b	108260	4x5 Lateral Pad (Each)	1	\$96.39
1-2, 3a, 4a-4d, 5-8	108249	Posterior - Fixed Lateral Right	1	\$69.62
1-2, 3b, 4a-4d, 5-8	108248	Posterior - Fixed Lateral Left	1	\$69.62
1-2, 3c, 4a-4d, 5-8	108252	Posterior Lateral - Fixed Lateral Right	1	\$69.62
1-2, 3d, 4a-4d, 5-8	108251	Posterior Lateral - Fixed Lateral Left	1	\$69.62