


## Ethos Wheels - Spinergy LX



Recommended to order the Spinergy LX Replacement Spoke Tool Kit (110119) when replacing spokes. Replacement Spinergy Spokes and Spinergy Wheels with Flex Rims will have an extended lead time. If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles may be needed.

| Item Number | Part Number | Description                        | Quantity |
|-------------|-------------|------------------------------------|----------|
| 1a          | 100886      | 22" x 1/2" Spinergy LX Black Spoke | 1        |
| 1b          | 100939      | 22" x 1/2" Spinergy LX Red Spoke   | 1        |
| 1c          | 100940      | 22" x 1/2" Spinergy LX White Spoke | 1        |

| Item Number | Part Number | Description                                                                                                                           | Quantity |
|-------------|-------------|---------------------------------------------------------------------------------------------------------------------------------------|----------|
| 1d          | 100938      | 22" x 1/2" Spinergy LX Blue Spoke                                                                                                     | 1        |
| 1j          | 100942      | 24" x 1/2" Spinergy LX Black Spoke                                                                                                    | 1        |
| 1k          | 100941      | 24" x 1/2" Spinergy LX Red Spoke                                                                                                      | 1        |
| 1l          | 100944      | 24" x 1/2" Spinergy LX White Spoke                                                                                                    | 1        |
| 1m          | 100943      | 24" x 1/2" Spinergy LX Blue Spoke                                                                                                     | 1        |
| 1s          | 118348      | 24" x 1/2" Spinergy LX Black Spoke w/<br>Pneumatic w/ Airless Insert Tire (Black)                                                     | 1        |
| 1t          | 118349      | 24" x 1/2" Spinergy LX Red Spoke w/ Pneumatic<br>w/ Airless Insert Tire (Black)                                                       | 1        |
| 1u          | 118350      | 24" x 1/2" Spinergy LX White Spoke w/<br>Pneumatic w/ Airless Insert Tire (Black)                                                     | 1        |
| 1v          | 118351      | 24" x 1/2" Spinergy LX Blue Spoke w/ Pneumatic<br>w/ Airless Insert Tire (Black)                                                      | 1        |
| 1ab         | 100994      | 25" x 1/2" Spinergy LX Black Spoke                                                                                                    | 1        |
| 1ac         | 100996      | 25" x 1/2" Spinergy LX Red Spoke                                                                                                      | 1        |
| 1ad         | 100998      | 25" x 1/2" Spinergy LX White Spoke                                                                                                    | 1        |
| 1ae         | 100995      | 25" x 1/2" Spinergy LX Blue Spoke                                                                                                     | 1        |
| 1ak         | 100921      | 26" x 1/2" Spinergy LX Black Spoke                                                                                                    | 1        |
| 1al         | 100946      | 26" x 1/2" Spinergy LX Red Spoke                                                                                                      | 1        |
| 1am         | 100945      | 26" x 1/2" Spinergy LX White Spoke                                                                                                    | 1        |
| 1an         | 100947      | 26" x 1/2" Spinergy LX Blue Spoke                                                                                                     | 1        |
| 2           | 102854      | Spinergy LX Spoke Wrench<br><br><i>No charge when ordered with Spinergy wheel or<br/>replacement spoke</i>                            | 1        |
|             | 110119      | Spinergy LX Replacement Spoke Tool Kit<br><br><i>This tool kit contains: Spinergy LX Spoke Wrench and<br/>Spinergy Spoke T-Wrench</i> | 1        |
| 3a          | 101929      | Replacement Spoke Spinergy LX 22", Black<br>(511291-8.75BLA)                                                                          | 1        |
| 3b          | 101930      | Replacement Spoke Spinergy LX 22", Red<br>(511291-8.75RED)                                                                            | 1        |
| 3c          | 101931      | Replacement Spoke Spinergy LX 22", White<br>(511291-8.75WHT)                                                                          | 1        |
| 3d          | 101932      | Replacement Spoke Spinergy LX 22", Blue<br>(511291-8.75BLU)                                                                           | 1        |
| 3i          | 101937      | Replacement Spoke Spinergy LX 24", Black<br>(511291-9.5BLA)                                                                           | 1        |
| 3j          | 101938      | Replacement Spoke Spinergy LX 24", Red<br>(511291-9.5RED)                                                                             | 1        |

| Item Number | Part Number | Description                                                | Quantity |
|-------------|-------------|------------------------------------------------------------|----------|
| 3k          | 101939      | Replacement Spoke Spinergy LX 24", White (511291-9.5WHT)   | 1        |
| 3l          | 101940      | Replacement Spoke Spinergy LX 24", Blue (511291-9.5BLU)    | 1        |
| 3q          | 101945      | Replacement Spoke Spinergy LX 25", Black (511291-9.84BLA)  | 1        |
| 3r          | 101946      | Replacement Spoke Spinergy LX 25", Red (511291-9.84RED)    | 1        |
| 3s          | 101947      | Replacement Spoke Spinergy LX 25", White (511291-9.84WHT)  | 1        |
| 3t          | 101948      | Replacement Spoke Spinergy LX 25", Blue (511291-9.84BLU)   | 1        |
| 3y          | 101953      | Replacement Spoke Spinergy LX 26", Black (511291-10.47BLA) | 1        |
| 3z          | 101954      | Replacement Spoke Spinergy LX 26", Red (511291-10.47RED)   | 1        |
| 3aa         | 101955      | Replacement Spoke Spinergy LX 26", White (511291-10.47WHT) | 1        |
| 3ab         | 101956      | Replacement Spoke Spinergy LX 26", Blue (511291-10.47BLU)  | 1        |

## Tire Component Chart

| Wheel             |                                 | 1 3/8"<br>Street Pneumatic |            |           | 1"<br>High Pressure |            |           | 1 3/8"<br>Euro Trek Knobby |            |           | 1"<br>Schwalbe |          |            |           | 1 3/8"<br>Schwalbe HS-127 |            |           | 2"<br>Mountain |            |           |
|-------------------|---------------------------------|----------------------------|------------|-----------|---------------------|------------|-----------|----------------------------|------------|-----------|----------------|----------|------------|-----------|---------------------------|------------|-----------|----------------|------------|-----------|
| Wheel             | Wheel Part #                    | Tire                       | Inner Tube | Rim Strip | Tire                | Inner Tube | Rim Strip | Tire                       | Inner Tube | Rim Strip | Mara-thon Tire | One Tire | Inner Tube | Rim Strip | Tire                      | Inner Tube | Rim Strip | Tire           | Inner Tube | Rim Strip |
| 22" Spoke         | 116730                          | 116735                     | 100588     | 101908    | 102961              | 100609     | 101908    |                            |            |           |                |          |            |           |                           |            |           | 101475         | 100588     | 101908    |
| 24" Spoke         | 116728                          | 116736                     | 100589     | 100596    | 102958              | 100580     | 100596    | 102982                     | 100589     | 100596    | 100864         | 103176   | 100580     | 100596    | 102920                    | 100589     | 100596    | 101476         | 101188     | 100596    |
| 18" MAXX Spoke    | 109285                          | 101587                     | 101588     | 101137    |                     |            |           |                            |            |           |                |          |            |           |                           |            |           |                |            |           |
| 20" MAXX Spoke    | 108244                          | 116843                     | 100912     | 101137    | 102957              | 102844     | 101137    |                            |            |           |                |          |            |           |                           |            |           |                |            |           |
| 22" MAXX Spoke    | 105135                          | 116735                     | 100588     | 101908    | 102961              | 100609     | 101908    |                            |            |           |                |          |            |           |                           |            |           | 101475         | 100588     | 101908    |
| 24" MAXX Spoke    | 105136                          | 116736                     | 100589     | 100596    | 102958              | 100580     | 100596    | 102982                     | 100589     | 100596    | 100864         | 103176   | 100580     | 100596    | 102920                    | 100589     | 100596    | 101476         | 101188     | 100596    |
| 25" MAXX Spoke    | 107436                          |                            |            |           | 102959              | 101457     | 101452    |                            |            |           | 101357         | 103177   | 101457     | 101452    |                           |            |           | 100166         | 100208     | 101452    |
| 26" MAXX Spoke    | 107437                          | 100636                     | 100590     | 100608    | 102960              | 100582     | 100597    |                            |            |           | 100865         |          | 100582     | 100597    |                           |            |           |                |            |           |
| 20" MAXX Mag      | 111853                          | 116843                     | 100912     | 101137    |                     |            |           |                            |            |           |                |          |            |           |                           |            |           |                |            |           |
| 22" MAXX Mag      | 111854                          | 116735                     | 100588     | 101908    |                     |            |           |                            |            |           |                |          |            |           |                           |            |           | 101475         | 100588     | 101908    |
| 24" MAXX Mag      | 111855                          | 116736                     | 100589     | 100596    |                     |            |           | 102982                     | 100589     | 100596    |                |          |            |           | 102920                    | 100589     | 100596    | 101476         | 101188     | 100596    |
| 24" Superlight    | 101159                          | 111818                     | 100589     | 100596    | 102958              | 100580     | 100596    | 102982                     | 100589     | 100596    | 100864         | 103176   | 100580     | 100596    | 102920                    | 100589     | 100596    |                |            |           |
| 25" Superlight    | 101158                          |                            |            |           | 102959              | 101457     | 101452    |                            |            |           | 101357         | 103177   | 101457     | 101452    |                           |            |           |                |            |           |
| 20" Spinergy Spox | See<br>Spinergy<br>Spox<br>Page | 116843                     | 100912     | 101137    | 102957              | 102844     | 101137    |                            |            |           |                |          |            |           |                           |            |           |                |            |           |
| 22" Spinergy Spox |                                 | 116735                     | 100588     | 101908    | 102961              | 100609     | 101908    |                            |            |           |                |          |            |           |                           |            |           | 101475         | 100588     | 101908    |
| 24" Spinergy Spox |                                 | 116736                     | 100589     | 100596    | 102958              | 100580     | 100596    | 102982                     | 100589     | 100596    | 100864         | 103176   | 100580     | 100596    | 102920                    | 100589     | 100596    | 101476         | 101188     | 100596    |
| 25" Spinergy Spox |                                 |                            |            |           | 102959              | 101457     | 101452    |                            |            |           | 101357         | 103177   | 101457     | 101452    |                           |            |           | 100166         | 100208     | 10145     |
| 26" Spinergy Spox |                                 | 100636                     | 100590     | 100608    | 102960              | 100582     | 100597    |                            |            |           | 100865         |          | 100582     | 100597    |                           |            |           |                |            |           |
| 22" Spingery LX   | See<br>Spinergy<br>LX<br>Page   | 116735                     | 100588     | 101908    | 102961              | 100609     | 101908    |                            |            |           |                |          |            |           |                           |            |           | 101475         | 100588     | 101908    |
| 24" Spingery LX   |                                 | 116736                     | 100589     | 100596    | 102958              | 100580     | 100596    | 102982                     | 100589     | 100596    | 100864         | 103176   | 100580     | 100596    | 102920                    | 100589     | 100596    | 101476         | 101188     | 100596    |
| 25" Spingery LX   |                                 |                            |            |           | 102959              | 101457     | 101452    |                            |            |           | 101357         | 103177   | 101457     | 101452    |                           |            |           | 100166         | 100208     | 101452    |
| 26" Spingery LX   |                                 | 100636                     | 100590     | 100608    | 102960              | 100582     | 100597    |                            |            |           | 100865         |          | 100582     | 100597    |                           |            |           |                |            |           |

See Page 2 for Discontinued Wheels.

### Tire Component Chart (Discontinued Wheels)

| Wheel              |              | 1 3/8"<br>Street Pneumatic |            |           | 1"<br>High Pressure |            |           | 1 3/8"<br>Euro Trek Knobby |            |           | 1"<br>Schwalbe |          |            |           | 1 3/8"<br>Schwalbe HS-127 |            |           | 2"<br>Mountain |            |           |
|--------------------|--------------|----------------------------|------------|-----------|---------------------|------------|-----------|----------------------------|------------|-----------|----------------|----------|------------|-----------|---------------------------|------------|-----------|----------------|------------|-----------|
| Wheel              | Wheel Part # | Tire                       | Inner Tube | Rim Strip | Tire                | Inner Tube | Rim Strip | Tire                       | Inner Tube | Rim Strip | Mara-thon Tire | One Tire | Inner Tube | Rim Strip | Tire                      | Inner Tube | Rim Strip | Tire           | Inner Tube | Rim Strip |
| 18" Ki Spoke       | 200529       | 101587                     | 101588     | 101137    |                     |            |           |                            |            |           |                |          |            |           |                           |            |           |                |            |           |
| 20" Ki Spoke       | 200530       | 100911                     | 100912     | 101137    | 102957              | 102844     | 101137    |                            |            |           |                |          |            |           |                           |            |           |                |            |           |
| 22" Ki Spoke       | 200531       | 100634                     | 100588     | 101908    | 102961              | 100609     | 101908    |                            |            |           |                |          |            |           |                           |            |           | 101475         | 100588     | 101908    |
| 24" Ki Spoke       | 200532       | 111818                     | 100589     | 100596    | 102958              | 100580     | 100596    | 102982                     | 100589     | 100596    | 100864         | 103176   | 100580     | 100596    | 102920                    | 100589     | 100596    | 101476         | 101188     | 100596    |
| 25" Ki Spoke       | 200533       |                            |            |           | 102959              | 101457     | 101452    |                            |            |           | 101357         | 103177   | 101457     | 101452    |                           |            |           | 100166         | 100208     | 101452    |
| 26" Ki Spoke       | 200534       | 100636                     | 100590     | 100608    | 102960              | 100582     | 100597    |                            |            |           | 100865         |          | 100582     | 100597    |                           |            |           |                |            |           |
| 20" 5-Spoke X Core | 101961       | 100911                     | 100912     | 101137    |                     |            |           |                            |            |           |                |          |            |           |                           |            |           |                |            |           |
| 22" 5-Spoke X Core | 101962       | 100634                     | 100588     | 101908    |                     |            |           |                            |            |           |                |          |            |           |                           |            |           |                |            |           |
| 24" 5-Spoke X Core | 100960       | 111818                     | 100589     | 100596    |                     |            |           | 102982                     | 100589     | 100596    |                |          |            |           | 102920                    | 100589     | 100596    |                |            |           |

| Tire Compatibility Chart |                         |                         |                          |                     |                                                    |                   |               |                  |                         |                      |                 |                        |             |
|--------------------------|-------------------------|-------------------------|--------------------------|---------------------|----------------------------------------------------|-------------------|---------------|------------------|-------------------------|----------------------|-----------------|------------------------|-------------|
| Wheel                    |                         | Tire                    |                          |                     |                                                    |                   |               |                  |                         |                      |                 |                        |             |
| Wheel                    | Wheel Part #            | 1 3/8" Street Pneumatic | 1 3/8" Full Profile Poly | 1" Low Profile Poly | 1 3/8" Pneumatic w/ Airless Insert (Tire / Insert) | 1" Solid Sentinel | 1" Solid Shox | 1" High Pressure | 1 3/8" Euro Trek Knobby | 1" Schwalbe Marathon | 1" Schwalbe One | 1 3/8" Schwalbe HS-127 | 2" Mountain |
| 22" Spoke                | 116730                  | 116735                  | 100621                   | 100632              | 116735 / 100503                                    | 103165            | 102718        | 102961           |                         |                      |                 |                        | 101475      |
| 24" Spoke                | 116728                  | 116736                  | 100622                   | 100633              | 116736 / 100504                                    | 102976            | 100640        | 102958           | 102982                  | 100864               | 103176          | 102920                 | 101476      |
| 18" MAXX Spoke           | 109285                  | 101587                  | 101992                   | 100918              | 101587 / 101092                                    |                   |               |                  |                         |                      |                 |                        |             |
| 20" MAXX Spoke           | 108244                  | 116843                  | 100917                   | 102085              | 116843 / 100913                                    | 102974            | 102717        | 102957           |                         |                      |                 |                        |             |
| 22" MAXX Spoke           | 105135                  | 116735                  | 100621                   | 100632              | 116735 / 100503                                    | 103165            | 102718        | 102961           |                         |                      |                 |                        | 101475      |
| 24" MAXX Spoke           | 105136                  | 116736                  | 100622                   | 100633              | 116736 / 100504                                    | 102976            | 100640        | 102958           | 102982                  | 100864               | 103176          | 102920                 | 101476      |
| 25" MAXX Spoke           | 107436                  |                         |                          |                     |                                                    | 102977            | 101453        | 102959           |                         | 101357               | 103177          |                        | 100166      |
| 26" MAXX Spoke           | 107437                  | 100636                  |                          |                     | 100636 / 100505                                    | 102978            | 100884        | 102960           |                         | 100865               |                 |                        |             |
| 12" Mag <sup>1</sup>     | N/A                     | 001424                  | 001426                   |                     | 001425                                             |                   |               |                  |                         |                      |                 |                        |             |
| 16" Mag <sup>1</sup>     | N/A                     | 001159                  | 001161                   | 003453              | 001160                                             |                   |               |                  |                         |                      |                 |                        |             |
| 20" MAXX Mag             | 111853                  | 116843                  | 100917                   | 102085              | 116843 / 100913                                    | 102974            |               |                  |                         |                      |                 |                        |             |
| 22" MAXX Mag             | 111584                  | 116735                  | 100621                   | 100632              | 116735 / 100503                                    | 103165            |               |                  |                         |                      |                 |                        | 101475      |
| 24" MAXX Mag             | 111855                  | 116736                  | 100622                   | 100633              | 116736 / 100504                                    | 102976            |               |                  | 102982                  |                      |                 | 102920                 | 101476      |
| 24" Superlight           | 101159                  | 100635                  |                          |                     |                                                    |                   | 100894        | 102958           | 102982                  | 100864               | 103176          | 102920                 |             |
| 25" Superlight           | 101158                  |                         |                          |                     |                                                    |                   | 101146        | 102959           |                         | 101357               | 103177          |                        |             |
| 20" Spinerger Spox       | See Spinerger Spox Page | 116843                  | 100917                   | 102085              | 116843 / 100913                                    | 102974            |               | 102957           |                         |                      |                 |                        |             |
| 22" Spinerger Spox       |                         | 116735                  | 100621                   | 100632              | 116735 / 100503                                    | 103165            |               | 102961           |                         |                      |                 |                        | 101475      |
| 24" Spinerger Spox       |                         | 116736                  | 100622                   | 100633              | 116736 / 100504                                    | 102976            | 100894        | 102958           | 102982                  | 100864               | 103176          | 102920                 | 101476      |
| 25" Spinerger Spox       |                         |                         |                          |                     |                                                    | 102977            | 101146        | 102959           |                         | 101357               | 103177          |                        | 100166      |
| 26" Spinerger Spox       |                         | 100636                  |                          |                     | 100636 / 100505                                    | 102978            | 100895        | 102960           |                         | 100865               |                 |                        |             |
| 22" Spinerger LX         | See Spinerger LX Page   | 116735                  | 100621                   | 100632              | 116735 / 100503                                    |                   |               | 102961           |                         |                      |                 |                        | 101475      |
| 24" Spinerger LX         |                         | 116736                  | 100622                   | 100633              | 116736 / 100504                                    | 102976            | 100894        | 102958           | 102982                  | 100864               | 103176          | 102920                 | 101476      |
| 25" Spinerger LX         |                         |                         |                          |                     |                                                    | 102977            | 101146        | 102959           |                         | 101357               | 103177          |                        | 100166      |
| 26" Spinerger LX         |                         | 100636                  |                          |                     | 100636 / 100505                                    | 102978            | 100895        | 102960           |                         | 100865               |                 |                        |             |

<sup>1</sup> 12" and 16" Mag Wheels are sold as a combination Wheel and Tire Assembly

See Page 2 for Discontinued Wheels

Effective 2/27/2023

| Tire Compatibility Chart (Discontinued Wheels) |              |                         |                          |                     |                                                    |                   |               |                  |                         |                      |                 |                        |             |
|------------------------------------------------|--------------|-------------------------|--------------------------|---------------------|----------------------------------------------------|-------------------|---------------|------------------|-------------------------|----------------------|-----------------|------------------------|-------------|
| Wheel                                          |              | Tire                    |                          |                     |                                                    |                   |               |                  |                         |                      |                 |                        |             |
| Wheel                                          | Wheel Part # | 1 3/8" Street Pneumatic | 1 3/8" Full Profile Poly | 1" Low Profile Poly | 1 3/8" Pneumatic w/ Airless Insert (Tire / Insert) | 1" Solid Sentinel | 1" Solid Shox | 1" High Pressure | 1 3/8" Euro Trek Knobby | 1" Schwalbe Marathon | 1" Schwalbe One | 1 3/8" Schwalbe HS-127 | 2" Mountain |
| 18" Ki Spoke                                   | 200529       | 101587                  | 101992                   | 100918              | 101587 / 101092                                    |                   |               |                  |                         |                      |                 |                        |             |
| 20" Ki Spoke                                   | 200530       | 100911                  | 100917                   | 102085              | 100911 / 100913                                    | 102974            | 102717        | 102957           |                         |                      |                 |                        |             |
| 22" Ki Spoke                                   | 200531       | 100634                  | 100621                   | 100632              | 100634 / 100503                                    | 103165            | 102718        | 102961           |                         |                      |                 |                        | 101475      |
| 24" Ki Spoke                                   | 200532       | 111818                  | 100622                   | 100633              | 111818 / 100504                                    | 102976            | 100640        | 102958           | 102982                  | 100864               | 103176          | 102920                 | 101476      |
| 25" Ki Spoke                                   | 200533       |                         |                          |                     |                                                    | 102977            | 101453        | 102959           |                         | 101357               | 103177          |                        | 100166      |
| 26" Ki Spoke                                   | 200534       | 100636                  |                          |                     | 100636 / 100505                                    | 102978            | 100884        | 102960           |                         | 100865               |                 |                        |             |
| 20" 5-Spoke X Core                             | 101961       | 100911                  | 100917                   | 102085              | 100911 / 100913                                    |                   | 102717        |                  |                         |                      |                 |                        |             |
| 22" 5-Spoke X Core                             | 101962       | 100634                  | 100621                   | 100632              | 100634 / 100503                                    |                   | 102718        |                  |                         |                      |                 |                        |             |
| 24" 5-Spoke X Core                             | 100960       | 111818                  | 100622                   | 100633              | 111818 / 100504                                    |                   | 100640        |                  | 102982                  |                      |                 | 102920                 |             |

Effective 2/27/2023

| Handrim Hardware Chart |                |                    |                                                              |                                              |                                                                  |                   |                |                  |                  |                    |                  |                      |                  |          |  |
|------------------------|----------------|--------------------|--------------------------------------------------------------|----------------------------------------------|------------------------------------------------------------------|-------------------|----------------|------------------|------------------|--------------------|------------------|----------------------|------------------|----------|--|
| Wheel                  |                | Handrim Connection |                                                              |                                              |                                                                  | Handrim           |                |                  |                  |                    |                  |                      |                  |          |  |
| Wheel                  | Wheel Part No. | Connection Points  | Rivet/Tab                                                    | Spacer*                                      | Screw*                                                           | Aluminum Anodized | Superlight     | Plastic Coated   | Projection       | Ergonomic Standard | Ergonomic LT     | Natural Fit Standard | Natural Fit LT   | Flex Rim |  |
| 22" Spoke              | 116730         | 5 or 6             | 116732 (Used with 5 Rivet)<br><br>100698 (Used with 6 Rivet) | 100653 (Not used w/ Superlight Handrim)      | Aluminum, Plastic Ctd, Projection<br>Ergonomic: 100654           | 101898 (5 Rivet)  |                | 101964 (5 Rivet) | 112824 (6 Rivet) | 113082 (6 Rivet)   | 113085 (6 Rivet) | 200538 (6 Rivet)     | 200201 (6 Rivet) |          |  |
| 24" Spoke              | 116728         | 5 or 6             | Not used w/ Superlight Handrims                              |                                              | Natural Fit: 100835<br><br>Superlight: 100669<br><br>Nut: 100657 | 100975 (5 Rivet)  | 101161 (6 Tab) | 100976 (5 Rivet) | 112825 (6 Rivet) | 113083 (6 Rivet)   | 113806 (6 Rivet) | 100793 (6 Rivet)     | 200202 (6 Rivet) |          |  |
| 18" MAXX Spoke         | 109285         | 3                  | 100698 (Not used w/ Superlight Handrim)                      | 100653 (Not used w/ Superlight Handrim)      | Aluminum, Plastic Ctd, Projection, Ergonomic: 100654             | 100206            |                |                  |                  |                    |                  |                      |                  |          |  |
| 20" MAXX Spoke         | 108244         | 6                  |                                                              |                                              | Natural Fit: 100835                                              | 200536            |                | 112819           |                  |                    |                  |                      |                  |          |  |
| 22" MAXX Spoke         | 105135         |                    |                                                              |                                              | Superlight: 100669                                               | 100560            |                | 112820           | 112824           | 113082             | 113085           | 200538               | 200201           |          |  |
| 24" MAXX Spoke         | 105136         |                    |                                                              |                                              | Nut: 100657                                                      | 200349            | 101161         | 112821           | 112825           | 113083             | 113086           | 100793               | 200202           |          |  |
| 25" MAXX Spoke         | 107436         |                    |                                                              |                                              | 200350                                                           |                   | 112822         | 200548           | 113084           | 113087             | 200539           | 200540               |                  |          |  |
| 26" MAXX Spoke         | 107437         |                    |                                                              |                                              | 200351                                                           |                   | 112823         | 200549           |                  |                    | 100907           | 101454               |                  |          |  |
| 20" MAXX Mag           | 111853         | 6                  |                                                              | Aluminum, Plastic Coated, Projection: 100629 | Aluminum, Plastic Coated, Projection, Ergonomic: 103545          | 200536            |                | 112819           |                  |                    |                  |                      |                  |          |  |
| 22" MAXX Mag           | 111854         |                    |                                                              | Natural Fit & Ergonomic: 101756              | Natural Fit: 101893                                              | 100560            |                | 112820           | 112824           | 113082             | 113085           | 200538               | 200201           |          |  |
| 24" MAXX Mag           | 111855         |                    |                                                              |                                              |                                                                  | 105137            |                | 112821           | 112825           | 113083             | 113086           | 100793               | 200202           |          |  |
| 24" Superlight         | 101159         | 6                  |                                                              |                                              | 100536                                                           | 100754            | 101161         | 100836           |                  | 113080             | 113081           | 100830               | 100828           |          |  |
| 25" Superlight         | 101160         |                    |                                                              |                                              |                                                                  |                   |                | 101160           | 101091           |                    |                  |                      | 101464           | 101460   |  |

\* Spacer and and screw part numbers listed in the chart are for standard handrim mounting. If using close mount handrim mounting on aluminum anodized, plastic coated or projection handrims, use part number 100792 for the spacer, part number 102616 for the screw on Ki Spoke and MAXX Spoke wheels, part number 100666 for the screw on 5-Spoke X Core wheels, or part number 100654 for the screw on MAXX Mag wheels.

Effective 11/27/2023

| Handrim Hardware Chart (Continued) |                         |                    |           |         |                    |                   |            |                |            |                    |              |                      |                |                       |
|------------------------------------|-------------------------|--------------------|-----------|---------|--------------------|-------------------|------------|----------------|------------|--------------------|--------------|----------------------|----------------|-----------------------|
| Wheel                              |                         | Handrim Connection |           |         |                    | Handrim           |            |                |            |                    |              |                      |                |                       |
| Wheel                              | Wheel Part No.          | Connection Points  | Rivet/Tab | Spacer* | Screw*             | Aluminum Anodized | Superlight | Plastic Coated | Projection | Ergonomic Standard | Ergonomic LT | Natural Fit Standard | Natural Fit LT | Flex Rim              |
| 20" Spinerger Spox                 | See Spinerger Spox Page | 6                  |           |         | Ergonomic: 100536  | 103125            |            | 103179         |            |                    |              |                      |                |                       |
| 22" Spinerger Spox                 |                         |                    |           |         | All Others: 100669 | 100827            |            | 100808         |            |                    |              | 100889               | 100888         |                       |
| 24" Spinerger Spox                 |                         |                    |           |         |                    | 100766            | 101161     | 100615         |            |                    |              | 100830               | 100828         |                       |
| 25" Spinerger Spox                 |                         |                    |           |         | Nut: 100657        | 100767            | 101160     | 100765         |            |                    |              | 101464               | 101460         |                       |
| 26" Spinerger Spox                 |                         |                    |           |         |                    | 101477            |            | 101148         |            |                    |              | 200200               | 100950         |                       |
| 22" Spinerger LX                   | See Spinerger LX Page   | 6                  |           |         | Ergonomic: 100536  | 100827            |            | 100808         |            |                    |              | 100889               | 100888         |                       |
| 24" Spinerger LX                   |                         |                    |           |         | All Others: 100669 | 100766            | 101161     | 100615         |            | 113080             | 113081       | 100830               | 100828         | See Spinerger LX Page |
| 25" Spinerger LX                   |                         |                    |           |         |                    | 100767            | 101160     | 100765         |            |                    |              | 101464               | 101460         |                       |
| 26" Spinerger LX                   |                         |                    |           |         | Nut: 100657        | 101477            |            | 101148         |            |                    |              | 200200               | 100950         |                       |

| Handrim Hardware Chart (Discontinued Wheels) |                |                    |           |                                           |                                           |                     |            |                |            |                      |                |          |  |
|----------------------------------------------|----------------|--------------------|-----------|-------------------------------------------|-------------------------------------------|---------------------|------------|----------------|------------|----------------------|----------------|----------|--|
| Wheel                                        |                | Handrim Connection |           |                                           |                                           | Handrim             |            |                |            |                      |                |          |  |
| Wheel                                        | Wheel Part No. | Connection Points  | Rivet/Tab | Spacer*                                   | Screw*                                    | Aluminum Anodized   | Superlight | Plastic Coated | Projection | Natural Fit Standard | Natural Fit LT | Flex Rim |  |
| 18" Ki Spoke                                 | 200529         | 3                  | 100698    | 100653                                    | Aluminum, Plastic Ctd, Projection: 100654 | 100206              |            | 101106         |            |                      |                |          |  |
| 20" Ki Spoke                                 | 200530         | 6                  |           |                                           |                                           | Natural Fit: 100835 | 200536     |                | 200542     |                      |                |          |  |
| 22" Ki Spoke                                 | 200531         |                    |           |                                           |                                           |                     | 100560     |                | 100576     | 100569               | 200538         | 200201   |  |
| 24" Ki Spoke                                 | 200532         |                    |           |                                           |                                           |                     | 200349     |                | 100577     | 200547               | 100793         | 200202   |  |
| 25" Ki Spoke                                 | 200533         |                    |           |                                           | 200350                                    |                     |            | 101870         | 200548     | 200539               | 200540         |          |  |
| 26" Ki Spoke                                 | 200534         |                    |           |                                           | 200351                                    |                     |            | 100578         | 200549     | 100907               | 101454         |          |  |
| 20" 5-Spoke X Core                           | 101961         | 5                  |           | Aluminum, Projection, Natural Fit: 100629 | 100724                                    | 101897              |            | 101963         |            |                      |                |          |  |
| 22" 5-Spoke X Core                           | 101962         |                    |           |                                           |                                           | 101898              |            | 101964         |            |                      |                |          |  |
| 24" 5-Spoke X Core                           | 100960         |                    |           |                                           |                                           | 100975              |            | 100976         | 200546     |                      |                |          |  |
| 24" 5-Spoke X Core                           | 100960         |                    |           | Plastic Ctd: 101756                       | 101893                                    |                     |            |                |            | 100768               | 100769         |          |  |

\* Spacer and screw part numbers listed in the chart are for standard handrim mounting. If using close mount handrim mounting on aluminum anodized, plastic coated or projection handrims, use part number 100792 for the spacer, part number 102616 for the screw on Ki Spoke and MAXX Spoke wheels, part number 100666 for the screw on 5-Spoke X Core wheels, or part number 100654 for the screw on MAXX Mag wheels.

Effective 11/27/2023