

# Clik Wheels - Mag



<p>If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles may be needed. <p>

If Drum Brake is present, see Mag Wheel for Drum Brake page

| Item Number | Part Number | Description  | Quantity | Retail   |
|-------------|-------------|--|----------|----------|
| 1a          | 101961      | Wheel X Core Mag 5-Spoke 20" Black   | 1        | \$147.08 |
| 1b          | 101962      | Wheel X Core Mag 5-Spoke 22" Black   | 1        | \$141.83 |
| 1c          | 100960      | Wheel X Core Mag 5-Spoke 24" Black   | 1        | \$136.58 |
| 3           | 100601      | Wheel Skyway Mag 8-Spoke 22"<br>(H6803XXXXAAA)   | 1        | \$125.02 |
| 4           | 102039      | Hub Cap for 22" 8-Spoke Mag - Discontinued   | 1        | \$3.35   |
| 5           | 200368      | Wheel Complete 24" 6-Spoke Mag/Treaded<br>Poly/Composite HR<br><i>Not available with Heavy Duty Option</i> | 1        | \$173.35 |
| 6a          | 102728      | Wheel X Core Mag 5-Spoke 20" Pink  | 1        | \$159.69 |
| 6b          | 102727      | Wheel X Core Mag 5-Spoke 20" Blue  | 1        | \$159.69 |
| 6c          | 102726      | Wheel X Core Mag 5-Spoke 20" White   | 1        | \$159.69 |
| 6d          | 102725      | Wheel X Core Mag 5-Spoke 20" Red   | 1        | \$159.69 |
| 6e          | 102732      | Wheel X Core Mag 5-Spoke 22" Pink  | 1        | \$170.20 |
| 6f          | 102731      | Wheel X Core Mag 5-Spoke 22" Blue  | 1        | \$170.20 |
| 6g          | 102730      | Wheel X Core Mag 5-Spoke 22" White   | 1        | \$170.20 |
| 6h          | 102729      | Wheel X Core Mag 5-Spoke 22" Red   | 1        | \$170.20 |

# Handrim Hardware Chart

| Wheel                     |                         | Handrim Connection |   |   |   | Handrim             |            |                |            |                      |                |                       |  |  |
|---------------------------|-------------------------|--------------------|---|---|---|---------------------|------------|----------------|------------|----------------------|----------------|-----------------------|--|--|
| Wheel                     | Wheel Part No.          | Connection Points  | Tab                                     | Spacer*   | Screw*                                    | Aluminum Anodized   | Superlight | Plastic Coated | Projection | Natural Fit Standard | Natural Fit LT | Flex Rim              |  |  |
| 18" Ki Spoke <sup>1</sup> | 200529                  | 3                  | 100698                                  | 100653  | Aluminum, Plastic Ctd, Projection: 100654 | 100206              |            | 101106         |            |                      |                |                       |  |  |
| 20" Ki Spoke <sup>1</sup> | 200530                  | 6                  |   |   |   | Natural Fit: 100835 | 200536     |                | 200542     |                      |                |                       |  |  |
| 22" Ki Spoke <sup>1</sup> | 200531                  |                    |   |   |   |                     | 100560     |                | 100576     | 100569               | 200538         | 200201                |  |  |
| 24" Ki Spoke <sup>1</sup> | 200532                  |                    |   |   | 200349                                    |                     |            | 100577         | 200547     | 100793               | 200202         |                       |  |  |
| 25" Ki Spoke <sup>1</sup> | 200533                  |                    |   |   | 200350                                    |                     | 101870     | 200548         | 200539     | 200540               |                |                       |  |  |
| 26" Ki Spoke <sup>1</sup> | 200534                  |                    |   |   | 200351                                    |                     | 100578     | 200549         | 100907     | 101454               |                |                       |  |  |
| 24" Superlight            | 101159                  |                    | 6                                       |   |   | 100536              | 100754     | 101161         | 100836     |                      | 100830         | 100828                |  |  |
| 25" Superlight            | 101160                  |                    |   |   |   |                     |            | 101160         | 101091     |                      | 101464         | 101460                |  |  |
| 20" 5-Spoke X Core        | 101961                  | 5                  |   | Aluminum, Projection, Natural Fit: 100629 Plastic Ctd: 101756 | 100724                                    | 101897              |            | 101963         |            |                      |                |                       |  |  |
| 22" 5-Spoke X Core        | 101962                  |                    |   |   |   | 101898              |            | 101964         |            |                      |                |                       |  |  |
| 24" 5-Spoke X Core        | 100960                  |                    |   |   |   | 100975              |            | 100976         | 200546     |                      |                |                       |  |  |
| 24" 5-Spoke X Core        | 100960                  |                    |   |   | 101893                                    |                     |            |                |            | 100768               | 100769         |                       |  |  |
| 20" Spinerger Spox        | See Spinerger Spox Page | 6                  |   |   | Screw: 100669 Nut: 100657                 | 103125              |            | 103179         |            |                      |                |                       |  |  |
| 22" Spinerger Spox        |                         |                    |   |   |   | 100827              |            | 100808         |            | 100889               | 100888         |                       |  |  |
| 24" Spinerger Spox        |                         |                    |   |   |   | 100766              | 101161     | 100615         |            | 100830               | 100828         |                       |  |  |
| 25" Spinerger Spox        |                         |                    |   |   |   | 100767              | 101160     | 100765         |            | 101464               | 101460         |                       |  |  |
| 26" Spinerger Spox        |                         |                    |   |   |   | 101477              |            | 101148         |            | 200200               | 100950         |                       |  |  |
| 22" Spinerger LX          | See Spinerger LX Page   | 6                  |   |   | Screw: 100669 Nut: 100657                 | 100827              |            | 100808         |            | 100889               | 100888         |                       |  |  |
| 24" Spinerger LX          |                         |                    |   |   |   | 100766              | 101161     | 100615         |            | 100830               | 100828         | See Spinerger LX Page |  |  |
| 25" Spinerger LX          |                         |                    |   |   |   | 100767              | 101160     | 100765         |            | 101464               | 101460         |                       |  |  |
| 26" Spinerger LX          |                         |                    |   |   |   | 101477              |            | 101148         |            | 200200               | 100950         |                       |  |  |
| 18" MAXX Spoke            | 109285                  | 3                  | 100698 (Not used w/ Superlight Handrim) | 100653 (Not used w/ Superlight Handrim)                       | Aluminum, Plastic Ctd, Projection: 100654 | 100206              |            | 101106         |            |                      |                |                       |  |  |
| 20" MAXX Spoke            | 108244                  | 6                  |   |   |   | Natural Fit: 100835 | 200536     |                | 200542     |                      |                |                       |  |  |
| 22" MAXX Spoke            | 105135                  |                    |   |   |   |                     | 100560     |                | 100576     | 100569               | 200538         | 200201                |  |  |
| 24" MAXX Spoke            | 105136                  |                    |   |   | 200349                                    |                     | 101161     | 100577         | 200547     | 100793               | 200202         |                       |  |  |
| 25" MAXX Spoke            | 107436                  |                    |   |   | 200350                                    |                     | 101870     | 200548         | 200539     | 200540               |                |                       |  |  |
| 26" MAXX Spoke            | 107437                  |                    |   |   | 200351                                    |                     | 100578     | 200549         | 100907     | 101454               |                |                       |  |  |

\* Spacer and screw part numbers listed in the chart are for standard handrim mounting. If using close mount handrim mounting on aluminum anodized, plastic coated or projection handrims, use part number 000240 for the spacer, part number 102616 for the screw on Ki Spoke and MAXX Spoke wheels or part number 100666 for the screw on 5-Spoke X Core wheels.

<sup>1</sup> Discontinued

## Tire Component Chart

| Wheel                     |                        | 1 3/8" Street Pneumatic |            |           | 1" High Pressure |            |           | 1 3/8" Euro Trek Knobby |            |           | 1" Schwalbe   |          |            |           | 1 3/8" Schwalbe HS-127 |            |           | 2" Mountain |            |           |        |
|---------------------------|------------------------|-------------------------|------------|-----------|------------------|------------|-----------|-------------------------|------------|-----------|---------------|----------|------------|-----------|------------------------|------------|-----------|-------------|------------|-----------|--------|
| Wheel                     | Wheel Part No.         | Tire                    | Inner Tube | Rim Strip | Tire             | Inner Tube | Rim Strip | Tire                    | Inner Tube | Rim Strip | Marathon Tire | One Tire | Inner Tube | Rim Strip | Tire                   | Inner Tube | Rim Strip | Tire        | Inner Tube | Rim Strip |        |
| 18" Ki Spoke <sup>1</sup> | 200529                 | 101587                  | 101588     | 101137    |                  |            |           |                         |            |           |               |          |            |           |                        |            |           |             |            |           |        |
| 20" Ki Spoke <sup>1</sup> | 200530                 | 100911                  | 100912     | 101137    | 102957           | 102844     | 101137    |                         |            |           |               |          |            |           |                        |            |           |             |            |           |        |
| 22" Ki Spoke <sup>1</sup> | 200531                 | 100634                  | 100588     | 101908    | 102961           | 100609     | 101908    |                         |            |           |               |          |            |           |                        |            |           | 101475      | 100588     | 101908    |        |
| 24" Ki Spoke <sup>1</sup> | 200532                 | 111818                  | 100589     | 100596    | 102958           | 100580     | 100596    | 102982                  | 100589     | 100596    | 100864        | 103176   | 100580     | 100596    | 102920                 | 100589     | 100596    | 101476      | 101188     | 100596    |        |
| 25" Ki Spoke <sup>1</sup> | 200533                 |                         |            |           | 102959           | 101457     | 101452    |                         |            |           | 101357        | 103177   | 101457     | 101452    |                        |            |           | 100166      | 100208     | 101452    |        |
| 26" Ki Spoke <sup>1</sup> | 200534                 | 100636                  | 100590     | 100608    | 102960           | 100582     | 100597    |                         |            |           | 100865        |          | 100582     | 100597    |                        |            |           |             |            |           |        |
| 24" Superlight            | 101159                 | 111818                  | 100589     | 100596    | 102958           | 100580     | 100596    | 102982                  | 100589     | 100596    | 100864        | 103176   | 100580     | 100596    | 102920                 | 100589     | 100596    |             |            |           |        |
| 25" Superlight            | 101158                 |                         |            |           | 102959           | 101457     | 101452    |                         |            |           | 101357        | 103177   | 101457     | 101452    |                        |            |           |             |            |           |        |
| 20" 5-Spoke X Core        | 101961                 | 100911                  | 100912     | 101137    |                  |            |           |                         |            |           |               |          |            |           |                        |            |           |             |            |           |        |
| 22" 5-Spoke X Core        | 101962                 | 100634                  | 100588     | 101908    |                  |            |           |                         |            |           |               |          |            |           |                        |            |           |             |            |           |        |
| 24" 5-Spoke X Core        | 100960                 | 111818                  | 100589     | 100596    |                  |            |           | 102982                  | 100589     | 100596    |               |          |            |           | 102920                 | 100589     | 100596    |             |            |           |        |
| 20" Spinergy Spox         | See Spinergy Spox Page | 100911                  | 100912     | 101137    | 102957           | 102844     | 101137    |                         |            |           |               |          |            |           |                        |            |           |             |            |           |        |
| 22" Spinergy Spox         |                        | 100634                  | 100588     | 101908    | 102961           | 100609     | 101908    |                         |            |           |               |          |            |           |                        |            |           | 101475      | 100588     | 101908    |        |
| 24" Spinergy Spox         |                        | 111818                  | 100589     | 100596    | 102958           | 100580     | 100596    | 102982                  | 100589     | 100596    | 100864        | 103176   | 100580     | 100596    | 102920                 | 100589     | 100596    | 101476      | 101188     | 100596    |        |
| 25" Spinergy Spox         |                        |                         |            |           | 102959           | 101457     | 101452    |                         |            |           | 101357        | 103177   | 101457     | 101452    |                        |            |           |             | 100166     | 100208    | 101452 |
| 26" Spinergy Spox         |                        | 100636                  | 100590     | 100608    | 102960           | 100582     | 100597    |                         |            |           | 100865        |          | 100582     | 100597    |                        |            |           |             |            |           |        |
| 22" Spinergy LX           | See Spinergy LX Page   | 100634                  | 100588     | 101908    | 102961           | 100609     | 101908    |                         |            |           |               |          |            |           |                        |            |           | 101475      | 100588     | 101908    |        |
| 24" Spinergy LX           |                        | 111818                  | 100589     | 100596    | 102958           | 100580     | 100596    | 102982                  | 100589     | 100596    | 100864        | 103176   | 100580     | 100596    | 102920                 | 100589     | 100596    | 101476      | 101188     | 100596    |        |
| 25" Spinergy LX           |                        |                         |            |           | 102959           | 101457     | 101452    |                         |            |           | 101357        | 103177   | 101457     | 101452    |                        |            |           |             | 100166     | 100208    | 101452 |
| 26" Spinergy LX           |                        | 100636                  | 100590     | 100608    | 102960           | 100582     | 100597    |                         |            |           | 100865        |          | 100582     | 100597    |                        |            |           |             |            |           |        |
| 18" MAXX Spoke            | 109285                 | 101587                  | 101588     | 101137    |                  |            |           |                         |            |           |               |          |            |           |                        |            |           |             |            |           |        |
| 20" MAXX Spoke            | 108244                 | 100911                  | 100912     | 101137    | 102957           | 102844     | 101137    |                         |            |           |               |          |            |           |                        |            |           |             |            |           |        |
| 22" MAXX Spoke            | 105135                 | 100634                  | 100588     | 101908    | 102961           | 100609     | 101908    |                         |            |           |               |          |            |           |                        |            |           | 101475      | 100588     | 101908    |        |
| 24" MAXX Spoke            | 105136                 | 111818                  | 100589     | 100596    | 102958           | 100580     | 100596    | 102982                  | 100589     | 100596    | 100864        | 103176   | 100580     | 100596    | 102920                 | 100589     | 100596    | 101476      | 101188     | 100596    |        |
| 25" MAXX Spoke            | 107436                 |                         |            |           | 102959           | 101457     | 101452    |                         |            |           | 101357        | 103177   | 101457     | 101452    |                        |            |           | 100166      | 100208     | 101452    |        |
| 26" MAXX Spoke            | 107437                 | 100636                  | 100590     | 100608    | 102960           | 100582     | 100597    |                         |            |           | 100865        |          | 100582     | 100597    |                        |            |           |             |            |           |        |

<sup>1</sup> Discontinued

## Tire Compatibility Chart

| Wheel                      |                        | Tire                    |                          |                     |  |                   |               |                  |                         |                      |                 |                        |             |
|----------------------------|------------------------|-------------------------|--------------------------|---------------------|--|-------------------|---------------|------------------|-------------------------|----------------------|-----------------|------------------------|-------------|
| Wheel                      | Wheel Part Number      | 1 3/8" Street Pneumatic | 1 3/8" Full Profile Poly | 1" Low Profile Poly | 1 3/8" Pneumatic w/ Airless Insert (Tire / Insert) | 1" Solid Sentinel | 1" Solid Shox | 1" High Pressure | 1 3/8" Euro Trek Knobby | 1" Schwalbe Marathon | 1" Schwalbe One | 1 3/8" Schwalbe HS-127 | 2" Mountain |
| 18" Ki Spoke               | 200529                 | 101587                  | 101992                   | 100918              | 101587 / 101092                                    |                   |               |                  |                         |                      |                 |                        |             |
| 20" Ki Spoke               | 200530                 | 100911                  | 100917                   | 102085              | 100911 / 100913                                    | 102974            | 102717        | 102957           |                         |                      |                 |                        |             |
| 22" Ki Spoke               | 200531                 | 100634                  | 100621                   | 100632              | 100634 / 100503                                    | 103165            | 102718        | 102961           |                         |                      |                 |                        | 101475      |
| 24" Ki Spoke               | 200532                 | 100635                  | 100622                   | 100633              | 100635 / 100504                                    | 102976            | 100640        | 102958           | 102982                  | 100864               | 103176          | 102920                 | 101476      |
| 25" Ki Spoke               | 200533                 |                         |                          |                     |  | 102977            | 101453        | 102959           |                         | 101357               | 103177          |                        | 100166      |
| 26" Ki Spoke               | 200534                 | 100636                  |                          |                     | 100636 / 100505                                    | 102978            | 100884        | 102960           |                         | 100865               |                 |                        |             |
| 24" Superlight             | 101159                 | 100635                  |                          |                     |  |                   | 100894        | 102958           | 102982                  | 100864               | 103176          | 102920                 |             |
| 25" Superlight             | 101158                 |                         |                          |                     |  |                   | 101146        | 102959           |                         | 101357               | 103177          |                        |             |
| 12" Mag <sup>1</sup>       | N/A                    | 001424                  | 001426                   |                     | 001425   |                   |               |                  |                         |                      |                 |                        |             |
| 16" Mag <sup>1</sup>       | N/A                    | 001159                  | 001161                   | 003453              | 001160   |                   |               |                  |                         |                      |                 |                        |             |
| 20" 5-Spoke X Core (Black) | 101961                 | 100911                  | 100917                   | 102085              | 100911 / 100913                                    |                   | 102717        |                  |                         |                      |                 |                        |             |
| 22" 5-Spoke X Core (Black) | 101962                 | 100634                  | 100621                   | 100632              | 100634 / 100503                                    |                   | 102718        |                  |                         |                      |                 |                        |             |
| 24" 5-Spoke X Core (Black) | 100960                 | 100635                  | 100622                   | 100633              | 100635 / 100504                                    |                   | 100640        |                  | 102982                  |                      |                 | 102920                 |             |
| 20" Spinergy Spox          | See Spinergy Spox Page | 100911                  | 100917                   | 102085              | 100911 / 100913                                    | 102974            |               | 102957           |                         |                      |                 |                        |             |
| 22" Spinergy Spox          |                        | 100634                  | 100621                   | 100632              | 100634 / 100503                                    | 103165            |               | 102961           |                         |                      |                 |                        | 101475      |
| 24" Spinergy Spox          |                        | 100635                  | 100622                   | 100633              | 100635 / 100504                                    | 102976            | 100894        | 102958           | 102982                  | 100864               | 103176          | 102920                 | 101476      |
| 25" Spinergy Spox          |                        |                         |                          |                     |  | 102977            | 101146        | 102959           |                         | 101357               | 103177          |                        | 100166      |
| 26" Spinergy Spox          |                        | 100636                  |                          |                     |  | 100636 / 100505   | 102978        | 100895           | 102960                  |                      | 100865          |                        |             |
| 22" Spinergy LX            | See Spinergy LX Page   | 100634                  | 100621                   | 100632              | 100634 / 100503                                    |                   |               | 102961           |                         |                      |                 |                        | 101475      |
| 24" Spinergy LX            |                        | 100635                  | 100622                   | 100633              | 100635 / 100504                                    | 102976            | 100894        | 102958           | 102982                  | 100864               | 103176          | 102920                 | 101476      |
| 25" Spinergy LX            |                        |                         |                          |                     |  | 102977            | 101146        | 102959           |                         | 101357               | 103177          |                        | 100166      |
| 26" Spinergy LX            |                        | 100636                  |                          |                     |  | 100636 / 100505   | 102978        | 100895           | 102960                  |                      | 100865          |                        |             |
| 18" MAXX Spoke             | 109285                 | 101587                  | 101992                   | 100918              | 101587 / 101092                                    |                   |               |                  |                         |                      |                 |                        |             |
| 20" MAXX Spoke             | 108244                 | 100911                  | 100917                   | 102085              | 100911 / 100913                                    | 102974            | 102717        | 102957           |                         |                      |                 |                        |             |
| 22" MAXX Spoke             | 105135                 | 100634                  | 100621                   | 100632              | 100634 / 100503                                    | 103165            | 102718        | 102961           |                         |                      |                 |                        | 101475      |
| 24" MAXX Spoke             | 105136                 | 100635                  | 100622                   | 100633              | 100635 / 100504                                    | 102976            | 100640        | 102958           | 102982                  | 100864               | 103176          | 102920                 | 101476      |
| 25" MAXX Spoke             | 107436                 |                         |                          |                     |  | 102977            | 101453        | 102959           |                         | 101357               | 103177          |                        | 100166      |
| 26" MAXX Spoke             | 107437                 | 100636                  |                          |                     | 100636 / 100505                                    | 102978            | 100884        | 102960           |                         | 100865               |                 |                        |             |

<sup>1</sup> 12" and 16" Mag Wheels are sold as a combination Wheel and Tire assembly