


## (Discontinued) Little Wave Transit



*This Little Wave Transit is discontinued and was used on serial numbers prior to TSTR004553 and LW0004958 (if ordered on chair).*

Item Number	Part Number	Description	Quantity
3	003232	M6 x 55 BHCS Blk Zc	4
4	100658	M6 Nylock Nut Blk Zn	4
5	101708	M6 x 14 BTN SHCS BLKZC w/Patch	4

Item Number	Part Number	Description	Quantity
6	000248	Saddle 1.375"	4