



ACCOUNT INFORMATION

P.O. Number : _____
 Quote Number : _____
 Contact Name: _____
 Email Address: _____
 Phone Number: _____ Fax Number: _____
 User Name or Number: _____
 Billing Name: _____
 Shipping Name: _____
 Ship Address: _____
 City: _____
 Province: _____ Postal Code: _____



CG AIR - THERAPEUTIC SEAT CUSHION

▶ 14" TO 20" WIDE

▼ SINGLE (1) VALVE				
CGA1VSCxxx-x	CG Air Single Valve (Includes a Standard 2-way Stretch Cover)	(E2622)		\$510
CGA1V3Dxxx-x	CG Air Single Valve (Includes a 3D Mesh Cover)	(E2622)		\$564
CGA1VICxxx-x	CG Air Single Valve (Includes an Incontinence Cover)	(E2622)		\$535
▼ DUAL (2) VALVE				
CGA2VSCxxx-x	CG Air Dual Valve (Includes a Standard 2-way Stretch Cover)	(E2624)		\$510
CGA2V3Dxxx-x	CG Air Dual Valve (Includes a 3D Mesh Cover)	(E2624)		\$564
CGA2VICxxx-x	CG Air Dual Valve (Includes an Incontinence Cover)	(E2624)		\$535
▼ SELECT CELL HEIGHT				
	2" (Only available in Single Valve)	3"	4"	
14" x 14" Qty: _____	16" x 16" Qty: _____	18" x 18" Qty: _____	20" x 18" Qty: _____	
14" x 16" Qty: _____	16" x 18" Qty: _____	18" x 20" Qty: _____	20" x 20" Qty: _____	
16" x 14" Qty: _____	18" x 16" Qty: _____	20" x 16" Qty: _____	Other, specify _____	Qty: _____

▶ 22" TO 24" WIDE

▼ SINGLE (1) VALVE				
CGA1VSCxxx-x	CG Air Single Valve (Includes a Standard 2-way Stretch Cover)	(E2623)		\$510
CGA1V3Dxxx-x	CG Air Single Valve (Includes a 3D Mesh Cover)	(E2623)		\$564
CGA1VICxxx-x	CG Air Single Valve (Includes an Incontinence Cover)	(E2623)		\$535
▼ DUAL (2) VALVE				
CGA2VSCxxx-x	CG Air Dual Valve (Includes a Standard 2-way Stretch Cover)	(E2625)		\$510
CGA2V3Dxxx-x	CG Air Dual Valve (Includes a 3D Mesh Cover)	(E2625)		\$564
CGA2VICxxx-x	CG Air Dual Valve (Includes an Incontinence Cover)	(E2625)		\$535
▼ SELECT CELL HEIGHT				
	2" (Only available in Single Valve)	3"	4"	
22" x 18" Qty: _____	22" x 22" Qty: _____	24" x 22" Qty: _____		
22" x 20" Qty: _____	24" x 20" Qty: _____	24" x 24" Qty: _____	Other, specify _____	Qty: _____

▶ 26" + WIDE

▼ SINGLE (1) VALVE				
CGA1VSCxxx-x	CG Air Single Valve (Includes a Standard 2-way Stretch Cover)	(E2623)		\$610
CGA1V3Dxxx-x	CG Air Single Valve (Includes a 3D Mesh Cover)	(E2623)		\$680
CGA1VICxxx-x	CG Air Single Valve (Includes an Incontinence Cover)	(E2623)		\$625
▼ DUAL (2) VALVE				
CGA2VSCxxx-x	CG Air Dual Valve (Includes a Standard 2-way Stretch Cover)	(E2625)		\$610
CGA2V3Dxxx-x	CG Air Dual Valve (Includes a 3D Mesh Cover)	(E2625)		\$680
CGA2VICxxx-x	CG Air Dual Valve (Includes an Incontinence Cover)	(E2625)		\$625
▼ SELECT CELL HEIGHT				
	2" (Only available in Single Valve)	3"	4"	
26" x 20" Qty: _____	28" x 20" Qty: _____	30" x 20" Qty: _____	Other, specify _____	Qty: _____
26" x 22" Qty: _____	28" x 22" Qty: _____	30" x 22" Qty: _____		
26" x 24" Qty: _____	28" x 24" Qty: _____	30" x 24" Qty: _____		



► **REPLACEMENT COMPONENTS**

CGASCxxx-x	Standard 2-way Stretch Cover. Specify Dimensions: W _____ " D _____ " H _____ "	(E2619)	\$130
CGA3Dxxx-x	3D Mesh Cover. Specify Dimensions: W _____ " D _____ " H _____ "	(E2619)	\$150
CGAICxxx-x	Incontinence Cover. Specify Dimensions: W _____ " D _____ " H _____ "	(E2619)	\$130
CGAPK3	CG Air Repair Kit - (Includes 3 alcohol prep pads and 1 small tube of Glue Gel adhesive)		\$31
CGAPK15	CG Air Repair Kit - (Includes 15 alcohol prep pads and 1 large tube of Glue Gel adhesive)		\$149
CUST-133627-01-OF	Custom Transfer Trellis, Black		\$67

CG AIR DELUXE - THERAPEUTIC CUSTOM SEAT CUSHION

► **14" TO 20" WIDE (BUILD AS PER SPECIFICATIONS ON THE NEXT PAGE)**

CGADLX10-20-SC	CG Air Deluxe (Includes a Standard 2-way Stretch Cover) - Below 22" Wide	(E2624)	\$510
CGADLX10-20-3D	CG Air Deluxe (Includes a 3D Mesh Cover) - Below 22" Wide	(E2624)	\$564
CGADLX10-20-IC	CG Air Deluxe (Includes an Incontinence Cover) - Below 22" Wide	(E2624)	\$535

► **22" TO 24" WIDE (BUILD AS PER SPECIFICATIONS ON THE NEXT PAGE)**

CGADLX22-24-SC	CG Air Deluxe (Includes a Standard 2-way Stretch Cover) - 22" to 24" Wide	(E2625)	\$510
CGADLX22-24-3D	CG Air Deluxe (Includes a 3D Mesh Cover) - 22" to 24" Wide	(E2625)	\$564
CGADLX22-24-IC	CG Air Deluxe (Includes an Incontinence Cover) - 22" to 24" Wide	(E2625)	\$535

► **26" + WIDE (BUILD AS PER SPECIFICATIONS ON THE NEXT PAGE)**

CGADLX26+-SC	CG Air Deluxe (Includes a Standard 2-way Stretch Cover) - 26" Wide or more	(E2625)	\$610
CGADLX26+-3D	CG Air Deluxe (Includes a 3D Mesh Cover) - 26" Wide or more	(E2625)	\$680
CGADLX26+-IC	CG Air Deluxe (Includes an Incontinence Cover) - 26" Wide or more	(E2625)	\$625

► **REPLACEMENT COMPONENTS**

CGASCxxx-x	Standard 2-way Stretch Cover. Specify Dimensions: W _____ " D _____ " H _____ "	(E2619)	\$130
CGA3Dxxx-x	3D Mesh Cover. Specify Dimensions: W _____ " D _____ " H _____ "	(E2619)	\$150
CGAICxxx-x	Incontinence Cover. Specify Dimensions: W _____ " D _____ " H _____ "	(E2619)	\$130
CGAPK3	CG Air Repair Kit - (Includes 3 alcohol prep pads and 1 small tube of Glue Gel adhesive)		\$31
CGAPK15	CG Air Repair Kit - (Includes 15 alcohol prep pads and 1 large tube of Glue Gel adhesive)		\$149
CUST-133627-01-OF	Custom Transfer Trellis, Black		\$67

CG AIR DELUXE - INSTRUCTIONS

REFERENCE GUIDE

Cell heights: 2 3 4

Cell not required / base required: X

Cell and base not required: #

Valve(s) and valve extension(s)

- Enter "V" in the corner box(es) for desired location of the valve(s).
- Make sure the number of compartment(s) matches the number of valve(s) indicated in Step 3 on the next page.
- For more than one valve on the same corner, enter "VV" in the box at that corner.



Note: For any valve required at a location other than a corner, specify location in the special instructions section.

Specifying number of compartments (the example opposite shows 2 compartments)

For more than one inflation compartment, use the following procedure:

- Click on the "Tool" tab
- Select the "Comment" tool in the list
- Click on "Drawing Tools" in the bar at the top, and select the "Line" tool
- Use the "Line" tool to identify each inflation compartment

EXAMPLE



Colors on this cushion are to show how a combination of cell heights can be selected to customize the CG Air.
(Blue = 2", Green = 3", Red = 4")

The **CG Air DELUXE** can be customized with a combination of cells (2", 3", and/or 4") to meet the specific needs of a user.

Use this form to order.



CG AIR DELUXE - CUSTOM ORDER MATRIX

		WIDTH - FRONT OF CUSHION														
		2	4	6	8	10	12	14	16	18	20	22	24	26	28	30
DEPTH	2															
	4															
	6															
	8															
	10															
	12															
	14															
	16															
	18															
	20															
	22															
	24															
	26															
	28															
	30															
		BACK OF CUSHION														

STEPS TO FOLLOW TO COMPLETE THE TABLE	Step 1:	Enter cushion dimensions: (Width x Depth) (W) _____" x (D) _____"	Step 5:	If more than one valve was entered in Step 3, make sure that each compartment is identified with lines (As entered in the previous procedure on number of compartments).
	Step 2:	Enter cell Height: (H) _____". For various heights, enter 2, 3, or 4 in the appropriate boxes.	Step 6:	For cells that are not required: a) Enter "X" in boxes where a cell is not required but the base is. b) Enter "#" where neither a cell or the base are required.
	Step 3:	Check the appropriate box for the number of valves: 1 2 3 4 5 6	Step 7:	Check the box for the type of cover required: Standard 2-way Stretch Mesh 3D Incontinence
	Step 4:	For each valve position, enter "V" in the desired corner.		



SIZING CHART

CG Air			OTHER MANUFACTURERS		
WIDTH OR DEPTH	# OF CELLS	SURFACE OF CELLS	WIDTH OR DEPTH	# OF CELLS	SURFACE OF CELLS
10"	5	10"	10"	5	8.75"
11"	5	10"	11"	6	10.50"
12"	6	12"	12"	6	10.50"
13"	6	12"	13"	7	12.25"
14"	7	14"	14"	8	14"
15"	7	14"	15"	8	14"
16"	8	16"	16"	9	15.75"
17"	8	16"	17"	9	15.75"
18"	9	18"	18"	10	17.5"
19"	9	18"	19"	10	17.5"
20"	10	20"	20"	11	19.25"
21"	10	20"	21"	12	21"
22"	11	22"	22"	12	21"
23"	11	22"	23"	13	22.75"
24"	12	24"	24"	13	22.75"
25"	12	24"	25"	14	24.5"

CG AIR BENEFITS

The CG Air cushion is made of medical grade PVC and provides the following benefits:

- Pressure relief on skin and soft tissue
- Better stability
- Better air retention
- Made from fire/flame retardant CAL TB 117-2013
- Better durability
- Best warranty – The warranty covers the cushion for a period of two (2) years and the cover for six (6) months from the date of purchase. See limited warranty on the Owner’s Manual for more details.



SPECIAL INSTRUCTIONS